

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Supplies								
Determine supply status	Power Processor Instructions For Use (IFU), Chapter 8	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Load supplies Aliquot Tips and Aliquot Tubes		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Replace printer labels		Power Processor Instructions For Use (IFU), Appendix D	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Samples								
Load tubes into a rack for samples processing	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Process routine patient samples		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Process priority patient samples		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
View sample information		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
View sample history		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Samples								
Retrieve samples from the Stockyard	Power Processor Instructions For Use (IFU), Chapter 9	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Remap samples	Power Processor Instructions For Use (IFU), Chapter 10	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Priority load samples directly on instrument	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Manage samples in the Error Rack	Power Processor Instructions For Use (IFU), Chapter 7	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Manage samples in the Aliquot Error and Aliquot Storage Rack	Power Processor Instructions For Use (IFU), Chapter 8	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Manage Samples in the Pending Rack	Power Processor Instructions For Use (IFU), Chapter 4	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Unload samples from Outlets	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Unload samples from the Stockyard	Power Processor Instructions For Use (IFU), Chapter 9	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Assign rack number for off line storage	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Maintenance								
Perform Daily Shutdown and Daily Startup, and Pre-Startup Checklist	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform Daily maintenance	Power Processor Instructions For Use (IFU), Appendix D	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform Weekly maintenance		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform Monthly maintenance		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform As-Needed maintenance		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Hardware								
Identify Basic Power Processor system and function	Power Processor Instructions For Use (IFU), Chapter 1	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify LX/DxC Connection module and function	Power Processor Instructions For Use (IFU), Chapter 4	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Instrument Dynamic Connection Modules and function	Power Processor Instructions For Use (IFU), Chapter 5	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Identify Generic Connection module and function	Power Processor Instructions For Use (IFU), Chapter 6	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Hematology module and function	Power Processor Instructions For Use (IFU), Chapter 7	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Aliquot module and function	Power Processor Instructions For Use (IFU), Chapter 8	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Recapper module and function	Power Processor Instructions For Use (IFU), Chapter 9	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Secondary Recapper module and function	Power Processor Instructions For Use (IFU), Chapter 9	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Stockyard module and function	Power Processor Instructions For Use (IFU), Chapter 9	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Error Recovery								
Identify Alarms	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Performs General Error Recovery Procedures	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Locate error code using the Module Error Code Table	Power Processor Instructions For Use (IFU), Appendix B	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Locate sensor using the Module Sensor Diagram	Power Processor Instructions For Use (IFU), Appendix C	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Resolve jammed sample tube carrier	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Utilize PrepLink Event Log	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Configuration								
Configure racks in PrepLink and Line Control	Power Processor Instructions For Use (IFU), Chapter 2	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Configure Test PrepLink and Line Control	Power Processor Instructions For Use (IFU), Chapter 2	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Configure Sample Purge Time in PrepLink	Power Processor Instructions For Use (IFU), Chapter 2	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Miscellaneous								
Bypass a module	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Set LX/DxC from manual mode to automation mode	Power Processor Instructions For Use (IFU), Chapter 4	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Manually centrifuge samples	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Power Processor™ Sample Processing System

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
System Shutdown/Power up	Power Processor Instructions For Use (IFU), Chapter 3	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Saving Setup configurations	Power Processor Instructions For Use (IFU), Chapter 2	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Recovery for LX/DxC Stop Button	Power Processor Instructions For Use (IFU), Appendix A	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			