

DEVELOPMENT OF STANDARDIZED FLOW CYTOMETRY TOOLS TO ASSESS IMMUNOPHENOTYPING IN IMMUNOTOXICOLOGY STUDIES



L. Roqueblave^{1,5}, T. Ramanarivo^{1,5}, C. Le Boulaire¹, H.-J. Egner², H.-W. Vohr³, M. Herbert⁴, F. Montero-Julian¹, C. Fornelli¹
¹Beckman Coulter, Cellular Research & Application Development Centre, Marseille, France; ²Beckman Coulter GmbH, Krefeld, Germany; ³Bayer HealthCare, Wuppertal-Elberfeld, Germany; ⁴Beckman Coulter, Cellular Analysis Business Group, Marseille, France; ⁵L.R. and T.R. have equally contributed to this study.

Abstract

Immunotoxicology evaluations in animal models are mandatory pre-clinical studies of the drug discovery process, as they determine potential adverse health effects that a new drug may have on the immune system of patients. For this reason, the testing strategies need to be reliable and reproducible to accurately assess safety. Among Different Immunotoxicology tests recommended by the ICH Harmonized Tripartite Guideline, Flow Cytometry-based Lymphocyte Phenotyping is a useful assay. To address this, we have developed standardized and automated tools to determine rapidly lymphocyte phenotype in a rat model. By using a cocktail of pre-calibrated labeled antibodies that give minimum emission spectral overlap (FITC/PC7/APC), we enable the addition of a wide variety of PE-labeled antibodies as well as the use of the dead cell exclusion 7-AAD marker. In addition to this flexible 5-color Lymphocyte Immunophenotyping assay, a no wash Red Blood Cell lysing system was optimized for enhanced accuracy. Furthermore, immunostaining and lysing steps can be processed on the Biomek Laboratory Automation Workstation and analysed on the FC 500 cytometer, with auto-setup panels and gating protocols provided. Various examples are shown to demonstrate the usefulness of our methodology. For example, the described technique and analysis were applied for the detection of activated T cells following mitogenic activation with Concanavalin A, or for the exclusion of dead cells in Bone Marrow samples. We have established and compared the results obtained on Rat Peripheral Blood samples processed both manually and on the automation platform, and shown no significant difference. In addition, the sensitivity and specificity of our method to know immuno-suppressive drugs and/or immunogenic agents has been evaluated.

Introduction

All new human pharmaceuticals should be evaluated for the potential to produce immunotoxicity, which means immunosuppression, immunoenhancement, allergenicity or drug-specific autoimmunity. Nonclinical / preclinical testing can help to reveal the first two of these adverse effects, except that, apart from the so-called "standard toxicity studies", no standardized practical approaches are currently available; methods to raise data in practice are not described! Here we take the flow cytometry-based lymphocyte subset phenotyping of rat biological samples as one example of tests recommended [by the International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use (ICH)] in the "Immunotoxicity Studies for Human Pharmaceuticals" ICH S8 harmonized tripartite guideline. Usually, these are performed by home brewing tests mixing reagent components and analysis instruments from various origins that may prove difficult to reconcile one to the other, and to validate properly. We describe a set of reagents and adapted pre analysis tools that have been designed and developed to fit with the Biomek NX Liquid Handling system and the FC 500 MPL Flow Cytometer from Beckman Coulter. By design, this methodology aims at streamlining the preparation steps, with easy procedures, whatever the sample to be used, and whatever the lymphocyte sub-family level we want to reach. It is also aimed at opening the process to automation, while giving the freedom to manually perform smaller testing series with the same level of confidence in the results. We performed different evaluations of this methodology to demonstrate its flexibility linked to high robustness.

Materials and Methods

Animals and Study Protocol

Tests were performed mainly with Wistar rats and more rarely with Sprague Dawley or Brown Norway strains. Female (F) and Male (M) rats were assigned to a control group and 3 experimental groups as follows: Group 1: 10 F & 10 M treated with PBS; Group 2: 10 F & 10 M treated with 10⁶ SRBCs per animal by i.v. on day 2; Group 3: 10 F & 10 M treated with Cyclophosphamide (CY 40 mg/kg by i.p. on day 1) + SRBCs (10⁶ per animal by i.v. on day 2); Group 4: 10 F & 10 M treated with Cyclophosphamide (CY 80 mg/kg by i.p. on day 1) + SRBCs (10⁶ per animal by i.v. on day 2). Hematological investigations were performed on blood and organ samples were taken at terminal sacrifice on day 7. Lymphocyte enumeration following rat treatment with CY is known to be reflective of immunosuppression.

Monoclonal Antibodies (mAbs)


The following mouse mAbs were used to develop the three color combinations.
 • IOtest[®] Anti-Rat CD3-FITC / CD4-PC7 / CD8-APC (PN A32909); anti-rat CD3, clone 1F4, IgM; anti-rat CD4, clone OX38, IgG2a; anti-rat CD8, clone OX8, IgG1.
 • IOtest[®] Anti-Rat CD3-FITC / CD45RA-PC7 / CD161-APC (PN A32910); anti-rat CD3, clone 1F4, IgM; anti-rat CD45RA, clone OX33, IgG1; anti-rat CD8, clone 10/78, IgG1.

Sample Preparation (manual)

The cell population of interest is stained with monoclonal antibodies and the red blood cells in each sample are lysed in a No wash "Fix-and-Lyse" procedure based on the combination of VersaLys[™] Lysing Solution (PN IM3648) and IOtest 3 Fixative Solution (PN IM3515). Briefly, 25 µL of biological sample (whole blood, or spleen cell preparation, or lymph node preparation, or broncho-alveolar lavage preparation) are incubated with 25 µL of mAbs combination for 20 minutes at room temperature (RT), protect from light. The Fix & Lyse mixture is then added to the samples and incubated for 10 minutes at RT, protect from light. The stained samples will be kept at 2-8°C until use (≤ 2hrs following the lysis).




Flow Cytometry

The flow cytometer must be equipped to detect Forward Scatter, Side Scatter and three fluorescence channels allowing the analysis of FITC-, PC7-, and APC-conjugated antibodies (respectively 504 – 541 / 750 – 810 / 650 – 680 nm maximal peak emission). APC conjugates require an exciting source of 633 nm (He-Ne laser) or 635 nm (Red diode laser).
 1- Prepared samples are analyzed on FC 500 / Multi Carousel Loader (MCL) or FC 500 / Multi Plate Loader (MPL) using respectively Beckman Coulter Cytomics™ CXP™ and MXP™ analysis software. In all experiments, a mixture of Flow-Set™ Fluorospheres (PN 6607007) and Flow-Set 675 Fluorospheres (PN 6607120) were used to set up the PMT values while compensations were set by using Cyto-Comp Cells (PN 6607023) stained separately with either CD45-FITC (PN IM0782), CD45-PC7 (PN IM3548) or CD45-APC (PN IM2473).
 2- Analysis: a large gate is drawn around lymphocytes on an FS/SS histogram and the CD3⁺ and CD3⁻ subpopulations were analyzed as described below.



Sample Preparation using the Biomek[®] NX Span-8

In some experiments, the sample preparations (blood cell suspensions + antibodies mixing in view of flow cytometric immunophenotyping) have been performed with a Beckman Coulter Biomek NX Span-8 Laboratory Automated Workstation, using a series of eight 1000-µL teflon™ needles to distribute biologic samples and reagents. The Span-8 option allows for liquid level sensing. We have optimized the method to provide a walk-away preparation system (follow-up of reagent and sample levels, mixing step optimized to avoid sedimentation of blood cells, specific aspiration process to eliminate risk of contaminations using blowout air and leading air gap, etc.). Mixing has been performed by multi aspiration / dispense of fluids in microtubes (Micro, Ref. M32022) arranged on a 96-socket carrier rack. After preparation, the rack has been directly transferred on the plate loader of an FC 500 MPL. A fully integrated installation of the FC 500 MPL plate loader can be proposed as optional feature, with a direct link with the Biomek NX Span-8 working area.



Discussion

The data presented here demonstrates that the methodology is suitable for the immunophenotyping of T-, B- and NK-cell populations in various rat biological samples. The main features of this methodology are:
 1- Robustness, since straightforward results can be obtained with whole blood samples, spleen or lymph node preparations, and bone marrow samples. 2- Flexibility, since you can perform the sample preparation into regular lysing tubes and into small-sized tubes that fit into rack with microplate format. Also, from an analytical point of view, you can add to the three color combinations a choice of viability dye (7-AAD), gating tool (CD45-PE), or cell-status / -activation marker (e.g. CD25-PE). You might even combine drug ins. would you need to perform at the same time dead cell exclusion and CD45-gating. 3- Easeiness since a no wash lysis procedure is involved and an automated preparation process is possible using the Biomek NX S-8, and can be fully synchronized with the flow cytometric acquisition on the FC 500 MPL. Further evaluations are still to be performed, one of which being to explore absolute counting of the lymphocyte subsets using Flow-Count™ Fluorospheres, a feature that would allow a direct determination of absolute counts, together with the full standardization of this methodology.

References

- Cooper-Hannan, R. et al., 1999, "The Principles of Good Laboratory Practices: Application to In Vitro Toxicology Studies. The Report and Recommendations of ECVAM Workshop 37", ATLA, 27, 539-577.
- European Medicines Agency, 2006, ICH S8, Immunotoxicity studies for human pharmaceuticals, note for guidance on Immunotoxicity studies for human pharmaceuticals, (EMEA/CHMP/167235/2004), 1-13.
- Criswell, K.A., Blewins, M.R., Zielinski, D., Zandes, J.C., "Comparison of flow cytometric and manual bone marrow differentials in Wistar rats", 1998, Cytochemistry, 32, 9-17.
- Dean, J.H., House, R.V., Luster, M.I., 2001, "Immunotoxicology: Effects of, and Response to, Drugs and Chemicals", in Principles and Methods of Toxicology, 4th Edition, Wallace Hayes, Ed., 1415-1450.
- Tanaka, T., Masuko, T., Yagita, H., Tamura, T., Hashimoto, Y., "Characterization of a CD3-like rat T cell surface antigen recognized by a monoclonal antibody", 1989, J. Immunol., 142, 2791-2795.
- Barclay, A.N., "The localization of populations of lymphocytes defined by monoclonal antibodies in rat lymphoid tissues", 1981, Immunology, 42, 593-600.
- Franch A, Castelleto, C, Pelegri C, Tolosa E, Castell M., "Blood B, T, CD4+ and CD8+ lymphocytes in female Wistar rats", 1993, Ann Hematol., 67, 115-118.
- Thomas, M.L., Green, J.R. "Molecular nature of the W3/25 and MRC OX-8 marker antigens for T lymphocytes: comparisons with mouse and human antigens", 1983, Eur. J. Immunol., 13, 855-858.
- Brideau, R.J., Carter, P.B., McMaster, W.R., Mason, D.W., Williams, A.F., "Two subsets of rat T lymphocytes defined with monoclonal antibodies", 1980, Eur. J. Immunol., 10, 609-615.