BECKMAN COULTER IMMUNOCHEMISTRY SYSTEMS**Beckman Coulter CAL 2 (Calibrator 2)****REF 663940****TABLE OF CONTENTS**

English	2
Français (France) (fr-fr)	5
Deutsche (de-de)	9
Italiano (it-it)	13
Español (es-es)	17
Português (Portugal) (pt-pt)	21
Dansk (da-dk)	25
Svenska (sv-se)	28
Norsk (nb-no)	31
Ελληνικά (el-gr)	35
日本語 (ja-jp)	39
中文 (中国) (zh-cn)	42
Lietuvių (lt-lt)	45
Magyar (hu-hu)	49
Polski (pl-pl)	53
Čeština (cs-cz)	57
Slovák (sk-sk)	61
한국어 (ko-kr)	65
Türkçe (tr-tr)	68
Română (ro-ro)	72
Українська (uk-ua)	76
Português (Brasil) (pt-br)	80
Nederlands (nl-nl)	84
REFERENCES	88

BECKMAN COULTER IMMUNOCHEMISTRY SYSTEMS

Beckman Coulter CAL 2 (Calibrator 2)

REF 663940

For *In Vitro* Diagnostic Use

For Professional Use Only

Rx Only

INTENDED USE

CAL 2 (Calibrator 2), when used in conjunction with Beckman Coulter alpha1-antitrypsin (AAT), ceruloplasmin (CER), properdin factor B (PFB), alpha2-macroglobulin (AMG), and antithrombin III (AT3) reagents, is intended for use on IMMAGE Systems for the calibration of these reagents.

SUMMARY

The calibrator is derived from fresh frozen human plasma that has been defibrinated and processed. Assay of CAL 2 provides a response value that is utilized for the adjustment of pre-programmed calibration curves from which AAT, CER, PFB, AMG, and AT3 concentration values are determined in test specimens.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

STORAGE AND STABILITY

The calibrator stored at +2°C to +8°C is stable until the expiration. After reconstitution, Calibrator 2 is stable for 4 months if stored at -20°C to -15°C and for 1 month if stored at +2°C to +8°C, or until expiration date if sooner.

NOTICE

Adverse storage conditions of this material may cause erroneous test results.

CONTENTS

IMMAGE Component:	QUANTITY
Beckman Coulter Calibrator 2 (CAL 2)	4 x 1 mL vials
Calibrator 2 Bar Code Card	1
Calibrator 2 Bar Code Strips	2
Value Assignment Sheet	1

NOTICE

Lot number on the CAL 2 vials and cards must be the same.

REACTIVE INGREDIENTS

CALIBRATOR 2, 1.0 mL: Human Based Material

CAUTION

Because this product is of human source origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HB_sAg. Because no test method can offer complete assurance that HIV, Hepatitis B virus, and Hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

⚠ CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

INSTRUCTIONS FOR USE

1. Allow contents of the bottle to come to room temperature.
2. To reconstitute, add 1.0 milliliter of distilled water and allow to stand at room temperature for 30 minutes.
3. Mix the contents by gently inverting the bottle.
4. Refer to the appropriate Operations Manuals and reagent Instructions for Use for the calibration procedures.

ADDITIONAL BAR CODE INSTRUCTIONS

1. Locate the calibrator bar code label provided with the calibrator.
2. Place the appropriate label on an empty 13 x 100 mm or 16 x 100 mm test tube. The labeled test tube should be saved for re-use.
3. Place the appropriate calibrator in a sample cup.
4. Place the appropriate cup in the labeled test tube.
5. Repeat steps 1-5 for any additional calibrators. Continue placing calibrators into the same sample rack.
6. Open the sample compartment lid.
7. Place the rack(s) containing calibrators in positions A and B.
8. Close the sample compartment lid.

ASSIGNED VALUES

IMMAGE values are on the value assignment sheet.

Multiple instruments using representative samples from this lot of calibrator established the assigned values. Measurand values are specific to the assay methodologies of the IMMAGE reagents. Values assigned by other Methodologies may be different. Such differences, if present, may be caused by inter-method bias.

LIMITATIONS

1. Beckman Coulter CAL 2 values are specific for Beckman Coulter Immunochemistry Analyzers only when using Beckman Coulter Reagents. Use of CAL 2 with other reagents or test procedures is not recommended.
2. Dust particles or other particulate matter in the reaction solution may result in extraneous light-scattering signals for the IMMAGE resulting in variable calibration rate units.

ADDITIONAL INFORMATION

Complete instructions for the use of the calibrator are contained in the Beckman Coulter Immunochemistry Systems Operating Manuals. Copies of these manuals are available from Beckman Coulter, Inc., Brea, CA 92821.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

Each measurand (analyte) in this calibrator is traceable to the appropriate reference materials as listed below. Traceability process is based on ISO 17511.

Measurands Higher order Reference Materials

AAT and AMG ERM-DA470k/IFCC^{1,2,3*}
PFB, CER and AT3 Manufacturer's working calibrator

*. ERM is a registered trademark of EC-JRC-IRMM.

REVISION HISTORY

Revision AE

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AF

Added new language requirement: Brazilian Portuguese.

Revision AG

Added new language requirement: Dutch, and Slovak. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.



Revision AH

Updates to comply with IVDR requirements per Beckman Coulter Global Labeling Policy.

Revision AJ

Added Ukrainian Language

Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	CE Mark		Consult Instructions for Use
	Authorized Representative in the European Community		Date of Manufacture
	Caution		Biological risks
	Calibrator, After Reconstitution		Bar Code Card
	Value Assignment Sheet		Bar Code Strip
			

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REFERENCES

1. Ingrid Zegers, Thomas Keller et al. "Characterization of the New Serum Protein Reference Material ERM-DA470k/IFCC: Value Assignment by Immunoassay," *Clin. Chem.*, 56:12 1880-1888(2010)
2. JRC-IRMM Certification Report "Certification of proteins in the human serum Certified Reference Material ERM-DA470k/IFCC."
3. Whicher, J. T., et al., "New International Reference Preparation for Proteins in Human Serum (RPPHS)," *Clin. Chem.*, 40/6, 934-938 (1994).

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Revisie AG
Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.