## SYNCHRON Systems

### Cal 1

**REF** 468405

#### TABLE OF CONTENTS

English . . . . .	2
Français (FR) . . . . .	4
Deutsch (DE) . . . . .	7
Italiano (IT) . . . . .	10
Español (ES) . . . . .	13
Português Portugal (PT-PT) . . . . .	16
Dansk (DA) . . . . .	19
Svenska (SV) . . . . .	21
Norsk (NO) . . . . .	23
Ελληνικά (EL) . . . . .	25
日本語 (JP) . . . . .	28
中文 (ZH-CN) . . . . .	30
Lietuviškai (LT) . . . . .	32
Magyar (HU) . . . . .	34
Polski (PL) . . . . .	37
Čeština (CZ) . . . . .	40
한국어 (KO) . . . . .	43
Türkçe (TR) . . . . .	45
Български (BG) . . . . .	47
Română (RO) . . . . .	50
Србија (SR) . . . . .	52
Tiếng Việt (VI) . . . . .	55
REFERENCES . . . . .	58

---

# SYNCHRON Systems

## Cal 1

**REF** 468405

For *In Vitro* Diagnostic Use

Rx Only

### INTENDED USE

The Beckman Coulter DxC Calibrator 1 (CAL 1), used in conjunction with DxC reagents, is intended for the calibration of immunoprotein tests on UniCel DxC SYNCHRON Systems.

### SUMMARY

Assay of DxC CAL 1 provides a response value that is utilized for the adjustment of pre-programmed calibration curves from which IgA, IgG, IgM, Transferrin, C3, C4, and Haptoglobin concentration values are determined in test specimens. The UniCel DxC SYNCHRON Systems utilize immunoturbidimetry for the quantitative measurement of these analytes. The calibrator is prepared from processed human serum stabilized with < 0.1% Sodium Azide (w/w).

### ⚠ CAUTION

Because this product is of human source origin, it should be handled as though capable of transmitting infectious disease. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for the antibodies to HIV and HCV and nonreactive for HB<sub>s</sub>Ag. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C or other infectious agents are absent, this material and all patient samples should be handled as though capable of transmitting infectious disease. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2 guidelines.

### GHS HAZARD CLASSIFICATION

Not classified as hazardous

**SDS** Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

### STORAGE AND STABILITY

CAL 1 is stable until the expiration date printed on the calibrator bottle, if stored capped in the original containers, at +2°C to +8°C. DO NOT FREEZE.

### CONTENTS

Each kit contains the following items:

COMPONENT	QUANTITY
3 mL of DxC CAL 1	4
Value Assignment Sheet	1

### ⚠ CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

### MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

Calibration Diskette P/N 468408 (lot specific)

## INSTRUCTIONS FOR USE

1. Gently mix the calibrator by inversion prior to opening.
2. Transfer the calibration serum into the sample cup.
3. Immediately recap the bottle tightly and return the container to the appropriate storage conditions.
4. For detailed calibration instructions, refer to the SYNCHRON DxC Clinical System Operating Manual.

## LIMITATIONS

DxC Calibrator 1 should be used only in conjunction with UniCel DxC SYNCHRON Systems and DxC reagents. Adverse storage conditions of CAL 1 may cause erroneous test results.

## ADDITIONAL INFORMATION

Complete instructions for the use of this calibrator are contained in the Beckman Coulter Clinical Systems Operating Manuals. Copies of these manuals are available from Beckman Coulter, Inc., Brea, CA 92821.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see [www.beckmancoulter.com/patents](http://www.beckmancoulter.com/patents).

## SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

## TRACEABILITY

Each measurand (analyte) in this calibrator is traceable to the IFCC reference preparation for plasma proteins, ERM-DA470k/IFCC.<sup>123</sup> Traceability process is based on prEN ISO 17511.

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the DxC reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by intermethod bias.

System specific values are found on the enclosed value assignment sheet.

## Revision History

### Revision AG

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.


### Revision AH

Added new language requirement: Bulgarian, Serbian, and Vietnamese. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

### Revision AJ

Moved Calibration Diskette information from CONTENTS to MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT.

## Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	CE Mark		Consult Instructions for Use
	Authorized Representative in the European Community		Date of Manufacture
	Caution		Biological risks
	Calibrator		Value Assignment Sheet
	Calibration Diskette		
			

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## REFERENCES

1. Ingrid Zegers, Thomas Keller et al. "Characterization of the New Serum Protein Reference Material ERM-DA470k/IFCC: Value Assignment by Immunoassay," *Clinical Chemistry*, 56:12 1880-1888(2010)
2. JCR-IRMM Certification Report "Certification of proteins in the human serum Certified Reference Material ERM-DA470k/IFCC."
3. Johnson AM. A new international reference preparation for proteins in human serum (RPPHS). *Arch Pathol Lab Med* 1993; 117:29-31.

Beckman Coulter Eurocenter SA  
22, rue Juste-Olivier, Case Postale 1044,  
CH-1260 Nyon 1, Switzerland.  
Telephone: +41 (0)22 365 36 11

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### Phiên bản AG
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.