Immunochemistry Systems CAL 3 (Calibrator 3)

REF 449730

TABLE OF CONTENTS

| | |
|--------------------------------------|----|
| English | 2 |
| Français (FR) | 5 |
| Deutsch (DE) | 8 |
| Italiano (IT) | 11 |
| Español (ES) | 14 |
| Português Portugal (PT-PT) | 17 |
| Dansk (DA) | 20 |
| Svenska (SV) | 23 |
| Norsk (NO) | 26 |
| Ελληνικά (EL) | 29 |
| 日本語 (JP) | 33 |
| 中文 (ZH-CN) | 36 |
| Lietuviškai (LT) | 39 |
| Magyar (HU) | 42 |
| Polski (PL) | 45 |
| Čeština (CZ) | 48 |
| Slovenčina (SK) | 51 |
| 한국어 (KO) | 54 |
| Türkçe (TR) | 57 |
| Română (RO) | 60 |
| Português Brasil (PT-BR) | 63 |
| Nederlands (NL) | 66 |
| REFERENCES | 69 |

Immunochemistry Systems CAL 3 (Calibrator 3)

REF 449730

For *In Vitro* Diagnostic Use
Rx Only

INTENDED USE

CAL 3 (Calibrator 3), when used in conjunction with Beckman Coulter albumin (ALB) and prealbumin (PAB) reagents, is intended for use on IMMAGE Systems for the calibration of these reagents.

SUMMARY

The calibrator is derived from fresh frozen human plasma that has been defibrinated and processed. Assay of CAL 3 provides a response value that is utilized for the adjustment of pre-programmed calibration curves from which ALB and PAB concentration values are determined in test specimens.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

STORAGE AND STABILITY

The calibrator is stable until the expiration date printed on the calibrator bottle if stored capped in the original container at +2°C to +8°C.

NOTICE

Adverse storage conditions of this material may cause erroneous test results.

CONTENTS

| | |
|--------------------------------------|-----------------|
| Beckman Coulter Calibrator 3 (CAL 3) | 4 x 3 mL vials |
| IMMAGE Component: | QUANTITY |
| Calibrator 3 Bar Code Card | 1 |
| Calibrator 3 Bar Code Strips | 2 |
| Value Assignment Sheet | 1 |

NOTICE

Lot number on the CAL 3 vials and cards must be the same.

CAUTION

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HB_sAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled as specified in the Centers for Disease Control's Biosafety Level 2 Guidelines.

⚠ CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

INSTRUCTIONS FOR USE

1. Allow contents of the bottle to come to room temperature.
2. Mix the contents by gently inverting the bottle.
3. Refer to the appropriate Operations Manuals, Chemistry Information Manuals, or Reagent Package Insert for the calibration procedures.

ADDITIONAL IMMAGE BAR CODE INSTRUCTIONS

1. Locate the calibrator bar code label provided with the calibrator.
2. Place the appropriate label on an empty 13 x 100 mm or 16 x 100 mm test tube. The labeled test tube should be saved for re-use.
3. Place the appropriate calibrator in a sample cup.
4. Place the appropriate cup in the labeled test tube.
5. Place the tube into a sample rack.
6. Repeat steps 1-5 for any additional calibrators. Continue placing calibrators into the same sample rack.
7. Open the sample compartment lid.
8. Place the rack(s) containing calibrators in positions A and B.
9. Close the sample compartment lid.

ASSIGNED VALUES

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the IMMAGE systems. Values established by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

IMMAGE values can be printed from the console.

LIMITATIONS

1. Beckman Coulter CAL 3 values are specific for Beckman Coulter Immunochemistry Analyzers only when using Beckman Coulter Reagents. Use of CAL 3 with other reagents or test procedures is not recommended.
2. Dust particles or other particulate matter in the reaction solution may result in extraneous light-scattering signals for the IMMAGE resulting in variable calibration rate units.

ADDITIONAL INFORMATION

Complete instructions for the use of the calibrator are contained in the Beckman Coulter Immunochemistry Systems Operating Manuals. Copies of these manuals are available from Beckman Coulter, Inc., Brea, CA 92821.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

Each measurand (analyte) in this calibrator is traceable to the IFCC reference preparation for plasma proteins, ERM-DA470k/IFCC.^{1,2,3*} Traceability process is based on prEN ISO 17511. The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the IMMAGE reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias. IMMAGE values can be printed from the console.

* ERM is a registered trademark of EC-JRC-IRMM.

Revision History

Revision AE

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AF

Added new language requirement: Brazilian Portuguese.

Revision AG

Added new language requirement: Dutch, and Slovak. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

Symbols Key

| | | | |
|---|---|--|--|
|  | Catalogue Number |  | In Vitro Diagnostic |
|  | Contents |  | Temperature limit |
|  | Manufacturer |  | Expiration Date |
|  | Batch code |  | Safety Data Sheet |
|  | CE Mark |  | Consult Instructions for Use |
|  | Authorized Representative in the European Community |  | Date of Manufacture |
|  | Caution |  | Biological risks |
|  | Calibrator |  | Value Assignment Sheet |
|  | Bar Code Strip |  | Bar Code Card |
|  | |  | Made in USA of US and Foreign Components |

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REFERENCES

1. Ingrid Zegers, Thomas Keller et al. "Characterization of the New Serum Protein Reference Material ERM-DA470k/IFCC: Value Assignment by Immunoassay," *Clin. Chem.*, 56:12 1880-1888(2010)
2. JCR-IRMM Certification Report "Certification of proteins in the human serum Certified Reference Material ERM-DA470k/IFCC."
3. Whicher, J. T., et al., "New International Reference Preparation for Proteins in Human Serum (RPPHS)," *Clin. Chem.*, 40/6, 934-938 (1994).

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Revisie AG
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.