LIQUID COMPREHENSIVE CHEMISTRY CONTROL SERUM LEVEL 1, 2, 3

REF 657365

TABLE OF CONTENTS

English	2
Français (FR)	6
Deutsch (DE)	10
Italiano (IT)	14
Español (ES)	18
Português Portugal (PT-PT)	22
Dansk (DA)	26
Svenska (SV)	29
Norsk (NO)	33
Ελληνικά (EL)	37
日本語 (JP)	41
中文 (ZH-CN)	44
Lietuviškai (LT)	47
Magyar (HU)	51
Polski (PL)	55
Čeština (CZ)	59
한국어 (KO)	63
Türkçe (TR)	66
Български (BG)	70
Română (RO)	74
Србија (SR)	78
Tiếng Việt (VI)	82
REFERENCES	86

LIQUID COMPREHENSIVE CHEMISTRY CONTROL SERUM LEVEL 1, 2, 3

REF 657365

**For In Vitro Diagnostic Use
Rx Only**

INTENDED USE

SYNCHRON Controls are made from stabilized human serum and are designed to monitor the performance of SYNCHRON Systems in the clinical laboratory. SYNCHRON Control Levels 1, 2, and 3 bear a quantitative relationship to one another; Level 2 is manufactured by combining equal quantities of Levels 1 and 3. The use of three levels enables the laboratorian to monitor changes in calibration and linearity along with analytical error and imprecision.

SUMMARY

SYNCHRON control serum is prepared from fresh frozen human plasma that has been defibrinated and, then, it is stabilized with ethylene glycol.^{1,2,3,4} The stabilizing effects of ethylene glycol are threefold. The high osmolality minimizes bacterial growth while its antioxidant property stabilizes oxygen-labile constituents. Additionally, the presence of ethylene glycol leads to freezing point depression, which allows the control material to remain in liquid state at normal freezer temperature; i.e., between -15°C and -20°C. The preparation of these stabilized controls in liquid form eliminates errors commonly associated with filling, drying, and the reconstitution of lyophilized control products.^{5,6,7,8,9}

Raw materials used in these products are from purified preparations. Amylase, Lipase, ALT, AST, and γ GT are from a porcine source. LDH is from an avian source, CK from a simian source, and alkaline phosphatase from a bovine source.

CAUTION

Because this product is of human source origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HB_sAg. Because no test method can offer complete assurance that HIV, Hepatitis B virus, and Hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

GHS HAZARD CLASSIFICATION

Comprehensive Chemistry
Control Serum (Human
Origin)
Level 1 WARNING

H302
P264

P301+P312

Harmful if swallowed.
Wash hands thoroughly after
handling.
IF SWALLOWED: Call a POISON
CENTER or doctor/physician if you
feel unwell.
Ethylene Glycol 30 - 40%

Comprehensive Chemistry
Control Serum (Human
Origin)
Level 2 WARNING

H302
EUH208
P264

P301+P312

Harmful if swallowed.
May produce an allergic reaction.
Wash hands thoroughly after handling.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
γ-Glutamyltransferase < 0.1%
Ethylene Glycol 30 - 40%

Comprehensive Chemistry
Control Serum (Human
Origin)
Level 3

WARNING

H302
EUH208
P264

P301+P312

Harmful if swallowed.
May produce an allergic reaction.
Wash hands thoroughly after handling.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
γ-Glutamyltransferase < 0.1%
Ethylene Glycol 30 - 40%

Safety Data Sheet is available at techdocs.beckmancoulter.com

STORAGE AND STABILITY

Analytes in unopened bottles will remain stable until the expiration date on the label when the control serum is stored between -15°C and -20°C, or, once opened, up to 20 days when stored between +2°C and +8°C. Storage below -20°C may result in product damage.

CONTENTS

PRODUCT	QUANTITY	REF
SYNCHRON Control Multi-Level	6 x 20 mL Bottles (2 Bottles each of Level 1, Level 2, and Level 3) 1 Value Assignment Sheet	657365

INSTRUCTIONS FOR USE

Store unopened bottles between -15°C and -20°C. Materials, once opened, to be used within twenty days may be stored between +2°C and +8°C for laboratory convenience. Storage below -20°C may result in product damage. As for any solution containing bilirubin, minimize exposure of the material to light, heat and air.

Prior to initial use, remove the plastic screw cap and open the sealed plug by cutting the tip. Mix the contents with a gentle swirling motion before dispensing. To prevent leakage into the cap, store upright and do not mix by inversion. To maximize stability, remove the bottle from the specified storage temperature and immediately dispense the control as required. Then immediately return the bottle to the specified storage temperature. The plug and screw cap should remain on the bottle while in storage.

In order to assess system performance, it is essential that procedures in the SYNCHRON Operation Manual be followed precisely while using these control materials.

ASSIGNED VALUES

Multiple instruments and laboratories using representative samples from this lot of control established the MEAN assigned values. Values are specific to the assay methodologies of the SYNCHRON DxC reagents and Calibrator

system. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

Performance range was obtained by combining estimates of variance as determined from participating laboratory data and other currently available studies. Average value obtained in the laboratory should fall within the performance range, although recovery may not be identical to the mean value listed. Those outside the range may indicate unsatisfactory performance. The causes for such a discrepancy may be variation in techniques, equipment, calibration etc., which should be evaluated by the laboratory.

Since assay values are published in these instructions for all three levels of this product, care should be taken in matching lot numbers.

LIMITATIONS

Please note that light, heat and air will decompose bilirubin in any solution. Bilirubin may degrade faster than other analytes. Therefore, the laboratorian is encouraged to review the recovery of bilirubin throughout the life of the product.

These control products are not compatible with direct reading (undiluted sample) ion selective electrodes and should not be used on these systems.

Since the viscosity of these control products is higher than that of normal serum, care should be taken when pipetting small volumes of sample. It is recommended that the pipet tip be held only slightly below the meniscus while aspirating to avoid adherence of excess amounts of control on the tip of the pipet.

SYNCHRON CONTROL is fortified with porcine pancreatic amylase. SYNCHRON pancreatic amylase (P-AMY) and SYNCHRON total amylase (AMY) reagents utilize different substrates. Amylase isoenzymes from different species have different affinities for specific substrates, therefore pancreatic amylase values in this control will recover higher than total amylase values with these reagents. SYNCHRON CONTROL does contain endogenous levels of human salivary amylase.

ADDITIONAL INFORMATION

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

Each analyte in this control is traceable to the calibrator as indicated in Table 1. Traceability process is based on prEN ISO 17511.

Revision History




Revision AF

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AG

Added new language requirement: Bulgarian, Serbian, and Vietnamese. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	CE Mark		Consult Instructions for Use
	Authorized Representative in the European Community		Date of Manufacture
	Warning		Biological risks
	Control		Value Assignment Sheet
	Level		
			Made in USA of US and Foreign Components

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REFERENCES

1. U.S. Patent No. 3,876,375.
2. Frajola, W.J., and J.A. Maurukas, A Stable Human Reference Serum, *Health Laboratory Sciences* 13:25 (1976).
3. Grannis, G.F., Interlaboratory Survey of Enzyme Analyses, *Amer. J. Clin. Path.* 66:206 (1976).
4. Tanishimana, K., Y. Minamikawa, N. Yokogawa, and M. Takeshita, Protective Effect of Glycerol Against the Increasing Alkaline Phosphatase Activity of Lyophilized Quality Control Serum, *Clin. Chem.* 23:1873 (1977).
5. Williams, G.Z., D.S. Young, M.R. Stein, and E. Cotlove, Biological and Analytic Components of Variation in Long-Term Studies of Serum Constituents in Normal Subject, *Clin.Chem.* 16:1016 (1970).
6. Helman, E.Z., I.M. Reingold, I.O. Gleason, Plea for Standardization of Commercial Calibration Materials for Automated Instruments, *Clin. Chem.* 17:1144 (1971).
7. Grannis, G.F., An Optimum System of Quality Control, *Clin. Chem.* 19:649 (1973).
8. Fraser, C.G., A.N. Fudge, and L. Penberthy, Evaluation of Precision Using Lyophilized Quality Control Materials. *Ann. of Clinical Biochemistry* 15:121 (1978).
9. Peacher, W. A., The Evaluation of a Stabilized Liquid Chemistry Control, *Clin Chem.* 25:1119 (1979).

Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale 1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.