Vigil Protein Control Levels 1, 2, and 3

REF 450120, 450130, 450125

TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	6
Deutsche (de-de)	10
Italiano (it-it)	14
Español (es-es)	18
Português (Portugal) (pt-pt)	22
Dansk (da-dk)	26
Svenska (sv-se)	30
Norsk (nb-no)	34
Suomi (fi-fi)	38
Ελληνικά (el-gr)	42
日本語 (ja-jp)	46
中文 (中国) (zh-cn)	49
Lietuvių (lt-lt)	52
Magyar (hu-hu)	56
Polski (pl-pl)	60
Čeština (cs-cz)	64
Slovák (sk-sk)	68
한국어 (ko-kr)	72
Türkçe (tr-tr)	75
Русский (ru-ru)	79
Eestlane (et-ee)	83
Hrvatski (hr-hr)	87
Български (bg-bg)	91
Română (ro-ro)	95
Slovenščina (sl-si)	99
Srpski (sr-sp)	103
Latviešu (lv-lv)	107
Українська (uk-ua)	111
Português (Brasil) (pt-br)	115
Nederlands (nl-nl)	119
Tiếng Việt (vi-vn)	123
Македонски (mk-mk)	127
Қазақша (kk-kz)	131
REFERENCES	135

Vigil Protein Control Levels 1, 2, and 3

REF 450120, 450130, 450125

For *In Vitro* Diagnostic Use

For Professional Use Only

Rx Only

INTENDED USE

Vigil Protein stabilized liquid controls are designed for monitoring the reliability and overall performance of specific protein test systems in the clinical laboratory. The use of three levels of control enables the laboratorian to monitor changes in calibration along with analytical error and imprecision.

SUMMARY

The Vigil Protein Controls are derived from fresh frozen human plasma that has been defibrinated and then stabilized with ethylene glycol.^{1,2,3,4} The stabilizing effects of ethylene glycol are threefold. The high osmolality minimizes bacterial growth while its antioxidant property stabilizes oxygen-labile constituents. Additionally, the presence of ethylene glycol leads to freezing point depression, which allows the controls to remain in the liquid state at normal freezer temperature; i.e., between -15°C and -20°C. The preparation of these stabilized controls in liquid form eliminates errors commonly associated with filling, drying, and the reconstitution of lyophilized control products.^{3,5,6,7,8}

CAUTION

Because this product is of human source origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HB_sAg. Because no test method can offer complete assurance that HIV, Hepatitis B virus, and Hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

GHS HAZARD CLASSIFICATION

Vigil Protein Control Level 1 WARNING

H302
P264

P301+P312

Harmful if swallowed.
Wash hands thoroughly after handling.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
Ethylene Glycol 30 - 40%

Vigil Protein Control Level 2 WARNING

H302
P264

P301+P312

Harmful if swallowed.
Wash hands thoroughly after handling.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
Ethylene Glycol 30 - 40%

Vigil Protein Control Level 3 WARNING

H302
P264

P301+P312

Harmful if swallowed.
Wash hands thoroughly after handling.
IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
Ethylene Glycol 30 - 40%

 Safety Data Sheet is available at beckmancoulter.com/techdocs

STORAGE AND STABILITY

Analytes in unopened bottles will remain stable until the expiration date on the label when the control is stored between -15°C and -20°C, or, once opened, up to 60 days when stored between +2°C and +8°C. Storage below -20°C may result in product damage. Due to variability of protein stability, individual laboratories may notice a decrease in ceruloplasmin (CER) values over the life of the product. Therefore, the laboratorian is encouraged to review recovery of CER over time.

CONTENTS

PRODUCT	QUANTITY	REF
Vigil Protein Control, Level 1	4 x 5 mL	450120
Vigil Protein Control, Level 2	4 x 5 mL	450125
Vigil Protein Control, Level 3	4 x 5 mL	450130

REACTIVE INGREDIENTS

Vigil Protein Control, Level 1, 5.0 mL: Human Based Material

Vigil Protein Control, Level 2, 5.0 mL: Human Based Material

Vigil Protein Control, Level 3, 5.0 mL: Human Based Material

INSTRUCTIONS FOR USE

Store unopened bottles between -15°C and -20°C. Once opened, controls to be used within 60 days may be stored between +2°C and +8°C for laboratory convenience. Storage below -20°C may result in product damage. Mix the contents with a gentle swirling motion before dispensing. To prevent leakage into the cap, store upright and do not mix by inversion. To maximize stability, remove the bottle from the specified storage temperature and immediately dispense the control as required. Then immediately return the bottle to the specified storage temperature. The plug and screw cap should remain on the bottle while in storage. In order to assess system performance, it is essential that procedures in the system operation manual be followed precisely while using these control materials.

LIMITATIONS

This control contains ethylene glycol as a preservative and cannot be used for the determination of Total Protein by the index of refraction method. This is due to the alteration of the colligative properties of the control by ethylene glycol. Since the viscosity of these control materials is higher than that of normal serum, care should be taken when pipetting small volumes of control. It is recommended that the pipet tip be held only slightly below the meniscus while aspirating control to avoid adherence of excess amounts of material on the tip of the pipet. For automated systems, a slower speed setting is recommended.

ADDITIONAL INFORMATION

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a

serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The Summary of Safety and Performance is available from the EUDAMED database ec.europa.eu/tools/eudamed.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

Values are traceable and specific to the assay methodologies of the AU/DxC AU, IMAGE, and Synchron reagents and Calibrator system. Multiple instruments and laboratories using representative samples from this lot of control established the MEAN assigned values. Performance range was obtained by combining estimates of variance as determined from participating laboratory data and other currently available studies.

Average value obtained in the laboratory should fall within the performance range, although recovery may not be identical to the mean value listed. Those outside the range may indicate unsatisfactory performance. The causes for such a discrepancy may be variation in techniques, equipment, calibration etc., which should be evaluated by the laboratory. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias. Since assay values are published for all three levels of this product, care should be taken in matching lot numbers.

REVISION HISTORY

Revision AE

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AF

Added new language requirement: Brazilian Portuguese.

Revision AG

Added new language requirement: Dutch. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AH

Added new language requirement: Slovak

Updates to family of instruments.

Revision AJ

Updates to comply with IVDR requirements per Beckman Coulter Global Labeling Policy.

Revision AK

Added Summary of Safety and Performance statement.

Added languages: Latvian, Croatian, Estonian, Slovenian, Bulgarian, Finnish, Kazakh, Macedonian, Russian, Serbian, Ukrainian.



Revision AL

Added Vietnamese.

Revision AM

Updated Vietnamese translation.

Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	CE Mark		Consult Instructions for Use
	Authorized Representative in the European Community		Date of Manufacture
	Warning		Biological risks
	Control		Value Assignment Sheet
			Made in USA of US and Foreign Components

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REFERENCES

1. U.S. Patent No. 3,876,375.
2. Frajola, H. J. and J. A. Maurukas, *A Stable Liquid Human Reference Serum*, *Health Laboratory Sciences* 13:25 (1976).
3. Grannis, G.F., Interlaboratory Survey of Enzyme Analyses, *Amer. J. Clin. Path.* 66:206 (1976).
4. Tanishimana, K., Y. Minamikawa, N. Yokogawa, and M. Takeshita, Protective Effect of Glycerol Against the Increasing Alkaline Phosphatase Activity of Lyophilized Quality Control Serum, *Clin. Chem.* 23:1873 (1977).
5. Williams, G.Z., D.S. Young, M.R. Stein, and E. Cotlove, Biological and Analytic Components of Variation in Long-Term Studies of Serum Constituents in Normal Subject, *Clin. Chem.* 16:1016 (1970).
6. Helman, E.Z., I.M. Reingold, I.O. Gleason, Plea for Standardization of Commercial Calibration Materials for Automated Instruments, *Clin. Chem.* 17:1144 (1971).
7. Grannis, G.F., An Optimum System of Quality Control, *Clin. Chem.* 19:649 (1973).
8. Fraser, C.G., A.N. Fudge, and L. Penberthy, Evaluation of Precision Using Lyophilized Quality Control Materials. *Ann. of Clinical Biochemistry* 15:121 (1978).

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.