BECKMAN COULTER IMMUNOCHEMISTRY SYSTEMS

Beckman Coulter CAL 1 (Calibrator 1)

REF 449560

TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	5
Deutsche (de-de)	9
Italiano (it-it)	13
Español (es-es)	17
Português (Portugal) (pt-pt)	21
Dansk (da-dk)	25
Svenska (sv-se)	28
Norsk (nb-no)	31
Suomi (fi-fi)	34
Ελληνικά (el-gr)	38
日本語 (ja-jp)	42
中文 (中国) (zh-cn)	45
Lietuvių (lt-lt)	48
Magyar (hu-hu)	52
Polski (pl-pl)	56
Čeština (cs-cz)	60
Slovák (sk-sk)	64
한국어 (ko-kr)	68
Türkçe (tr-tr)	71
Eestlane (et-ee)	74
Hrvatski (hr-hr)	78
Български (bg-bg)	82
Română (ro-ro)	86
Slovenščina (sl-si)	90
Latviešu (lv-lv)	94
Українська (uk-ua)	97
Português (Brasil) (pt-br)	101
Nederlands (nl-nl)	105
REFERENCES	109

BECKMAN COULTER IMMUNOCHEMISTRY SYSTEMS

Beckman Coulter CAL 1 (Calibrator 1)

REF 449560

For *In Vitro* Diagnostic Use

For Professional Use Only

Rx Only

INTENDED USE

CAL 1 (Calibrator 1), when used in conjunction with Beckman Coulter alpha1-acid glycoprotein (AAG), complement C3 (C3), complement C4 (C4), haptoglobin (HPT), immunoglobulin A (IgA), immunoglobulin G (IgG), immunoglobulin M (IgM), kappa light chain (KAP), lambda light chain (LAM), and transferrin (TRF) reagents, is intended for use on IMMAGE Systems for the calibration of these reagents.

SUMMARY

The calibrator is derived from fresh frozen human plasma that has been defibrinated and processed. Assay of CAL 1 provides a response value that is utilized for the adjustment of pre-programmed calibration curves from which alpha1-acid glycoprotein, complement c3, complement c4, haptoglobin, immunoglobulin A, immunoglobulin G, immunoglobulin M, kappa light chain, lambda light chain, and transferrin concentration values are determined in test specimens.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

STORAGE AND STABILITY

Once opened, the calibrator is stable until expiration date printed on the calibrator bottle if stored capped in the original container at +2°C to +8°C.

NOTICE

Adverse storage conditions of this material may cause erroneous test results.

CONTENTS

COMPONENT	IMMAGE Component:
Beckman Coulter CAL 1 (Calibrator 1)	4 x 3 mL vials
Calibrator 1 Bar Code Card	1
Calibrator 1 Bar Code Strips	2
Value Assignment Sheet	1

NOTICE

Lot number on the Calibrator 1 vials and cards must be the same.

REACTIVE INGREDIENTS

CALIBRATOR 1, 3.0 mL: Human Based Material

⚠ CAUTION

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HB_sAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled as specified in the Centers for Disease Control's Biosafety Level 2 Guidelines.

⚠ CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

INSTRUCTIONS FOR USE

1. Allow contents of the bottle to come to room temperature.
2. Mix the contents by gently inverting the bottle.
3. Refer to the appropriate Operations Manuals and reagent Instructions for Use for the calibration procedures.
4. Locate the calibrator bar code label provided with the calibrator.
5. Place the appropriate label on an empty 13 x 100 mm or 16 x 100 mm test tube. The labeled test tube should be saved for re-use.
6. Place the appropriate calibrator in a sample cup.
7. Place the appropriate cup in the labeled test tube.
8. Place the tube into a sample rack.
9. Repeat steps 1-5 for any additional calibrators. Continue placing calibrators into the same sample rack.
10. Open the sample compartment lid.
11. Place the rack(s) containing calibrators in positions A and B.
12. Close the sample compartment lid.

ASSIGNED VALUES

The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the IMMAGE reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by intermethod bias.

IMMAGE values can be printed from the console.

LIMITATIONS

1. CAL 1 values are specific for Beckman Coulter Immunochemistry Analyzers only when using Beckman Coulter Reagents. Use of CAL 1 with other reagents or test procedures is not recommended.
2. Dust particles or other particulate matter in the reaction solution may result in extraneous light-scattering signals for the Immunochemistry system resulting in variable calibration rate units.

ADDITIONAL INFORMATION

Complete instructions for the use of this calibrator are contained in the Beckman Coulter Immunochemistry Systems Operating Manuals. Copies of these manuals are available from Beckman Coulter, Inc., Brea, California 92821.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The Summary of Safety and Performance is available from the EUDAMED database ec.europa.eu/tools/eudamed.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

Each measurand (analyte) in this calibrator is traceable to the IFCC reference preparation for plasma proteins, ERM-DA470/IFCC*.^{1,2,3} The traceability process is based on ISO 17511.

REVISION HISTORY

Revision AD

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AE

Added new language requirement: Brazilian Portuguese.

Revision AF

Added new language requirement: Dutch, and Slovak. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AG

Updates to comply with IVDR requirements per Beckman Coulter Global Labeling Policy.

Revision AH





Added Summary of Safety and Performance statement.

Added new languages: Latvian, Croatian, Estonian, Slovenian, Bulgarian, Finnish.

Revision AJ

Added Ukrainian Language

Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	CE Mark		Consult Instructions for Use
	Authorized Representative in the European Community		Date of Manufacture
	Caution		Biological risks
	Calibrator		Bar Code Card
	Value Assignment Sheet		Bar Code Strip
		Made in USA of US and Foreign Components	

*ERM is a registered trademark of EC-JRC-IRMM

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REFERENCES

1. Ingrid Zegers, Thomas Keller et al. Characterization of the New Serum Protein Reference Material ERM-DA470k/IFCC: Value Assignment by Immunoassay *Clinical Chemistry* 56:12 1880-1888(2010)
2. JRC-IRMM Certification Report "Certification of proteins in the human serum Certified Reference Material ERM - DA470k/IFCC
3. Johnson AM. A new international reference preparation for proteins in human serum (RPPHS). *Arch Pathol Lab Med* 1993; 117:29-31.

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.