DxH 500 Series Control

REF B36872

TABLE OF CONTENTS

English	2
Français (FR)	5
Deutsche (DE)	8
Italiano (IT)	11
Español (ES)	14
Português Portugal (PT-PT)	17
Dansk (DK)	20
Svenska (SE)	23
Norsk (NO)	26
Ελληνικά (GR)	29
日本語 (JP)	32
中文 (中国) (ZH-CN)	35
Lietuvių (LT)	38
Magyar (HU)	41
Polski (PL)	44
Čeština (CZ)	47
Slovák (SK)	50
한국어 (KR)	53
Türkçe (TR)	56
Русский (RU)	59
Hrvatski (HR)	62
Български (BG)	65
中文 (台灣) (ZH-TW)	68
Română (RO)	71
Srpski (Latinski) (SP)	74
Українська (UA)	77
Português Brasil (PT-BR)	80
Nederlands (NL)	83
Tiếng Việt (VN)	86

DxH 500 Series Control

REF B36872

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE

DxH 500 Series Control is an assayed whole blood control designed to monitor values on multi parameter hematology cell counters. Refer to the Table of Expected Results for specific instrument models.

SUMMARY AND PRINCIPLE

It is an established laboratory practice to use a stable control to monitor the performance of diagnostic tests. This control is composed of stable materials that provide a means of monitoring the performance of hematology blood cell counters. It is sampled in the same manner as a patient specimen.

REAGENTS

DxH 500 Series Control is composed of human erythrocytes, mammalian leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL. Each human donor/unit used in the preparation of this product has been tested by a FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures. DxH 500 Series Control is intended for use only by trained personnel.

GHS HAZARD CLASSIFICATION			
DxH 500 Series Control	WARNING	Causes skin irritation.	
			Causes serious eye irritation.
			Wear protective gloves, protective clothing and eye/face protection.
			IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
			If skin irritation occurs: Get medical advice/attention.
			If eye irritation persists: Get medical advice/attention.
			Potassium Hydroxide <1%

STORAGE, STABILITY, AND DISPOSAL

Store DxH 500 Series Control upright at 2-8°C when not in use. **Protect tubes from overheating and freezing.** Unopened tubes are stable until the expiration date, shown on the TABLE OF EXPECTED RESULTS. For open vial stability, refer to the TABLE OF EXPECTED RESULTS for your system. Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

After mixing, product should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. **Do not use the product if deterioration is suspected.**

INSTRUCTIONS FOR USE

1. Remove tubes from the refrigerator and verify that the lot number on the tubes matches the lot numbers on the Table of Expected Results.
2. Allow the tubes to warm to room temperature (15-30°C or 59-86°F) for 15 minutes before mixing.
3. To mix, hold a tube horizontally between the palms of the hands. **Do not pre-mix on a mechanical mixer.**
 - A. Roll the tube back and forth for 20-30 seconds; occasionally invert the tube.
 - B. Continue to mix in this manner until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.

- C. Gently invert the tube 8-10 times immediately before sampling.
4. Set Up and analyze the controls according to the procedure in the Instructions for Use for your instrument.
5. After sampling:
 - A. Clean residual material from the cap and rim if the tube is opened for sampling. Replace the cap tightly.
 - B. Return tubes to refrigerator within 30 minutes of use.

ASSIGNED VALUES AND EXPECTED RESULTS

Reagent differences, maintenance, operating technique, and calibration may contribute to inter-laboratory variation. Refer to the Table of Expected Results for the corresponding lot number of the control tube.

PERFORMANCE LIMITS

The performance of this product is assured only if it is properly stored and used as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

PERFORMANCE CHARACTERISTICS

Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.

Assay values on a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument is in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.

For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the Assay Sheet, if the control is suitable for the method.

TROUBLESHOOTING PROCEDURE

The Instructions for Use of the DxH 500 Series instrument contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other control lots currently in use.
- Previous control lot numbers and the data for these previous lots.
- Data from current repeatability study (N=10) using a fresh whole blood specimen and performed according to the Instructions for Use.
- Data from the IQAP and data from your last instrument calibration where appropriate.

PRODUCT AVAILABILITY

DxH 500 Series Control

REF B36872 - 2 x 2.3 mL Abnormal Low, 2 x 2.3 mL Normal, 2 x 2.3 mL Abnormal High

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

Revision History

Revision AH, 11/2018

Changes were made to:

- PRODUCT AVAILABILITY - Removed Product Part Number C21335
- REVISION HISTORY - Removed Revision AE

Revision AJ, 12/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch language

Revision AK, 10/2021

- Updated STORAGE, STABILITY, AND DISPOSAL section
- Updated website address
- Updated Legal Manufacturer phone number
- Updated copyright year

Revision AL, 02/2023

- Updated STORAGE, STABILITY, AND DISPOSAL section
-

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Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland +(353) (0) 65 683 1100
Manufactured for Beckman Coulter, Inc.