UniCel DxC Synchron Systems

HDL Calibrator

REF B23634

TABLE OF CONTENTS

English	2
Français (FR)	5
Deutsch (DE)	8
Italiano (IT)	11
Español (ES)	14
Português Portugal (PT-PT)	17
Dansk (DA)	20
Svenska (SV)	23
Norsk (NO)	26
Ελληνικά (EL)	29
日本語 (JP)	32
中文 (ZH-CN)	35
Lietuviškai (LT)	38
Magyar (HU)	41
Polski (PL)	44
Čeština (CZ)	47
한국어 (KO)	50
Türkçe (TR)	53
Български (BG)	56
Română (RO)	59
Србија (SR)	62
Tiếng Việt (VI)	65
REFERENCES	68

UniCel DxC Synchron Systems HDL Calibrator

REF B23634

For *In Vitro* Diagnostic Use
Rx Only

INTENDED USE

The UniCel DxC SYNCHRON Systems HDL Calibrator is designed to provide suitable calibration levels for Beckman Coulter UniCel DxC 600/800 SYNCHRON Systems employing the quantitative UniCel DxC SYNCHRON Systems HDL Cholesterol reagent (HDL).

SUMMARY

The UniCel DxC SYNCHRON Systems HDL Calibrator is a lyophilized human serum containing lipoproteins from the various lipoprotein classes including high-density lipoproteins.

REACTIVE INGREDIENTS

Human serum	
HDL Cholesterol	Variable
Sodium Azide	0.30%

After contact with skin, wash immediately with plenty of water. Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment. Avoid release to the environment. Refer to Safety Data Sheets.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use.

Do not ingest. Harmful if swallowed.

CAUTION

BIOHAZARD: This material and/or container must be disposed of as hazardous waste. Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum or plasma donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV-Ag, HIV-1/2 and HCV and nonreactive for HbsAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples to be handled as specified in Centers for Disease Control's Biosafety Level 2 guidelines.¹

GHS HAZARD CLASSIFICATION

UniCel DxC SYNCHRON Systems HDL Calibrator	H412	Harmful to aquatic life with long lasting effects.
	P273	Avoid release to the environment. Sodium Azide < 0.5%

 Safety Data Sheet is available at techdocs.beckmancoulter.com

STORAGE AND STABILITY

If unopened, HDL Calibrator should be stored at +2°C to +8°C until the expiration date printed on the calibrator bottle. Reconstituted calibrators are stable for 14 days at +2°C to +8°C or for 30 days at ≤ -20°C unless the expiration date is exceeded. Frozen calibrator should be thawed only once. Visible signs of microbial growth, gross turbidity, or precipitate in the calibrator may indicate degradation and warrant discontinuation of use.

CALIBRATOR PREPARATION

Use deionized water as the level one calibrator. Remove the metal cap around the calibrator bottle. Gently tap the bottle on the table to remove powder at the top of the stopper. Open the HDL calibrator bottle carefully, avoiding

loss of lyophilizate. Add exactly 1.00 mL of deionized water to the bottle of calibrator. Replace the stopper and let stand for 5 minutes at room temperature.

Gently invert until the contents are dissolved avoiding the formation of foam. DO NOT SHAKE. For detailed calibration instructions, refer to the calibration section of the appropriate UniCel DxC Synchron Clinical System Instructions for Use (IFU).

CONTENTS

COMPONENT	QUANTITY
HDL Calibrator	3 x 1 mL Bottles
Value Assignment Sheet	1

CAUTION

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

Calibration Diskette P/N B32236 (lot specific)

LIMITATIONS

The calibrator can only be used in conjunction with the UniCel DxC 600/800 Systems HDL Cholesterol reagent. This calibrator has not been tested with any other chemistry system. The results obtained with this calibrator are dependent on several factors, including proper storage of the calibrator, proper techniques and use of Beckman Coulter UniCel DxC 600/800 Systems and their respective reagents. Dust particles or other particulate matter in the reaction solution may result in extraneous light-scattering signals, resulting in variable sample analysis.

ADDITIONAL INFORMATION

For more detailed information on the UniCel DxC SYNCHRON Systems, refer to the appropriate system manual.

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

SHIPPING DAMAGE

If damaged product is received, notify your Beckman Coulter Clinical Support Center.

TRACEABILITY

HDL Cholesterol measurand (analyte) in this calibrator is traceable to the US CDC HDL Cholesterol reference method. Any instrument or reagent modification may invalidate the assigned value.

Revision History

Revision AD

Updates to comply with requirements per Beckman Coulter Global Labeling Policy.



Revision AE

Added new language requirement: Bulgarian, Serbian, and Vietnamese. Additional changes to comply with requirements per Beckman Coulter Global Labeling Policy.

Revision AF

Moved Calibration Diskette information from CONTENTS to MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT and revised CALIBRATOR PREPARATION section to align with content updates.

Symbols Key

	Catalogue Number		In Vitro Diagnostic
	Contents		Temperature limit
	Manufacturer		Expiration Date
	Batch code		Safety Data Sheet
	Caution		Consult Instructions for Use
	Calibration Diskette		Made in Japan
	Authorized Representative in the European Community		Date of Manufacture
	CE Mark		Biological risks
	Value Assignment Sheet		Calibrator

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.
Manufactured for Beckman Coulter, Inc.