
For DxC 500 AU customers only

Liquid Urine Chemistry Calibrator

REF DR0090

LOT 2209052B

## INTENDED USE

BECKMAN COULTER LIQUID URINE CHEMISTRY CALIBRATOR is intended for use as a urine chemistry calibrator for methods run on the BECKMAN COULTER AU® series of Clinical Chemistry Systems (refer to "Use Limitations" section below).

## CONSTITUENTS

LIQUID URINE CHEMISTRY CALIBRATOR is a human urine based product assayed for the following 8 urine constituents: Urine Urea Nitrogen (UUN), Calcium (Arsenazo), Calcium (oCPC), Glucose, Inorganic Phosphorus, Magnesium, Total Protein, and Uric Acid.

## PRECAUTIONS

- 1 FOR IN VITRO DIAGNOSTIC USE
- 2 WARNING: POTENTIAL BIOHAZARDOUS MATERIAL

LIQUID URINE CHEMISTRY CALIBRATOR is manufactured from human urine. This calibrator should be handled at the biosafety Level 2 as recommended for any potentially infectious human serum, urine or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 1984

- 3 DO NOT INGEST-Toxicity has not been established.

## PREPARATION

LIQUID URINE CHEMISTRY CALIBRATOR is ready to use. No preparation is needed.

## STORAGE AND STABILITY

Store unopened vials of LIQUID URINE CHEMISTRY CALIBRATOR at 2 to 8°C. When stored as directed, the calibrator is stable until the expiration date printed on the label.

Once opened, vials of LIQUID URINE CHEMISTRY CALIBRATOR are stable for 30 days.

## USE LIMITATIONS

This calibrator has not been tested for use with any other Chemistry Systems other than those listed on the Assay Value section.

For best results avoid prolonged exposure of the samples to air. Withdraw sample and replace cap immediately. Maintain the vial at 2 to 8°C. Assay the sample without delay.

The results obtained using this calibrator are dependent upon several factors, including proper storage of the calibrator, proper technique, and use of AU Clinical Chemistry Analyzers and their respective reagents.

The DxC 700 AU is not available in all geographies.

## INDICATIONS OF DETERIORATION

Bacterial contamination produces an increase in turbidity, a characteristic odor, and/or a decrease in the glucose level. Discard the vial if any of these phenomena are observed.

## RECOMMENDED PROCEDURE

- 1 Gently invert the vial several times to ensure a homogeneous mixture.
- 2 Transfer sufficient volume of the calibrator to sample cups. Handle this calibrator with the same care used for patient samples.
- 3 Replace the cap immediately and store unused calibrator at 2 to 8°C.
- 4 Refer to the appropriate Instrument User's Guide for System Calibration Information.
- 5 Good Quality Control Practices should be observed to assure proper System performance.
- 6 Do not pour previously dispensed calibrator back into original bottle.
- 7 The 2 dimensional bar code is intended for use with the DxC 700 AU analyzer only.

## VALUE ASSIGNMENT

The assay values for each constituent have been determined in accordance with BECKMAN COULTER test protocol (adopted from CLSI) and the values apply only to this lot of calibrator. All values for each BECKMAN COULTER AU Clinical Chemistry Analyzer have been assigned using their respective reagents. Any instrument or reagent modifications may invalidate these assay values.

**PLEASE ENSURE THAT THE LOT NUMBER ON THE CALIBRATOR VIAL IS THE SAME AS THE ONE LISTED IN THE ASSAY VALUE SECTION OF THIS PACKAGE INSERT.**

The assay values for the constituents are traceable to the National Institute of Standards and Technology (NIST).

## TECHNICAL ASSISTANCE

For technical assistance contact the BECKMAN COULTER Technical Services Department at 1 (800) 223-0130.

## REFERENCES

1. National Institute of Standards and Technology, Office of Standard Reference Material, Room 8322, Chemistry Building, Gaithersburg, MD 20899.

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## REFERENCIAS

- 1 National Institute of Standards and Technology, Office of Standard Reference Material, Room 8322, Chemistry Building, Gaithersburg, MD 20899.

**ORDERING INFORMATION:**

To place an order for LIQUID URINE CHEMISTRY CALIBRATOR, contact the Customer Service Department at 1 (800) 223-0125.

Catalog No. : DR0090  
 Description: URINE CHEMISTRY CALIBRATOR  
 Packaging: 6 X 5 mL

**ASSAY VALUES****LOT: 2209052B****EXP. DATE: 11/30/24****Beckman Coulter AU® Series Analyzers**

Test Name	Constituents	Units			
CAZ1U	Calcium (Arsenazo)	mmol/L	2.41	mg/dL	9.64
CAO1U	Calcium (oCPC)	mmol/L	2.7	mg/dL	10.8
GLU1U	Glucose	mmol/L	6.78	mg/dL	122
PHO1U	Inorganic Phosphorus	mmol/L	16.8	mg/dL	52.1
MG-1U	Magnesium	mEq/L	2.64	mg/dL	3.17
BUN1U	Urea Nitrogen	mmol/L	212	mg/dL	593
UA-1U	Uric Acid	µmol/L	1690	mg/dL	28.4

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
Form No. 32037B.02  
 Issued 12/15/16

# LIQUID URINE CHEMISTRY CALIBRATOR

DR0090

6 X 5 mL

Manufactured for:  
 Beckman Coulter, Inc.  
 250 S. Kraemer Blvd.  
 Brea, CA 92821  
 Made in USA