



AU/DxC AU US

Instructions For Use

© 2025 Beckman Coulter, Inc. All rights reserved.

HbA1c

HbA1c (Hemoglobin A1c)

REF

B00389 2 x 37.5 mL HbA1c R1

2 x 7.5 mL HbA1c R2

2 x 34.5 mL Total Hemoglobin R1

5 x 2 mL LYO HbA1c Calibrator (Levels 2-6)

For *In Vitro* Diagnostic Use

For Rx use only

PRINCIPLE

INTENDED USE

The HbA1c (Hemoglobin A1c) reagent, when used in conjunction with Beckman Coulter Systems, HbA1c Calibrators, and SYNCHRON and AU Hemolyzing Reagent, is intended for the quantitative determination of hemoglobin A1c concentration in human whole blood.

The absolute HbA1c and Total Hemoglobin (THb) values generated as part of the HbA1c assay are intended for use in the calculation of the HbA1c/Total Hemoglobin ratio, and must not be used individually for diagnostic purposes.

SUMMARY AND EXPLANATION

Reference^{1,2,3,4,5,6,7,8}

Measurement of hemoglobin A1c is accepted as a method to measure long-term glucose control in patients with diabetes mellitus (a chronic disorder associated with disturbances in carbohydrate, fat and protein metabolism and characterized by hyperglycemia). Determination of HbA1c provides an important tool for monitoring the efficiency of dietary control and therapy during treatment of diabetes mellitus. Long term treatment of the disease emphasizes control of blood glucose levels in preventing the acute complications of ketosis and hyperglycemia. In addition, long term complications such as retinopathy, neuropathy and cardiovascular disease can be minimized if blood glucose levels are effectively controlled.

The process of conversion from hemoglobin A to hemoglobin A1c depends on the blood glucose concentration. Since the average life of a red blood cell is 120 days, measurement of hemoglobin A1c can reflect the mean daily blood glucose concentration over the preceding two to three months and provides a much better indication of glycemic control than blood or urinary glucose determinations.

METHODOLOGY

The HbA1c assay (B00389) involves the use of four reagents: Total Hemoglobin R1, HbA1c R1, HbA1c R2, and Hemolyzing Reagent (sold separately as Cat. No. 472137). In a pre-treatment step, whole blood is mixed with the Hemolyzing Reagent in a 1:100 dilution and the resultant hemolysate is used. Tetradecyltrimethylammonium bromide (TTAB) in the hemolyzing reagent eliminates interference from leukocytes.

The concentrations of both HbA1c and Total Hemoglobin are determined. The HbA1c/Total Hemoglobin ratio is expressed either as mmol/mol (IFCC) or %HbA1c (DCCT/NGSP).

Total Hemoglobin Reagent is used to measure total hemoglobin concentration by a colorimetric method. Change in absorbance is measured at 570/660 nm.

HbA1c reagent is used to measure hemoglobin A1c concentration by a turbidimetric immunoinhibition method. In the reaction, hemoglobin A1c antibodies combine with HbA1c from the sample to form soluble antigen-antibody complexes. Polyhapten from the reagent then bind with the excess antibodies and the resulting agglutinated complex is measured turbidimetrically. Change in absorbance is measured at 340/700 nm.

SPECIMEN

SPECIMEN STORAGE AND STABILITY

Reference^{7,9}

Samples (non-pretreated) are stable up to 8 hours when stored at 25°C, 7 days when stored at 2...8°C and up to 3 months when frozen at -20°C. Whole blood samples are stable for 18 months at -70°C. Frozen samples should be thawed only once.

Hemolyzed (pre-treated) samples are stable up to 4 hours when stored at 15...25°C, up to 24 hours when stored at 2...8°C, if stored in a sealed container.

Note: All hemolyzed samples should be mixed thoroughly immediately prior to assay.

Each laboratory should evaluate sample handling procedures to avoid variable results.

Specimen storage and stability information provides guidance to the laboratory. Based on specific needs, each laboratory may establish alternative storage and stability information according to good laboratory practice or from alternative reference documentation.

Additional handling conditions as designated by this laboratory:

SPECIMEN COLLECTION AND PREPARATION

K₂-EDTA, K₃-EDTA, Li-Heparin or Na-Heparin whole blood (freshly drawn blood treated with EDTA is the preferred specimen). Draw a sufficient amount of blood to yield the necessary tube volume. (It is important to follow the tube manufacturer's recommendations).

It is recommended to aliquot the required volume of Hemolyzing Reagent from the primary container, immediately returning the primary container to recommended storage conditions. Allow the aliquot of the Hemolyzing Reagent to reach room temperature (15...25°C) before use, to allow for reproducible accurate dispensing. If the aliquot has not been brought to room temperature, it is less efficient and may require a longer period of time (e.g., > 1min) to fully hemolyze the whole blood sample.

Pre-treat samples and controls by dilution of whole blood with Hemolyzing Reagent. 1 part sample or control with 100 parts Hemolyzing Reagent (for example 10 µL of sample or control plus 1000 µL of Hemolyzing Reagent Cat. No. 472137). To ensure the highest level of pipetting accuracy when pipetting very low volumes, always choose the pipette with the lowest nominal volume possible and the smallest tip. Select a pipette whereby the volume required lies in the centre of the pipette capacity range. The pipette should be calibrated and maintained correctly.

1. Thoroughly mix (on a roller if possible) the whole blood sample before the aliquot is taken for pre-treatment. Settling of the red cells before the aliquot is taken for pre-treatment may cause altered total hemoglobin (THb) and glycated hemoglobin (A1c) results. Avoid the formation of foam. This could lead to unpredictable sampling or localized dried blood spots around the inside of the sample tube.
2. Ensure Pipette is measuring accurately to ensure quantity of blood is dispensed accurately. Ensure outside of the pipette tip is free of excess blood by careful cleaning with a lint free wipe taking care not to come in contact with the pipette opening as this will decrease the amount of whole blood being dispensed.
3. Use pipette mixing (carefully aspirating and dispensing several times) to ensure all the whole blood is dispensed into the hemolyzing solution. Mix thoroughly (on a roller if possible), avoid foaming and assay the hemolysate after

hemolysis is complete (allow at least 1 minute for Hemolysis). This is indicated by a colour change from red to brown-green. Incomplete Hemolysis can manifest as significantly lower levels of Total Hemoglobin.

4. Do not pipette directly from the Hemolyzing Reagent bottle; use a disposable tube. This will avoid potential contamination.

Please note that only SYNCHRON and AU Systems Hemolyzing Reagent Cat. No. 472137 can be used with this method.

If these steps are **not** carried out as outlined, **imprecision issues may ensue.**

Additional instructions for patient sample preparation as designated by this laboratory:

Additional type conditions as designated by this laboratory:

REAGENTS

CONTENTS

HbA1c Reagent

Reagent storage location in this laboratory:

WARNING AND PRECAUTIONS

Reference¹⁰

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL.

The calibrator is manufactured from human material; each donor used in the preparation of this material was tested by an FDA approved method for the presence of the antibody to HIV-1/2 and HCV as well as for hepatitis B surface antigen and was not repeatedly reactive. Because no test method can offer complete assurance that HIV-1/2, HCV, hepatitis B virus or other infectious agents are absent from biological materials, this product should be handled at the Biosafety Level 2 as recommended for any infectious human serum or blood specimen in the Centers for Disease Control and Prevention/National Institutes of Health manual, *Biosafety in Microbiological and Biomedical Laboratories*.

Exercise the normal precautions required for handling all laboratory reagents.

Dispose of all waste material in accordance with local guidelines.

This product contains material of animal origin. The product should be considered as potentially capable of transmitting infectious diseases.

REACTIVE INGREDIENTS

HbA1c Calibrator
Hemolysate (human and sheep)
0.9% tetradecyltrimethylammonium bromide
Non-reactive chemicals necessary for optimum system performance

HbA1c R1 Antibody Reagent		HbA1c R2 Polyhaptan Reagent	
Anti-human HbA1c Antibody (sheep)	≥ 0.5mg/mL	HbA1c Polyhaptan	≥ 8 µg/mL
MES (2-morpholino-ethanesulphonic acid) Buffer	0.025 mol/L	MES (2-morpholino-ethanesulphonic acid) Buffer	0.025 mol/L
TRIS (tris(hydroxymethyl)aminomethane) Buffer (pH 6.2)	0.015 mol/L	TRIS (tris(hydroxymethyl)aminomethane) Buffer (pH 6.2)	0.015 mol/L
Non-reactive chemicals necessary for optimum system performance			

Total Hemoglobin

Total Hemoglobin R1	
Phosphate Buffer, pH 7.4	0.02 mol/L
Non-reactive chemicals necessary for optimum system performance	

The concentrations of the reactive components of the reagents shown on the kit label are the actual concentrations in the individual R1/R2 vials. The reagent composition which is shown in the Instructions For Use is the final concentration of these components in the reaction cuvette after addition of R1, Sample, and R2.

GHS HAZARD CLASSIFICATION

HbA1c R1

WARNING



H317	May cause an allergic skin reaction.
H402	Harmful to aquatic life.
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing and eye/face protection.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before use.
	2-Methyl-4-isothiazolin-3-one < 0.06%

HbA1c R2

WARNING



H317 May cause an allergic skin reaction.
H402 Harmful to aquatic life.
P273 Avoid release to the environment.
P280 Wear protective gloves, protective clothing and eye/face protection.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before use.
2-Methyl-4-isothiazolin-3-one < 0.06%

HbA1c THb R1

WARNING



H317 May cause an allergic skin reaction.
H402 Harmful to aquatic life.
P273 Avoid release to the environment.
P280 Wear protective gloves, protective clothing and eye/face protection.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before use.
2-Methyl-4-isothiazolin-3-one < 0.06%

HbA1c Calibrator (Levels 2 - 6) WARNING



H315 Causes skin irritation.
H317 May cause an allergic skin reaction.
H319 Causes serious eye irritation.
H402 Harmful to aquatic life.
P273 Avoid release to the environment.
P280 Wear protective gloves, protective clothing and eye/face protection.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P337+P313 If eye irritation persists: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before use.
2-Methyl-4-isothiazolin-3-one < 0.05%
Tetradecyltrimethylammonium bromide 1 - 2%

MATERIALS NEEDED BUT NOT SUPPLIED WITH REAGENT KIT

SYNCHRON and AU Systems Hemolyzing Reagent (for use in sample preparation); Cat. No. 472137.

Storage location of Hemolyzing Reagent, 472137, in this laboratory:

0.9% Saline solution.

Storage location of the saline solution in this laboratory:

At least two levels of control material.

0.1M NaOH (Analar Grade (AR)).

Storage location of 0.1M NaOH in this laboratory:**EQUIPMENT AND MATERIALS**

For AU400/AU480, AU640/640^e/680, AU2700 /AU5400/AU5800 and DxC 700 AU Beckman Coulter Analyzers.

Storage location of test tubes or sample cups in this laboratory:**REAGENT PREPARATION**

Total Hemoglobin R1 and HbA1c R1 and R2 are ready for use, and can be placed directly on board the analyzer.

Bring Hemolyzing Reagent to room temperature prior to use.

REAGENT STORAGE AND STABILITY

1. The unopened reagents are stable, up to the stated expiration date when stored at 2...8°C.
2. Opened bottles of reagent are stable for 30 days when stored in the refrigerated compartment of the analyzer.

DO NOT FREEZE. Once opened, the SYNCHRON and AU Systems Hemolyzing Reagent is stable until the expiration date printed on the bottle label when stored and capped at +2°C to +8°C.

INDICATIONS OF DETERIORATION

Visible signs of microbial growth, gross turbidity, precipitate, or change in color in the HbA1c or Total Hemoglobin reagents may indicate degradation and warrant discontinuation of use.

Additional storage requirements as designated by this laboratory:

CALIBRATION

CALIBRATOR PREPARATION

1. For each individual level of calibrator, carefully remove the cap and stopper, avoiding any loss of lyophilized material. Ensure each level of calibrator is reconstituted individually to remove the possibility of cap and stopper mix up.
2. Add 2.00 mL of sterile deionized water at 15 to 25°C (room temperature) to the lyophilized material using a calibrated pipette or gravimetrically using an analytical balance. Ensure the deionized water being used for reconstitution is at room temperature for accurate aspiration and dispensing into the calibrator bottle. Failure to accurately measure 2.00 mL for each level of calibrator will lead to variability in test results.

Note: Gravimetric preparation is recommended over volumetric preparation where possible.

3. Replace the cap and stopper.
4. Invert to ensure any material on the stopper is in solution. Dissolve the contents completely by gently mixing for 30 minutes. Vortex each bottle for 5 seconds at medium speed. Avoid foaming. Ensure that all lyophilized material is reconstituted.
5. Record the date the calibrator was reconstituted on the bottle label.

CALIBRATOR STORAGE AND STABILITY

Unopened, the calibrators should be stored at 2 ... 8°C until the expiration date printed on the calibrator bottle.

Reconstituted calibrators are stable for 8 hours stored at 15...25°C or 30 hours stored at 2...8°C until the expiration date is exceeded.

Calibrators that are aliquoted immediately after reconstitution and stored at -20°C are stable for 60 days, as long as the expiration date of the un-reconstituted calibrator is not exceeded. Frozen calibrators should be thawed only once. After thawing, vortex each bottle for 5 seconds at medium speed. Avoid foaming.

The calibrator should be aliquoted to avoid multiple freeze thaw cycles.

Mix the calibrator thoroughly before aliquotting.

Ensure Calibrators are placed at 2...8°C or -20°C when not in use to limit degradation of the calibrator. Adhere to calibrator storage and stability periods as outlined. If any of these conditions have been breached please ensure a new calibrator is prepared.

These calibrators should only be used in conjunction with the reagent system described here on Beckman Coulter AU Systems.

Calibration Stability

Recalibrate the assay every 14 days or when the following occur:

1. Change in reagent lot or significant shift in control values.
2. Major preventative maintenance was performed on the analyzer or a critical part was replaced.

Following calibration, the resulting curve should be visually reviewed on the Beckman Coulter AU Analyzer for acceptability using the software options to access the Calibration Monitor.

Quality control procedures should be undertaken immediately following calibration in accordance with good laboratory practice.

CALIBRATION INFORMATION

Calibrators are included in the kit.

CALIBRATOR ASSIGNED VALUES

Refer to table of assigned values included in the kit.

NOTE: Calibrators are lot-specific and should not be interchanged. Calibrators DO NOT require pre-treatment with hemolyzing reagent prior to assay.

For the level 1 calibrator, 0.9% saline should be used.

Total Hemoglobin: Two Point Calibration

Level 1 (0.9% Saline) and HbA1c Calibrator Level 3 are used for calibration of the Total Hemoglobin assay.

HbA1c: Multi-point Calibration.

HbA1c Calibrator Levels 2 to 6 are used for calibration of the HbA1c assay.

TRACEABILITY

Reference^{11,12,13,14,15,16,17}

The calibrator HbA1c values are traceable to the IFCC HbA1c reference method via IFCC HbA1c reference material. The relationship between results from the NGSP network (DCCT aligned) and the IFCC network has been evaluated and a Master Equation has been developed for interconversion of results from IFCC (mmol/mol) to NGSP (%) units.

MASTER EQUATION

NGSP = (0.0915 x IFCC (mmol/mol)) + 2.15

The definition of the relationship between the two networks links IFCC traceable results to clinically meaningful HbA1c results from the DCCT and the United Kingdom Prospective Diabetes Study (UKPDS). The Master Equation also provides these DCCT results with traceability to a higher order reference method.

Results

To report %HbA1c in NGSP units, this has to be defined as a CALCULATED TEST in the INTER TEST or CALCULATED TEST menu. Enter the formula $(A/B)*a + b$, where A=HbA1c, B=THb, a=91.5, b=2.15.

To report HbA1c in IFCC units (mmol/mol), this has to be defined as a CALCULATED TEST in the INTER TEST or CALCULATED TEST menu. Enter the formula $(A/B)*a$, where A=HbA1c, B=THb, a=1000.

Note: HbA1c and THb must be selected in the same units.

QUALITY CONTROL

At least 2 levels of control material should be analyzed. Each laboratory should establish its own control frequency however good laboratory practice suggests that controls be tested each day patient samples are tested and each time calibration is performed.

The results obtained by any individual laboratory may vary from the given mean value. It is therefore recommended that each laboratory generates analyte specific control target values and intervals based on multiple runs according to their requirements. These target values should fall within the corresponding acceptable ranges given in the relevant product literature.

If any trends or sudden shifts in values are detected, review all operating parameters.

Each laboratory should establish guidelines for corrective action to be taken if controls do not recover within the specified limits.

Location of controls used at this laboratory.

--

CONTROL NAME	SAMPLE TYPE	STORAGE

TESTING PROCEDURE(S)

Total Hemoglobin and HbA1c tests must be performed on each pre-treated sample and control. Refer to the appropriate Beckman Coulter AU analyzer User Guide/Instructions For Use (IFU) and the accompanying AU Instrument Setting Sheets for analyser-specific assay instructions.

REPORTING RESULTS

REFERENCE INTERVALS

Reference^{18,19,20,21}

Adults: 4.0 – 6.0% (NGSP)
20 – 42 mmol/mol (IFCC units)

Reference Intervals shown above were taken from the literature. Expected values may vary with age, sex, sample type, diet and geographical location. Each laboratory should verify the transferability of the expected values to its own

population, and if necessary determine its own reference interval according to good laboratory practice. Results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Data contained within this section is representative of performance on Beckman Coulter systems. Data obtained in your laboratory may differ from these values.

DYNAMIC RANGE / ANALYTICAL MEASURING RANGE

Total Hemoglobin

The analytical range for Total Hemoglobin is 3.7-13.0 mmol/L (6-21 g/dL). When the result for Total Hemoglobin is outside the analytical range an "F" or "G" flag is generated and the calculated % HbA1c should not be reported. Settling of the red cells before the aliquot is taken for pre-treatment may cause an elevated Total Hemoglobin result. Samples exceeding the upper limit of Total Hemoglobin may be mixed well and the analyses repeated on a freshly hemolyzed sample.

HbA1c

The analytical range of this assay extends from 0.19 mmol/L (0.3 g/dL) to the concentration of Calibrator 6. If the HbA1c concentration is outside the analytical range an "F" or "G" flag is generated and the calculated % HbA1c should not be reported. Samples exceeding the upper limit of the analytical range for HbA1c should not be diluted, but instead should be reported as "% HbA1c > 15%" or "HbA1c>140 mmol/mol".

% HbA1c

The reportable range for the calculated HbA1c is 20 – 140 mmol/mol HbA1c (IFCC) and 4 – 15% HbA1c (NGSP), based on a typical Calibrator 6 value of 1.36 mmol/L (2.19g/dL), at a total hemoglobin of 9.6 mmol/L (15.5 g/dL). Note: a calculated result may be outside the reportable range based on either the Total Hemoglobin or the HbA1c result being outside of their respective analytical ranges.

Expected reference ranges in this laboratory:

INTERVALS	SAMPLE TYPE	UNITS

Additional reporting information as designated by this laboratory:

PROCEDURAL NOTES

ANTICOAGULANT TEST RESULTS

The following anticoagulants were assessed by Deming regression analysis with paired EDTA and heparin whole blood samples. Values of K2-EDTA (X), ranging from 4.7% HbA1c to 14.7% HbA1c, were compared with values for K3-EDTA, lithium heparin whole blood and sodium heparin whole blood (Y), yielding the following results:

ANTICOAGULANT	CONCENTRATION OF ANTICOAGULANT TESTED	DEMING REGRESSION ANALYSIS
K ₃ -EDTA	1.8 mg/mL	Y = 0.995X + 0.023 ; r = 0.999

Lithium Heparin	17 I.U./mL	$Y = 0.990X + 0.039$; $r = 0.999$
Sodium Heparin	17 I.U./mL	$Y = 0.998X + 0.007$; $r = 0.999$

LIMITATIONS

Reference^{18,22,23,24,25}

1. This assay is designed only for the measurement of mmol/mol HbA1c (IFCC) and %HbA1c (NGSP). Individual results for Hb and HbA1c concentration should not be reported.
2. Do not use this test for the diagnosis of diabetes mellitus. Performance characteristics for this use have not been determined.
3. This assay is not useful in judging day-to-day glucose control and should not be used to replace daily home testing of glucose.
4. Shortened red cell survival time will reduce the exposure of red cells to glucose, with a resultant decrease in HbA1c values. Causes of reduced red cell survival time include hemolytic anemia, or other hemolytic diseases, significant blood loss, blood transfusions, iron deficiency and pregnancy. Caution should be exercised when interpreting the HbA1c results from patients with these or other conditions affecting red cell survival time, and when the total hemoglobin is 5.6 mmol/L (<9 g/dL).
5. Sample carryover from the hemolysate can affect subsequent Urinary/CSF Protein (OSR6x70) results, and reagent carryover from HbA1c (B00389) can adversely affect Urine Amphetamines (OSR6323) and Urine THC (OSR6322) assays. Please refer to AU Contamination Avoidance Parameters for guidance.
6. Prior to weekly maintenance a mandatory additional W2 cycle using 0.1M NaOH is required when running HbA1c from standby, in batch mode or random access mode. When running the HbA1c assay in **batch mode from standby**, a mandatory W2 using 0.1M NaOH followed by a W2 using 1M HCl and photocal, is required after every 4th batch of samples, where the analyzer is returned to standby between batches. This represents, potentially, 4 occupations of the same cuvette. Batch sample capacity is based on the cuvette capacity of the specific analyzer in use.

When running the HbA1c assay in **random access mode** (i.e. together with other assays) additional parameters are necessary due to increased risk of cuvette coating by this assay. Please refer to AU Contamination Avoidance Parameters for guidance. If unacceptable drift or imprecision is observed in any Quality Control results or any calibration failures are observed, an additional W2 cuvette cleaning using 0.1M NaOH followed by 1M HCl and photocal, is recommended.

Note 1: Weekly maintenance should be carried out according to Beckman Coulter AU analyzer User Guide/Instructions For Use (IFU).

Note 2: As with any chemical reaction, users should be aware of the possible effect on results due to unknown interferences from medication or endogenous substances.

INTERFERENCES

Reference^{26,27,28}

Results of studies conducted to evaluate the susceptibility of the method to interference were as follows:

Icterus: Interference ≤ 6% up to 513 μmol/L (30 mg/dL) Bilirubin
Lipemia: Interference ≤ 7% up to 500 mg/dL Intralipid^{-a}

Ascorbic Acid Interference \leq 6% up to 50 mg/dL Ascorbic Acid
Rheumatoid factor Interference \leq 6% up to 1,000 IU/mL Rheumatoid Factor (RF):

Intralipid[®] is a registered trademark of Fresenius Kabi AB., Uppsala, Sweden.

Refer to references for further information on interfering substances.

SPECIFICITY

Reference^{26,9,29,30,23,24,31}

The HbA1c test shows no cross-reactivity with HbA0, HbA1a, HbA1b, acetylated hemoglobin, carbamylated hemoglobin and glycated albumin.

No significant effect of HbS, HbD, HbE, HbC, and HbF up to 10% was observed with this assay. Glycated HbF is not detected by the HbA1c assay as it does not contain the glycated β -chain. However, HbF is measured in the THb assay.

Samples containing $>10\%$ HbF may result in lower than expected HbA1c results.

No significant effect ($\leq 10\%$) of labile glycated hemoglobin (up to 2,000 mg/dL, 5 hours at $+37^{\circ}\text{C}$) was observed with this assay.

Laboratory specific procedure notes:

PERFORMANCE CHARACTERISTICS

SENSITIVITY

Reference³²

The limit of blank (LoB) and limit of detection (LoD) were determined in accordance with the CLSI EP17-A guideline. The LoB is calculated from $n \geq 60$ measurements of an analyte-free sample, and corresponds to the concentration below which analyte-free samples are found with 95% confidence. The Limit of Detection (LoD) corresponds to the sample concentration above the LoB which is detectable with 95% confidence.

Total Hemoglobin

LoB = 0.05 mmol/L (0.09 g/dL)

LoD = 0.10 mmol/L (0.16 g/dL)

HbA1c

LoB = 0.12 mmol/L (0.19 g/dL)

LoD = 0.13 mmol/L (0.22 g/dL)

METHODS COMPARISON

Reference³³

Patient samples were run in singlicate to compare this HbA1c (B00389) assay against OSR6192 on the AU680 analyzer. Results of Deming regression analysis were as follows:

Slope = 1.036	Intercept = -0.3821	r = 0.9968	n = 116	Sample range = 4.6 – 12.0% HbA1c
---------------	------------------------	------------	---------	----------------------------------

Patient samples were run in singlicate to compare this HbA1c (B00389) assay on the AU680 analyzer against HbA1c (650262) on the DXC800 analyzer. Results of Deming regression analysis were as follows:

Slope = 0.901	Intercept = 0.3140	r = 0.9941	n = 130	Sample range = 4.9 – 14.2% HbA1c
---------------	-----------------------	------------	---------	----------------------------------

PRECISION

Reference²²

Correctly operating AU Systems should exhibit precision values less than or equal to the following:

TYPE OF IMPRECISION	SAMPLE TYPE	%HbA1c (NGSP) % CV
Within-run	Whole Blood Hemolysate	4.0
Total	Whole Blood Hemolysate	4.0

Estimates of imprecision, based on CLSI recommendations are consistent with typical performance. Comparative performance data for the AU2700 system evaluated using the CLSI approved guideline EP5-A2 appears in the table below. Each laboratory should characterize their own instrument performance for comparison purposes.

The following data was obtained using 3 pools analyzed over 20 days.

TYPE OF IMPRECISION	SAMPLE TYPE	No. of Data Points	Mean %HbA1c	CV%	SD (%HbA1c)
Within-run	Hemolysate Control Pool 1	80	5.3	1.44	0.08
	Hemolysate Control Pool 2	80	7.4	1.03	0.08
	Hemolysate Control Pool 3	80	9.4	1.03	0.10
Total	Hemolysate Control Pool 1	80	5.3	2.07	0.11
	Hemolysate Control Pool 2	80	7.4	1.84	0.14
	Hemolysate Control Pool 3	80	9.4	1.68	0.16

ADDITIONAL INFORMATION

DxC 700 AU has a standard format for abbreviated Test Name for each reagent application. Refer to the table below for the abbreviated Test Name assigned to each application for this assay.

Test Name	Description
A1C1G	HbA1c (serum)
THB1G	Total Hemoglobin (Serum)

Setting Sheet Footnotes

§ For the Level 1 (zero) calibrator, saline should be used.

User defined

† HbA1c Calibrators included in the kit.

* Values set for working in mEq/L (mmol/L)

** Concentration of 0.19 mmol/L and Calibrator 6.

Note: HbA1c and THb must be selected in the same units.

For AU400, AU640/640^e and AU2700/AU5400 Beckman Coulter Analyzers.

¥ For determination of %HbA1c (NGSP units), the THb and HbA1c tests are used. A third test, %HbA1c must also be entered in the general tests and the calculated tests of the test section menu (no settings required). Set this test as CALCULATED TEST in the INTER TEST menu. Enter the formula $(A/B)*a+b$, where A=HbA1c, B=THb, a=91.5 and b=2.15. To report results as mmol/mol HbA1c (IFCC units), enter the formula $(A/B)*a$, where A=HbA1c, B=THb, a=1000.

For AU480, AU680, AU5800 and DxC 700 AU Beckman Coulter Analyzers.

¥ For determination of %HbA1c (NGSP units), the THb and HbA1c tests are used. A third test, %HbA1c must also be entered in the general tests and the calculated tests of the test section menu (no settings required). Set this test in common Test Parameter test name as CALCULATED TEST. In Specific Test Parameters Calculated Test enter the formula $(A/B)*a+b$, where A=HbA1c, B=THb, a=91.5 and b=2.15. To report results as mmol/mol HbA1c (IFCC units), enter the formula $(A/B)*a$, where A=HbA1c, B=THb, a=1000

Φ AU680

REVISION HISTORY

| Updated REAGENTS section


Preceding version revision history

Revised GHS section

REFERENCES

1. Sperling M, ed. Physician's Guide to Insulin-dependent (Type 1) Diabetes: Diagnosis and Treatment, American Diabetes Association, Inc., Alexandria, VA, 1998
2. Rochman H. Hemoglobin A1c and Diabetes Mellitus. *Ann Clin Lab Sci* 1980; 10(2): 111-115.
3. The Diabetes Control and Complications Trial Research Group. The Effect of Intensive Treatment of Diabetes on the development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. *New England Journal of Medicine* 1993; 329: 977-986.
4. Heinze, E, Kohne E, Meissner C, Beischer W, Teller WM, Kleihauer E. Hemoglobin A1c (HbA1c) in Children with Long Standing and Newly Diagnosed Diabetes Mellitus. *Acta Paediatr Scand.* 1979;68(4):609-12.
5. Mortensen HB. Glycated Hemoglobin. *Dan Med Bull* 1985; 32(6); 309-328.
6. Lehmann P. Homogeneous Immunospectrometric Assay for Hemoglobin A1c Adaptable for Most Clinical Chemistry Analyzers: A new Concept in the Care of Diabetic Patients. AACC 45th National Meeting, 1993.
7. Burtis CA, Ashwood ER, Bruns DE, eds. Carbohydrates In: *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*, 4th Edn. St. Louis, MO: Elsevier Saunders 2005.
8. Goldstein, D.E. et al. ADA; Tests of Glycemia in Diabetes. *DIABETES CARE*, VOLUME 27, SUPPLEMENT 1, JANUARY 2003: S106-S109.
9. Data on file at Beckman Coulter.
10. CDC-NIH manual, Biosafety in Microbiological and Biomedical Laboratories, 5th Edition (CDC 21-1112), U.S. Government Printing Office, Washington, D.C. (2009).
11. Jeppsson J-O, Kobold U, Barr J, Finke A, Hoebel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umemoto M, Weykamp C. Approved IFCC reference method for the Measurement of HbA1c in Human Blood. *Clin Chem Lab Med* 2002; 40(1): 78-89.
12. Hoebel W, Weykamp C, Jeppsson J-O, Miedema, Barr JR, Goodall I, Hoshino T, John, WG, Kobold U, Little R, Mosca A, Mauri P, Paroni R, Susanto F, Takei I, Theinpoint L, Umemoto M, Wiedmeyer H-M. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardisation Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. *Clin Chem* 2004; 50(1): 166-174.
13. Weykamp, C. et al. The IFCC Reference Measurement System for HbA1c: A 6-Year Progress Report *Clinical Chemistry* 54:2: 240-248 (2008).
14. Andrea Geistanger, Sabine Arends, Christoph Berding, Tadao Hoshino, Jan-Olof Jeppsson, Randie Little, Carla Siebelder and Cas Weykamp on behalf of the IFCC Working Group on Standardization of HbA1c: Statistical Methods for Monitoring the Relationship between the IFCC Reference Measurement Procedure for Hemoglobin A1c and the Designated Comparison Methods in the United States, Japan and Sweden. *Clin Chem* 2008, 54 (8): 1379-85.
15. Andrea Mosca, Ian Goodall, Tadao Hoshino, Jan O. Jeppsson, W. Garry John, Randie R. Little, Kor Miedema, Gary L. Myers, Hans Reinauer, David B. Sacks and Cas W. Weykamp. Global standardization of glycated hemoglobin measurement: the position of the IFCC Working Group. *Clin Chem Lab Med* 2007, 45(8): 1077-1080
16. Little RR, Rohlfing CL, Sacks DB. Status of HbA1c measurement and goals for improvement: From chaos to order for improving diabetes care. *Clin Chem* 2011; 57:205-214.

17. National Glycohemoglobin Standardization Program. NGSP protocol. <http://www.ngsp.org/docs/Protocol.pdf> (Accessed December 2011).
18. Panteghini M, John WG. Implementation of Hemoglobin A1c results traceable to the IFCC reference system: the way forward. *Clin Chem Lab Med* 2007; 45(8): 942-944.
19. Clinical and Laboratory Standards Institute (CLSI), Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory, Approved Guideline, 3rd Edition, CLSI document C28-A3 (ISBN 1-56238-682-4), Wayne, PA (2008).
20. Wu, A., ed., Tietz Clinical Guide to Laboratory Tests, 4th Edition (ISBN 0-7216-0189-8), Saunders Elsevier, St. Louis, MO (2006).
21. McPherson, R.A., Pincus, M.R., Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd Edition (ISBN 978-1-4377-0974-2), Saunders Elsevier, Philadelphia, PA (2011).
22. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Quantitative Measurement Methods Approved Guideline - Second Edition, CLSI Document EP5-A2, Wayne PA 2004.
23. Greiling H, Gressner AM, Lehrbuch der Klinischen Chemie und Pathobiochemie, 2nd ed., Stuttgart New York: Schattaner-Verlag: 1989: 208.
24. Bessis M. Living Blood Cells and Their Ultrastructure, Berlin, Springer-Verlag 1973.
25. Sacks, D.B., A1C Versus Glucose Testing: A Comparison, *Diabetes Care*, 34(2), pp 518-523 (2011).
26. Young DS. Effects of Drugs on Clinical Laboratory Tests, AACC, 5th ed. AACC Press, 2000.
27. Young, D. S. and Friedman, R. B., Effects of Disease on Clinical Laboratory Tests, 4th Edition (ISBN 1-890883-45-X), AACC Press, Washington, D.C. (2001).
28. Young, D. S. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd Edition (ISBN 978-1-59425-068-2), AACC Press, Washington, D.C. (2007).
29. Chang J, Hoke C, Ettinger B, Penerian G. Evaluation and Interference Study of Hemoglobin A1c Measured by Turbidimetric Inhibition Immunoassay. *Am J Clin Pathol* 1998; 109: 274-278.
30. Rohlfing C, Connolly S, England J, Little R. Effect of Elevated Fetal Hemoglobin on HbA1c measurements: four common assay methods compared to the IFCC reference method. *Clin Chem* 2006; 52 Suppl 6: A108.
31. Miedema K., "Influence of Hemoglobin Variants on the Determination of Glycated Hemoglobin", *Klinisches Labor*, No. 39, pp 1029-1032 (1993).
32. Clinical and Laboratory Standards Institute (CLSI), Protocols for Determination of Limits of Detection and Limits of Quantitation, Approved Guideline, CLSI document EP17-A (ISBN 1-56238-551-8), Wayne, PA (2004).
33. Clinical and Laboratory Standards Institute (CLSI), Method Comparison and Bias Estimation Using Patient Samples, Approved Guideline, Second Edition (Interim Revision), CLSI Document EP09-A2-IR (ISBN 1-56238-731-6), Wayne, PA (2010).

 Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.
+(1) 800-854-3633
www.beckmancoulter.com