

# ALBUMIN, AU480/AU680

System Reagent: OSR6x02

Reagent ID: 002

Specific Test Parameters	
General	Range
Test Name: ALB1U	Type: Serum
Sample Volume: 1.6 $\mu$ L	Dilution: 0 $\mu$ L
Pre-Dilution Rate: 1	OD Limit: Min. OD Max. OD
Reagents Volume: R1(R1-1) 46 $\mu$ L	Dilution: 194 $\mu$ L
R2 (R2-1) 0 $\mu$ L	Dilution: 0 $\mu$ L
Wavelength: Pri. 600 nm	Sec. 800 nm
Method: END	Dynamic Range Low: 1.5* High: 6*
Reaction slope: +	Correlation Factor A: 1 B: 0
Measuring Point 1: First 0 Last 1	Factor for Maker A: 1 B: 0
Measuring Point 2: First Last	Onboard Stability: 90 Day 0 Hour
Linearity: %	LIH Influence Check: #
No Lag Time:	Lipemia: +++++
	Icterus: +++++
	Hemolysis: +++++

Specific Test Parameters	
General	Range
Test Name: ALB1U	Type: Serum
Value/Flag: #	Level L: # Level H: #
Specific Ranges:	
From	To
Sex	Year
1. #	#
2. #	#
3. #	#
4. #	#
5. #	#
6. #	#
7. No demographics	
8. Not within expected values	
Unit: g/dL*	Decimal Places: #
Panic Value	
Low	High
#	#

Calibration Specific	
General	ISE
Test Name: ALB1U	Type: Serum
Calibration Type: AB	Formula: Y=AX+B
Counts: #	Slope Check: None
<Calibrator Parameters>	
Calibrator $\uparrow$	OD
Point 1: #	Conc $\uparrow$
Point 2:	Factor Range
Point 3:	Low 4.5* High 8.5*
Point 4:	
Point 5:	
Point 6:	
Point 7:	
Point 8:	
Point 9:	
Point 10:	
<Point Cal. For Master Curve>	
Calibrator	No. of Correction Points
Point 1:	OD
Point 2:	Conc
MB Type Factor:	1-Point Calibration Point
Use Master Curve	Use Lot Calibration
Stability	Reagent Blanks
Calibration	30 Day 0 Hour
	30 Day 0 Hour
Operation	Interval (RB/ACAL)
Yes	Lot / Lot

# User defined  
 $\uparrow$  Beckman Coulter System Calibrator Cat. No. DR0070.  
 \* Values set for working in g/dL. To work in SI units (g/L) multiply by 10.

# ALBUMIN, AU5800

System Reagent: OSR6x02

Reagent ID: 002

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c		Calculated Test	Range
Test Name: <input type="text" value="ALB1U"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>						
Sample Volume	<input type="text" value="1.1"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	OD Limit		
Pre-Dilution Rate	<input type="text" value="1"/>	Diluent Bottle	<input type="text" value="#"/>	Min.OD	<input type="text" value=""/>	
Rgt. Volume	R1(R1-1) <input type="text" value="32"/> $\mu\text{L}$	Dilution	<input type="text" value="133"/> $\mu\text{L}$	Reagent OD Limit		
	R1-2 <input type="text" value=""/>	Dilution	<input type="text" value=""/>	1 <sup>st</sup> .	Low	<input type="text" value="-0.1"/> High <input type="text" value="0.5"/>
				Last	Low	<input type="text" value="-0.1"/> High <input type="text" value="0.5"/>
	R2(R2-1) <input type="text" value="0"/> $\mu\text{L}$	Dilution	<input type="text" value="0"/> $\mu\text{L}$	Dynamic Range Low <input type="text" value="1.5*"/> High <input type="text" value="6*"/>		
Common Rgt. Type	<input type="text" value="None"/>	Name	<input type="text" value="None"/>		Correlation Factor A	<input type="text" value="1"/>
Wavelength	Pri <input type="text" value="600"/> $\text{nm}$	Sec.	<input type="text" value="800"/> $\text{nm}$		Factor for Maker A	<input type="text" value="1"/>
Method	<input type="text" value="END"/>					
Reaction Slope	<input type="text" value="+"/>	Onboard Stability Period		<input type="text" value="90"/> Day	<input type="text" value="0"/> Hour	
Measuring Point1 1 <sup>st</sup>	<input type="text" value="0"/>	Last	<input type="text" value="1"/>		LIH Influence Check	<input type="text" value="#"/>
Measuring Point2 1 <sup>st</sup>	<input type="text" value=""/>	Last	<input type="text" value=""/>		Lipemia	<input type="text" value="++++"/>
Linearity Limit	<input type="text" value=""/>			Icterus	<input type="text" value="++++"/>	
Lag Time Check	<input type="text" value=""/>			Hemolysis	<input type="text" value="+++"/>	

Parameters		Specific Test Parameters				
General	LIH	ISE	HbA1c		Calculated Test	Range
Test Name: <input type="text" value="ALB1U"/> < > Type: <input type="text" value="Serum"/>						
Value/Flag: <input type="text" value="#"/>						
Level			Low	<input type="text" value="#"/>		High <input type="text" value="#"/>
Specific Ranges:						
	Sex	Year	Month	Year	Month	High
<input type="checkbox"/> 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>
7. Standard demographics						
8. Not within expected values						
Panic Value	Low	<input type="text" value="#"/>		High	<input type="text" value="#"/>	Unit <input type="text" value="g/dL*"/> Decimal Places <input type="text" value="#"/>

Parameters		Calibration Parameters				
Calibrators	Calibration Specific					
General	ISE					
Test Name: <input type="text" value="ALB1U"/> < > Type <input type="text" value="Serum"/> Cuvette . <input type="text" value=""/>						
<input type="checkbox"/> Use Serum Cal.						
Calibration Type: <input type="text" value="AB"/>		Formula: <input type="text" value="Y=AX+B"/>		Counts: <input type="text" value="#"/>		
<Calibrator Parameters>						
	Calibrator	OD	Conc	Range		
Point 1:	<input type="text" value="#"/>	<input type="text" value=""/>	<input type="text" value="†"/>	Low	High	
Point 2:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value="5.4*"/>	<input type="text" value="9.6*"/>	
Point 3:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 4:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 5:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 6:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 7:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
<Point Cal. For		No. of Correction Points	<input type="text" value=""/>	Use Master Curve	<input type="text" value=""/>	
Master Curve>		OD Range		<input type="checkbox"/> Lot Calibration		
	Calibrator	OD	Conc	Low	High	
Point-1	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	
MB Type Factor: <input type="text" value=""/>		1-Point Calibration Point <input type="text" value=""/>		<input type="checkbox"/> with Conc-0		
Slope Check		<input type="text" value=""/>				
Allowance Range Check						
<input type="checkbox"/> Reagent Blank		<input type="text" value=""/>				
<input type="checkbox"/> Calibration		<input type="text" value=""/>				
Advanced Calibration Operation		<input type="text" value="Yes"/>				
Interval (RB/ACAL)		<input type="text" value="Lot/Lot"/>				
Stability						
Reagent Blank	<input type="text" value="30"/> Day	<input type="text" value="0"/> Hour				
Calibration	<input type="text" value="30"/> Day	<input type="text" value="0"/> Hour				

# User defined  
 † Beckman Coulter System Calibrator Cat. No.: DR0070  
 \* Values set for working in g/dL. To work in SI units (g/L) multiply by 10

# Name: ALBUMIN, DxC 700 AU Serum Application

System Reagent: OSR6x02

Test name: ALB1U

Reagent ID: 002

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> ALB1U <input type="text"/> <b>Test No</b> <input type="text"/> <b>Type:</b> Serum <input type="text"/> <b>Operation</b> Yes <input type="text"/>				
Sample Volume <input type="text" value="1.1"/> $\mu\text{L}$	Dilution <input type="text" value="0"/> $\mu\text{L}$	OD Limit <input type="text"/>	Min. OD <input type="text"/>	Max OD <input type="text"/>
Pre-Dilution Rate <input type="text" value="1"/>		Reagent OD Limit <input type="text"/>	1 <sup>st</sup> Low <input type="text" value="-0.1000"/>	High <input type="text" value="0.5000"/>
Reagent Volume R1 (R1-1) <input type="text" value="32"/> $\mu\text{L}$	Dilution <input type="text" value="133"/> $\mu\text{L}$	Last Low <input type="text" value="-0.1000"/>	High <input type="text" value="0.5000"/>	
R1-2 <input type="text"/> $\mu\text{L}$	Dilution <input type="text"/> $\mu\text{L}$	Analytical Measuring Range Low <input type="text" value="1.50*"/>	High <input type="text" value="6.00*"/>	
R2 (R2-1) <input type="text" value="0"/> $\mu\text{L}$	Dilution <input type="text" value="0"/> $\mu\text{L}$	Correlation Factor A <input type="text" value="1"/>	B <input type="text" value="0"/>	
Common Reagent Type <input type="text" value="None"/>	Name <input type="text" value="None"/>	Manufacturer Factor A <input type="text" value="1"/>	B <input type="text" value="0"/>	
Wavelength Pri <input type="text" value="600"/> nm	Sec <input type="text" value="800"/> nm	Onboard Stability Period <input type="text" value="90"/> Day	<input type="text" value="0"/> Hour	
Method <input type="text" value="END"/>		LIH Influence Check <input type="text" value="#"/>		
Reaction Slope <input type="text" value="+"/>		Lipemia <input type="text" value="++++"/>		
Measuring Point-1 1st <input type="text" value="0"/>	Last <input type="text" value="1"/>	Icterus <input type="text" value="++++"/>		
Measuring Point-2 1st <input type="text"/>	Last <input type="text"/>	Hemolysis <input type="text" value="+++"/>		
Linearity Limit <input type="text"/> %				
Lag Time Check <input type="text"/>				

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> ALB1U <input type="text"/> <b>Test No</b> <input type="text"/> <b>Type:</b> Serum <input type="text"/>				
Value/Flag <input type="text" value="#"/>	Level Low <input type="text" value="#"/>	High <input type="text" value="#"/>		
<b>Specific Ranges</b>				
	From	To	Other Type	Low High
<input type="checkbox"/> 1: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 2: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6: Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	<input type="text" value="None"/>	<input type="text" value="#"/> <input type="text" value="#"/>
7: Standard demographics				<input type="text" value="#"/> <input type="text" value="#"/>
8: Not within expected values				<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="g/dL*"/>	<input type="text" value="Select"/>	Decimal Places <input type="text" value="#"/>

Calibrators	General	ISE	Calculated Test	Range
<b>Test Name:</b> ALB1U <input type="text"/> <b>Type:</b> Serum <input type="text"/>				
	<input type="checkbox"/> Use Serum Cal.			
Calibration Type: <input type="text" value="AB"/>	Formula: <input type="text" value="Y=AX+B"/>	Counts: <input type="text" value="#"/>		
<Calibrator Parameters>				Slope Check <input type="text"/>
	Calibrator	OD	Conc	Range Low High
Point-1	# <input type="text"/>		†*	5.4* 9.6*
Point-2	<input type="text"/>			
Point-3	<input type="text"/>			
Point-4	<input type="text"/>			
Point-5	<input type="text"/>			
Point-6	<input type="text"/>			
Point-7	<input type="text"/>			
MB Type Factor <input type="text"/>	1-Point Calibration Point <input type="text"/>	<input type="checkbox"/> with Conc-0		Interval (RB) <input type="text" value="Lot"/>
				Interval (ACAL) <input type="text" value="Lot"/>
		Stability		
		Reagent Blank <input type="text" value="30"/> Day <input type="text" value="0"/> Hour		
		Calibration <input type="text" value="30"/> Day <input type="text" value="0"/> Hour		

# User Defined

† Beckman Coulter System Calibrator Cat. No: DR0070

\* Values set for working in g/dL. To work in SI units (g/L) multiply by 10