iQ Control/Focus Set

REF 800-3104

TABLE OF CONTENTS

English	2
Français (FR)	4
Deutsch (DE)	7
Italiano (IT)	10
Español (ES)	13
Português Portugal (PT-PT)	16
Dansk (DA)	19
Svenska (SV)	21
Norsk (NO)	23
Ελληνικά (EL)	25
日本語 (JP)	28
中文 (ZH-CN)	31
Lietuviškai (LT)	33
Magyar (HU)	36
Polski (PL)	39
Čeština (CZ)	42
Slovenčina (SK)	45
한국어 (KO)	48
Türkçe (TR)	50
Русский (RU)	53
eesti keel (ET)	56
Hrvatski (HR)	59
Български (BG)	62
中文 (ZH-TW)	65
Română (RO)	67
Србија (SR)	69
Latviešu (LV)	72
Українська (UK)	74
Português Brasil (PT-BR)	77
Tiếng Việt (VI)	80

iQ Control/Focus Set

REF 800-3104

Rx Only in the U.S.A.

INTENDED USE

iQ Control/Focus Set is intended for *in vitro* diagnostic use with iQ 200 and DxU Microscopy Series Analyzers.

PRINCIPLE

iQ Positive Control and iQ Focus are suspensions of fixed human red blood cells in a buffered, isotonic balanced solution. They are used to control and focus the iQ200 and DxU Microscopy Series Analyzers. iQ Positive Control is used as an abnormal, and iQ Negative Control is used as a normal control to verify accurate counting by the instrument. iQ Focus is used to check light level and focus the instrument.

Barcode labels that contain Lot number, Expiration date and Expected value are provided for each product. An accurate concentration count of these cells is certified and the expected counting reference value (accurate to $\pm 10\%$) is printed on each label.

WARNING AND PRECAUTIONS

- The iQ Control/Focus has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

GHS HAZARD CLASSIFICATION

iQ Positive Control,
iQ Focus

EUH208

May produce an allergic reaction.
reaction mass of:
5-chloro-2-methyl-4-isothiazolin
-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC#
220-239-6](3:1) < 0.05%

Safety Data Sheet is available at beckmancoulter.com/techdocs

NON-REACTIVE INGREDIENTS

iQ Negative Control: Water, salts and preservatives.

iQ Positive Control/iQ Focus: Human erythrocytes, water, salts, dyes and preservatives. The stabilized human erythrocytes are <1% in an aqueous suspension and are non-reactive for HbsAg, HCV 3.0 Ab, HIV 1 / 2, HIV-1 p24 Ag, RPR and ALT <82 IU/mL. Stabilized human erythrocytes are chemically fixed, which minimizes any hazardous characteristics.

STORAGE AND STABILITY

- Store iQ Focus and iQ Positive Control at 2° C to 8° C (35.6° F to 46.4° F)
- Store iQ Negative Control at 2° C to 28° C (35° F to 82.4° F)
- iQ Control/Focus is stable for eight months from the date of manufacture
- Open expiration: 30 days from bottle open date
- Do not freeze
- Do not store in direct sunlight
- Bring to room temperature before use
- Do not save leftover material or mix Control/Focus material from different bottles
- Do not use iQ Control/Focus after the expiration date printed on the bottle labels

MATERIALS PROVIDED

- 2 X 125 mL iQ Focus
- 1 X 125 mL iQ Positive Control
- 1 X 125 mL iQ Negative Control
- 1 X 40 iQ Positive Control Barcode Labels
- 1 X 40 iQ Negative Control Barcode Labels
- 1 X 50 iQ Focus Control Barcode Labels

MATERIALS REQUIRED, BUT NOT PROVIDED

- Approved tubes listed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Do not use preservative tubes.
- iQ200 or DxU Microscopy Series Analyzer
- Personal Protective Equipment

PROCEDURE

Follow the directions for control rack, focus, and daily control as described in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Shake the Positive Control and Focus prior to use. **DO NOT SHAKE THE NEGATIVE CONTROL.** The included barcodes inform the instrument of the lot number, expiration date, and expected value. If the results are within the allowable range, the date and time of the QC will appear on the “Last Micro QC” field of the Instrument screen. Daily control results and statistics can be printed from the iQ QC Review and QC Statistics sections.

iQ CONTROL/FOCUS SET NOMINAL ASSAY VALUES AND RANGES

CONTROL	NOMINAL ASSAY VALUE	RANGE
POS	1000 (particles/ μ L)	\pm 100
FOCUS	1100 (particles/ μ L)	\pm 100
NEG	0 (particles/ μ L)	+ 20

Important Note: The actual assay values for the iQ Control/Focus Set may change from LOT number to LOT number.

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

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Updated copyright year
Phiên bản AB, 08/2021
Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland +(353) (0) 65 683 1100