iQ Calibrator

REF 800-3103

TABLE OF CONTENTS

English	2
Français (FR)	4
Deutsch (DE)	6
Italiano (IT)	8
Español (ES)	10
Português Portugal (PT-PT)	12
Dansk (DA)	14
Svenska (SV)	16
Norsk (NO)	18
Ελληνικά (EL)	20
日本語 (JP)	22
中文 (ZH-CN)	24
Lietuviškai (LT)	26
Magyar (HU)	28
Polski (PL)	30
Čeština (CZ)	32
Slovenčina (SK)	34
한국어 (KO)	36
Türkçe (TR)	38
Русский (RU)	40
eesti keel (ET)	43
Hrvatski (HR)	45
Български (BG)	47
中文 (ZH-TW)	49
Română (RO)	51
Србија (SR)	53
Latviešu (LV)	55
Українська (UK)	57
Português Brasil (PT-BR)	59
Tiếng Việt (VI)	61

iQ Calibrator

REF 800-3103

Rx Only in the U.S.A.

INTENDED USE

iQ Calibrator is intended for *in vitro* diagnostic use with iQ 200 and DxU Microscopy Series Analyzers.

PRINCIPLE

iQ Calibrator is a suspension of fixed human red blood cells in a particulate-free solution. It is used to automatically calibrate the iQ200 and DxU Microscopy Series Analyzers. An accurate count is certified and the expected value (accurate to approximately 8.3%) is printed on each label. Barcodes are provided that contain Lot number, Expiration date and Expected value.

WARNING AND PRECAUTIONS

- The iQ Calibrator has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910.1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

GHS HAZARD CLASSIFICATION

iQ Calibrator

EUH208

May produce an allergic reaction.
reaction mass of:
5-chloro-2-methyl-4-isothiazolin
-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC#
220-239-6](3:1) < 0.05%

Safety Data Sheet is available at beckmancoulter.com/techdocs

NON-REACTIVE INGREDIENTS

iQ Calibrator: Human erythrocytes, water, salts, dyes and preservatives. The concentration of stabilized human erythrocytes is <1% in an aqueous suspension. Non-reactive for HbsAg, HCV 3.0 Ab, HIV 1 / 2, HIV-1 p24 Ag, RPR and ALT <82 IU/mL. The stabilized human erythrocytes are chemically fixed, which thus minimizes any biohazardous risks.

STORAGE AND STABILITY

- Store at 2° C to 8° C (35.6° F to 46.4° F)
- iQ Calibrator is stable for eight months from date of manufacture
- Open expiration: 24 hours from bottle open time
- Do not freeze
- Do not store in direct sunlight
- Each bottle is expected to yield single calibration and should not be used again
- Do not mix the calibrator from different bottles
- Bring to room temperature before use
- Do not use the solution after the expiration date printed on the bottle label

MATERIALS PROVIDED

- 4 X 125 mL iQ Calibrator
- 1 X 40 iQ Calibrator Barcode Labels

MATERIALS REQUIRED, BUT NOT PROVIDED

- Approved tubes listed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use. Do not use preservative tubes.
- iQ200 or DxU Microscopy Series Analyzer
- Personal Protective Equipment.

LIMITATION OF THE PROCEDURE

No interferences are reported.

PROCEDURE

Shake the bottle as directed in the iQ200 or DxU Microscopy Series Analyzer Instructions for Use before use and follow the directions for calibration. The included barcodes provide the lot number, expiration date, and expected value to the instrument. The passing results are displayed on screen as soon as they are available.

iQ CALIBRATOR NOMINAL ASSAY VALUES AND RANGES

CONTROL	NOMINAL ASSAY VALUE	RANGE
iQ CALIBRATOR	1200 (MEAN; particles/ μ L)	\pm 100

Important Note: The actual assay values for the iQ Calibrator may change from LOT number to LOT number.

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

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Updated copyright year
Phiên bản AB, 08/2021
Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland +(353) (0) 65 683 1100