COULTER Body Fluid Control

REF 628030
TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	5
Deutsche (de-de)	8
Italiano (it-it)	11
Español (es-es)	14
Português (Portugal) (pt-pt)	17
Dansk (da-dk)	20
Svenska (sv-se)	23
Norsk (nb-no)	26
Suomi (fi-fi)	29
Ελληνικά (el-gr)	32
日本語 (ja-jp)	35
中文 (中国) (zh-cn)	38
Lietuvių (lt-lt)	41
Magyar (hu-hu)	44
Polski (pl-pl)	47
Čeština (cs-cz)	50
Slovák (sk-sk)	53
한국어 (ko-kr)	56
Türkçe (tr-tr)	59
Русский (ru-ru)	62
Eestlane (et-ee)	66
Hrvatski (hr-hr)	69
Български (bg-bg)	72
中文 (台湾) (zh-tw)	75
Română (ro-ro)	78
Slovenščina (sl-si)	81
Srpski (sr-rs)	84
Latviešu (lv-lv)	87
Українська (uk-ua)	90
Português (Brasil) (pt-br)	93
Nederlands (nl-nl)	96
Tiếng Việt (vi-vn)	99
Қазақша (kk-kz)	102
Bahasa Indonesia (id-id)	105
REFERENCES	108

COULTER Body Fluid Control

REF 628030

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

Body Fluid Control is a hematology quality control material used to monitor the performance and verify the measuring range of the body fluid cycle of UniCel DxH Coulter Cellular Analysis Systems listed in the TABLE OF EXPECTED RESULTS, in conjunction with specific COULTER reagents. Refer to the System HELP or Instructions for Use.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor performance, and verify the measurable range. This product can also be used to establish your own laboratory mean.

INTENDED USER

This device is intended to be used by a laboratory professional.

SUMMARY AND PRINCIPLE

The Body Fluid Control is a reference product prepared from stabilized human erythrocytes to confirm and monitor instrument accuracy and precision performance by providing quantitative measurements during counting.

The COULTER Body Fluid Control TNC and RBC concentrations span the instrument measurable range. Results for each concentration are compared to the established expected range, to assess the instrument performance and to verify the measurable range.

REACTIVE INGREDIENTS

Stabilized and/or fixed erythrocytes from:

Human < 55% w/w

Avian < 1% w/w

NON REACTIVE INGREDIENTS

Buffers, Stabilizers and Preservatives

REAGENTS

Body Fluid Control consists of treated, stabilized human erythrocytes in a buffered, isotonic medium. The control also contains fixed erythrocytes to simulate leukocytes.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced material: Human, Avian and Ungulate.

Each human donor unit used in preparation of this material was tested by an FDA approved method for the presence of the antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive).

Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen.¹

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines.²

STORAGE, STABILITY, AND DISPOSAL

The Body Fluid Control is shipped in a thermally insulated container designed to keep it cool. Reagents in unopened packs are stable for up to 120 days (until the expiration date stated on the label) when stored according to requirements. It is not recommended to use the reagent after shelf life expiration indicated on the label. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. Opened vials are stable for 16 days assuming that the Instructions for Use section of the package insert is performed a maximum of 18 times within the 16 days. Storage of product with cap down (inverted) position may require additional mixing for complete resuspension of cellular components.

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly colored supernatant) is indicative of product deterioration. A slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

INSTRUCTIONS FOR USE: DAILY QUALITY CONTROL

Before analyzing Body Fluid Control, review the following to ensure accurate results are obtained:

MATERIALS PROVIDED

- Body Fluid Control.
- Table of EXPECTED RESULTS containing Assigned Values and EXPECTED RANGES.

MATERIALS REQUIRED, BUT NOT PROVIDED

- Sufficient COULTER reagents to complete analysis. For the reagents used on the system refer to the IFU specific for the UniCel DxH Coulter Cellular Analysis Systems in use.
- Online System HELP or Instructions for Use.

INSTRUMENT SET UP

Prepare your instrument: Follow the instructions for Setting Up Controls in System HELP or the Instructions for Use when a new lot of controls is received.

Analyzing Body Fluid Controls:

1. Remove Body Fluid Control tube(s) from refrigerator and warm at AMBIENT TEMPERATURE for 10 to 15 minutes.
2. After warming, mix by hand as follows:
 - A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

3. Repeat steps 2A to 2C.

NOTE: DO NOT use any other type of mechanical mixer.

4. Run the control using the Single Tube Presentation on the instrument. Run Body Fluid Control tubes in consecutive order beginning with Level 1 and ending with Level 3.
5. Return the tube(s) to the refrigerator **within 30 minutes**.
6. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

ASSIGNED VALUES AND EXPECTED RESULTS

Refer to TABLE OF EXPECTED RESULTS.

Mean values should recover within the ranges. The ranges established by Beckman Coulter should be used as a guideline. Each laboratory must establish its own criteria for acceptable results.

ASSIGNED VALUES are determined on validated systems using specific COULTER reagents. The values are confirmed by multiple analyses of the control.

EXPECTED RANGES include variation between instruments, and represent 95% confidence limits for well maintained instrument systems.

- ◆ DxH Reagents for use on DxH 800 / 600 / 900 / 690T
- ♻ DxH ECO Reagents for use on DxH 900 / 690T

Whether you use the Beckman Coulter Assigned Value or your own laboratory mean, the instrument is considered well maintained and operating correctly (in control) if:

- ≥95% of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend OUTSIDE the EXPECTED RANGE.

Before your current cell control lot(s) expire, perform the following on your new lot(s):

- Confirm that recovered values are within the TABLE OF EXPECTED RESULTS.
- OR
- Establish your own laboratory mean.

IMPORTANT:

Flags, Codes and Messages

- DO NOT ignore System messages (R).
- Level 3 has a TNC high enough to trigger the TNC Carryover System Message if the TNC of the following sample is less than or equal to 200 cell/ μ L or the WBC Carryover System Message if the WBC of the following sample is less than or equal to 0.20×10^9 cell/ μ L. Samples giving carryover system messages should be repeated.

PERFORMANCE LIMITS

Individual laboratories can expect better repeatability than that shown in the EXPECTED RANGE column. Refer to your System HELP and Instructions for Use for precision performance specifications for your UniCel DxH Coulter Cellular Analysis Systems.

The report generated by IQAP/eIQAP will include an assessment of the measuring range.

TROUBLESHOOTING PROCEDURE

The System HELP and Instructions for Use contain troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other Body Fluid Control and COULTER 6C Control lots currently in use.
- Previous Body Fluid Control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON CP-X Control analyses.
- Data from the last calibration.
- Data from current repeatability study (N=10) using a fresh whole blood specimen and performed according to System HELP and Instructions for Use.

- Data from the IQAP survey as well as data from your last instrument calibration, where appropriate.

PRODUCT AVAILABILITY

COULTER Body Fluid Control:

REF 628030

4 x 3.5 mL Level 1

4 x 3.5 mL Level 2

4 x 3.5 mL Level 3

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

Glossary of Symbols is available at beckmancoulter.com/techdocs by searching for PN C07163.

May be covered by one or more pat. - see www.beckmancoulter.com/patents

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

The Summary of Safety and Performance is available from the EUDAMED database: <https://ec.europa.eu/tools/eudamed>.

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	February 2025	Initial release to include documentation updates for compliance to European In Vitro Diagnostic Regulations (IVDR) (2017/746). No changes were made to the product.
AB	June 2025	Updated wording in EXPECTED RANGES

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.