β 2-microglobulin RIA KIT

REF IM1113

TABLE OF CONTENTS

English	2	Čeština	29
Français	5	Slovenčina	32
Deutsch	8	한국어	35
Italiano	11	Türkçe	37
Español	14	Русский	39
Ελληνικά	17	中文 ZH-TW	42
Lietuviškai	20	APPENDIX	44
Magyar	23		
Polski	26		

β2-microglobulin RIA KIT

REF IM1113

RADIOIMMUNOASSAY FOR THE IN VITRO DETERMINATION OF β2-MICROGLOBULIN IN HUMAN SERUM, PLASMA AND URINE For *in vitro* diagnostic use.

PRINCIPLE

The radioimmunoassay of β2-microglobulin is a competition assay. Samples and calibrators are incubated with 125I-labeled β2-microglobulin, as tracer, in antibody-coated tubes. After incubation, the liquid content of tubes is aspirated and the bound radioactivity is determined in a gamma counter. A standard curve is constructed and unknown values are obtained from the curve by interpolation.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the immunoassay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Tracer

DANGER

H360

May damage fertility or the unborn child.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned: Get medical advice/attention.

P405

Store locked up.

P501

Dispose of contents/container in accordance with local/national regulations
Boric Acid 1 - 2%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

Serum and plasma

- Collect blood in dry tubes or in tubes with EDTA.
- Separate serum or plasma from cells by centrifugation within 3 hours after collection.
- Serum and plasma samples may be stored at 2-8°C, if the assay is performed within 24 hours. For longer storage keep frozen at <-20°C (6 months maximum) after aliquoting so as to avoid repeated freezing and thawing. Wait until the samples are completely thawed and homogenize them before the assay. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator they must be diluted by zero calibrator or assayed using smaller sample volume, 20 µL.

Serum and EDTA plasma values for 15 samples (serum values ranging from 1.18 to 1.91 mg/L) were compared using the β2-microglobulin RIA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.9927[\text{serum}] - 0.1426$$

$$R = 0.9557$$

Urine

- Collect urine in glass or plastic vessels. Add borate (10-20 mM) for storage.
- Urine samples may be stored at 2-8°C, if the assay is performed within 24 hours. For longer storage keep frozen at <-20°C (1 year maximum) after aliquoting so as to avoid repeated freezing and thawing. Wait until the samples are completely thawed and homogenize them before the assay. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator they must be diluted by zero calibrator.
- If samples have concentrations lower than 0.25 mg/L they should be assayed using 100 or 200 µL sample volume.

Note: At acid pH (below 5), β2-microglobulin is denatured and can no longer be assayed.

MATERIALS PROVIDED

All reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8 °C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution are indicated in paragraph Procedure.

Anti-β2-microglobulin monoclonal antibody-coated tubes: 2 x 50 tubes (ready-to-use)

125I-labeled β2-microglobulin tracer: one 55 mL vial (ready-to-use)

The vial contains 148 kBq at the date of manufacture of 125I-labelled β2-microglobulin in buffer with bovine serum albumin, sodium azide (<0.1 %) and a dye.

Calibrators: one 2 mL vial, five 0.5 mL vials (ready-to-use)

The calibrator vials contain from 0 to approximately 30 mg/L of β2-microglobulin in buffer with bovine serum albumin and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated using the international standard WHO, 1st IS 1985. 1 IU corresponds to 14 ng.

Control serum: one vial (lyophilised)

The vial contains β2-microglobulin lyophilised in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipets (50 µL).
- Semi-automatic pipets (500 µL).
- Vortex type mixer.
- Horizontal or orbital shaker.
- Aspiration system.
- Gamma counter set for ¹²⁵I

RESULTS

The results in the package insert were calculated using a logit-log curve fit (weighted cubic regression) with B/T (%) or B/B₀ (%) on vertical axis and the β2-microglobulin concentration of the calibrators on the horizontal axis (mg/L). Other data reduction methods may give slightly different results.

Standard curve

Total activity: 51,671 cpm				
Calibrators	β2-m (mg/L)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0	42,192	81.6	100
1	0.29	34,240	66.3	81.2
2	0.76	20,614	39.9	48.9
3	1.90	10,463	20.2	24.8
4	7.60	3,520	6.81	8.34
5	28.5	1,384	2.68	3.28

(Example of standard curve, do not use for calculation)

Samples

Locate for each sample the ratio B/T (%) or B/B₀ (%) on the vertical axis of the standard curve and read-off the corresponding β2-microglobulin concentration of the sample on the horizontal axis in mg/L. The concentrations of diluted samples must be corrected by the dilution factor. When different volume than 50 µL was taken, correction must be made on sample results, too.

EXPECTED VALUES

It is suggested that each laboratory establish its own normal values. The following values obtained with healthy subjects are indicative only.

Plasma and serum:

1.0 to 2.4 mg/L (95% of the normal population); mean 1.20 mg/L. The serum level increases slightly with age but it is not dependent on sex.

Urine

Up to 0.37 mg/L

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature.

Reconstitution of control serum

The content of the vial reconstitute with the volume of distilled water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solution at 2-8°C for 3 days or aliquoted at <-18 °C until the expiry date of the kit.

Immunoassay procedure

Step 1 Additions	Step 2 Incubation	Step 3 Counting
To antibody coated tubes, add successively: 50 µL of calibrator, control or sample and 500 µL of tracer* Mix.	Incubate 90 minutes at (18-25 °C) with shaking (>280 rpm).	Aspirate carefully the content of tubes (except of the 2 tubes «total cpm»): Count bound cpm (B) and total cpm (T) for 1 min.

*Add 500 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.06 mg/L

Functional sensitivity: 0.185 mg/L

Specificity

No cross-reactivity with human IgG was detected in the assay.

Precision

Intra-assay

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 6.8 % for serum and below or equal to 6.4 % for urine.

Inter-assay

Samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 10.8 % for serum and below or equal to 12.8 % for urine.

Accuracy

Dilution test

High-concentration samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 85.3 % and 119 % for serum and between 81.0 % and 116 % for urine.

Recovery test

Samples were spiked with known quantities of β2-microglobulin. The recovery percentages were obtained between 97.6 % and 119 % for serum and between 80.7 % and 97.3 % for urine.

Measurement range (from analytical sensitivity to highest calibrator):

0.06 to approximately 30 mg/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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