## HISTAMARINE ENZYME IMMUNOASSAY KIT IM2369

The “Histamarine” enzyme immunoassay kit is intended for the quantitative determination of histamine concentrations in fresh, frozen or canned fish by fish processors or quality control laboratories.

It has been shown by the AOAC Research Institute-independent laboratory evaluation that the Histamarine test kit detects the presence of histamine in tuna, canned tuna and mahi-mahi.

The intended user should be familiar with enzyme immunoassays.

### PRINCIPLE OF THE ASSAY

Histamine contained in a fish sample is determined by use of high affinity specific monoclonal antibodies coated to the wells in a microtiter plate. The amount of histamine bound to antibodies is measured by use of enzyme conjugate which catalyzes a color change in the plate wells. The intensity of color is inversely proportional to the concentration of histamine in the sample. The concentration is calculated by use of a calibration curve obtained with calibrators supplied with the kit.

### MATERIALS REQUIRED BUT NOT PROVIDED

Precision balance, blender, microtiter plate reader (405-414 nm).

Variable volume precision micropipets (20 µL, 200 µL, 1000 µL), 25 mL and 1 L graduate cylinders.

Plastic tubes and distilled water.

### REAGENTS PROVIDED IN THE KIT

Let the components of the kit **equilibrate** to room temperature before use.

All reagents in the unopened kit are stable at 2-8 °C until the expiration date of the kit.

Reagents of the kit may be labelled with IVD symbol, this is due to manufacturer internal purposes only. Do not take the IVD symbol into account. This kit is intended for the quantitative determination of histamine in fresh, frozen or canned fish by fish processors and quality control laboratories, not for use in diagnostic procedures.

Reagent	Description	Reagent preparation	Storage after use*	Precautions
<b>Plate</b>	96 wells	ready-to-use	store at 2-8 °C in self-lock bag	
<b>Calibrators</b>	5 x 1 mL	ready-to-use	store at 2-8 °C	DANGER H316 H360 
<b>Acylation reagent</b>	2 vials (powder)	<b>dissolve</b> with DMSO (volume on label)	store at -20 °C	<b>Homogenize</b> prior use.
<b>DMSO</b>	2 x 3 mL	ready-to-use	store at 2-8 °C, <b>melt</b> at room temperature	WARNING H227 H315 H319 H335 
<b>Acylation buffer</b>	4 x 5 mL	ready-to-use	store at 2-8 °C	DANGER H360 
<b>Conjugate</b>	2 vials lyophilized	<b>dissolve</b> with d.w.** (volume on label) Once dissolved, pour both vials into Vial for conjugate and mix.	<b>See Vial for conjugate</b>	When reconstituted <b>wait</b> 20 min before mixing. DANGER H316 H360 
<b>Wash solution</b>	1 x 50 mL	<b>dilute</b> with 950 mL of d.w.**	store 1 month at 2-8 °C (beyond at -20 °C)	DANGER H360 

<b>Substrate buffer</b>	1 x 30 mL	ready-to-use	store at 2-8 °C	DANGER H316 H318 H373  
-------------------------	-----------	--------------	-----------------	--

Reagent	Description	Reagent preparation	Storage after use*	Precautions
<b>Substrate</b>	2 tablets of pNPP	<b>dissolve</b> 1 tablet in 15 mL of substrate buffer	store 24 h at 2-8 °C (beyond at -20 °C)	Handle with forceps.
<b>Stop solution</b>	1 x 6 mL of NaOH	ready-to-use	store at 2-8 °C	DANGER H314 
<b>Vial for conjugate</b>	1 vial	Mix 2 vials of reconstituted conjugate	store 1 week at 2-8 °C (beyond at -20 °C)	DANGER H316 H360 

\* until the expiration date      \*\* distilled water

*Note : Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.*

#### ASSAY PROCEDURE

Conduct fish sampling per local regulatory requirements.

Extraction	Samples	Immunoassay	Results
<p>Weigh precisely the fish sample (between 1 and 10 grams).</p> <p>Crush to homogeneity in 8 mL of water per gram of fish in a blender.</p> <p>Centrifuge (10,000 x g 5 min) or filter (paper filter) or decant (carefully collect supernatant).</p>	<p>In plastic tubes containing 180 µL of acylation buffer dispense 20 µL of extraction supernatant.</p> <p>Add 50 µL of acylation reagent, mix.</p> <p>At this step, samples may be stored 48 h at 2-8 °C.</p>	<p>To plate wells add 50 µL of sample or of ready to use calibrator.</p> <p>Add 200 µL of conjugate.</p> <p>Incubate 30 min at 18-25 °C.</p> <p>Wash wells carefully.</p>	<p>To plate wells add 200 µL of substrate and incubate 30 min at 18-25 °C.</p> <p>Add 50 µL of stop solution.</p> <p>Read plate at 405-414 nm.</p>

## RESULTS

Results are obtained by interpolation from calibration curve.

Calibrator	Histamine ( $\mu\text{M}$ )	Absorbance	B/B0 (%)
1	0.1	1.754	100
2	1.0	0.907	47.7
3	5.0	0.384	15.4
4	11.0	0.275	8.7
5	50.0	0.165	1.9

*Example of calibration curve, not to be used in calculation.*

Using the absorbance of the 5 calibrators provided with the kit, draw a calibration curve on semi-log paper (or use appropriate software), plotting on the horizontal axis (log scale) the histamine concentration of the 5 calibrators and on the vertical axis (linear scale) the corresponding absorbance value obtained for each calibrator.

When the points have been plotted for each histamine calibrator, connect points resulting in a calibration curve. An example of how the curve should appear is presented in the table above.

Using the calibration curve, locate absorbance on vertical axis and read off corresponding histamine concentrations on horizontal axis.

A concentration of 1  $\mu\text{M}$  read on the horizontal axis corresponds to a histamine concentration of 10 ppm (10 mg/kg) in the original tissue sample. This correspondency takes into account dilutions performed during sample preparation.

## PERFORMANCE

The sensitivity, defined as the lowest concentration of acylated histamine significantly different from zero with a probability of 95 %, is 0.03  $\mu\text{M}$  (0.3 ppm).

At the 1  $\mu\text{M}$  concentration, 1-methyl histamine, 3-methyl histamine, putrescine and histidine do not show cross-reactivity higher than 0.05 % in the assay.

Recovery experiments performed on fish samples gave an average recovery of 103.5 %.

Intra-assay

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 8.91 %.

Inter-assay

Samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 17.16 %.

## TECHNICAL ASSISTANCE

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: [imunochem@beckman.com](mailto:imunochem@beckman.com).

**For the first test run please refer to our detailed Instructions (upon request)**

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