

SYNCHRON LX[®]i Maintenance Log

Serial #: _____ System ID: _____ Month: _____ Year: _____

LX[®]i DAILY Access 2

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Check Zone Temperatures																																
Check System Supplies																																
Check Liquid Waste Container																																
System Backup Successful?																																
Inspect Fluidic Module																																
Clean Probe Exteriors																																
Prime Substrate																																
Run Daily Clean System																																
(If Necessary) Perform Special Clean																																
Tech Initials																																

TWICE WEEKLY LX Maintenance

	Week 1 <i>Date / Initial</i>	Week 2 <i>Date / Initial</i>	Week 3 <i>Date / Initial</i>	Week 4 <i>Date / Initial</i>
Clean				
• Sample & reagent probes, mixers, EIC and flow cell ² using Clenz Solution and Sodium Hypochlorite.	_____	_____	_____	_____
• Probes ⁰ (Exterior)	_____	_____	_____	_____
• CAP Piercer Blades ⁰ (only systems with 4-blade foil CTS option)	_____	_____	_____	_____

SYNCHRON LX[®]i Maintenance Log

Serial #: _____ System ID: _____ Month: _____ Year: _____

LX[®]i WEEKLY Access 2

	Week 1	Week 2	Week 3	Week 4	Week 5
Today's Date					
Clean Instrument Exterior					
Check Waste Filter Bottle					
Clean Primary Probe					
Remove/Clean Aspirate Probes					
Run Daily Maintenance (Initials above)					
Run System Check					
Tech Initials					

Enter the System Check Results in the boxes. Initial the Tech Initials box. Access 2

SYSTEM CHECK RESULTS	Week 1	Week 2	Week 3	Week 4	Week 5
Today's Date					
Washed RLU/%CV					
Substrate RLU/%CV					
Unwashed RLU/%CV					
Wash Efficiency PPM					
Substrate Ratio					
Substrate: Washed Ratio					
Tech Initials					

System Check Expected Results	
Washed Check	
RLU mean	5,000 - 20,000
% CV	≤ 12.00
Substrate Check	
RLU mean	5,000 - 9,000
% CV	≤ 5.00
Substrate Ratio	0 - 1.40
Substrate: Washed Ratio	0 - 1.00*
Unwashed Check	
RLU mean	4 - 10 million*
% CV	≤ 2.00
Wash Efficiency	
PPM	0 - 5.00
<p style="font-size: small;">* The Substrate: Washed Ratio and the Unwashed Check RLU mean result are not system specifications. They are only reference guidelines.</p>	

LX[®]i WEEKLY SYNCHRON LX[®]

	Week 1 <i>Date / Initial</i>	Week 2 <i>Date / Initial</i>	Week 3 <i>Date / Initial</i>	Week 4 <i>Date / Initial</i>	Week 5 <i>Date / Initial</i>
Check					
• Level of Wash Concentrate II, No Foam and CTS Auto-Gloss	_____	_____	_____	_____	_____
• Replace cuvette wiper ¹	_____	_____	_____	_____	_____
Clean					
• Total Protein and Albumin Cup Modules ²	_____	_____	_____	_____	_____

SYNCHRON LX®i Maintenance Log

Instrument Serial #: _____ System ID: _____ Month: _____ Year: _____

**LX®i TWICE-MONTHLY
SYNCHRON Maintenance - For CTA Refer to IFU**

	<i>Date/Initial</i>	<i>Date/Initial</i>
Clean • CTA Piercing Probe, Aliquot Probe, Auto-Gloss Dispense Station, Probe Wash Station ⁰	_____	_____

**LX®i MONTHLY
SYNCHRON Maintenance**

	<i>Date/Initial</i>	<i>Date/Initial</i>
Replace • Alkaline buffer reagent and in-line filter, when completed, check and adjust damper fluid level ²	_____	_____
Clean • Albumin, BUNm/UREAm (including electrode), Creatinine, Glucose, Phosphorus, and Total Protein cup modules ² • CC sample and reagent mixers ⁰ • Flow Cell as per Flow Cell Maintenance Procedure using Clenz Solution • Chloride Electrode	_____ _____ _____	_____ _____ _____
Calibrate • CUPs Lamp/Sensor ²	_____	_____

**LX®i TWO-MONTH
SYNCHRON Maintenance**

	<i>Date</i>	<i>Initial</i>	<i>Date</i>	<i>Initial</i>	<i>Date</i>	<i>Initial</i>	<i>Date</i>	<i>Initial</i>	<i>Date</i>	<i>Initial</i>	<i>Date</i>	<i>Initial</i>
Replace • CTS Blade/Wick ⁰ (only systems with 1-Blade Thick CTS)												
• Blade ⁰ (Only Systems with Sarstedt S-Monovette [®] CTS Option)												
Change/clean • Air Filters ⁰												

SYNCHRON LX[®]i Maintenance Log

Serial #: _____ System ID: _____ Month: _____ Year: _____

**LX[®]i THREE-MONTH
SYNCHRON Maintenance Manual Reference Chapter 6**

	Date	Initial	Date	Initial	Date	Initial	Date	Initial
Replace • Sample (100µL) and reagent (500µL) syringe plunger rod with attached tips ¹								

FOUR-MONTH MAINTENANCE MANUAL REFERENCE CHAPTER 7

	Date	Initial
Replace • Blades ⁰ (Only Systems with 4-Blade Foil CTS Option)		
Clean • MC Reagent Lines, Cups and Stir Bars ²		

**ELECTRODES/SENSOR/PARTS REPLACEMENT
SIX-MONTH SYNCHRON Maintenance**

	Date / Initial	Electrodes/Sensor Number	Date / Initial	Electrodes/Sensor Number
Replace • Potassium electrode tip (P/N 669117) ² • Calcium electrode tip (P/N 467769) ² • AccuSense™ glucose sensor (P/N 467432) ²	____	____-____-____	____	____-____-____
	____	____-____-____	____	____-____-____
	____	____-____-____	____	____-____-____

**LX[®]i ONE-YEAR
SYNCHRON Maintenance - For CTA Refer to IFU**

	Date	Initial	Date	Initial	Date	Initial	Date	Initial
Replace • CTA Aliquot Probe and Piercing Probe ⁰ • CTA Syringes ⁰								

Following specific maintenance procedures, calibration and controls may be required. A code is displayed for each procedure to show the minimum requirements following maintenance: 0 = No action, 1 = Run Controls**, 2 = Calibrate and Run Controls.

*Maintenance procedure frequencies may vary depending upon individual laboratory workload volumes. These frequencies are based on the number of samples per day an average laboratory may process, approximately 400 samples/day.

**It is assumed that if controls are out, a Level 2 will occur.

