
DOC026.98.80406, Edition 7
April 2026

HIAC HRLD+

© 2026 Beckman Coulter, Inc.

Specifications

Specifications are subject to change without notice. For sensor configurations, refer to Table 1.

Specification	Details
Collection optics	Light-obscuration
Light source	Class III B, 780 nm, 5 mW laser diode
Sample temperature	65 °C (150 °F) maximum
Liquid sample types	Samples with fire points above 90 °C[a]
Power requirements	±15 VDC (furnished from the counter)
Counter control and signals	Cable provided with sensor must be used
Operating environment	7 to 52 °C (44.6 to 125.6 °F)
	30 to 95% relative humidity, non-condensing
Storage environment	0 to 71 °C (32 to 159.8 °F)
	0 to 98% relative humidity, non-condensing
Dimension (W x L x H)	203.2 x 38.1 x 38.1 mm (8 x 1.5 x 1.5 in.)
Weight	0.5 kg (1.1 lb)
Certifications	HRLD Laser Module Classification: Class 1 Laser Product
	US FDA-CDRH Laser Accession No. 8921539-005
	European Standard: EN 60825-1 and EN 61010-1
	CE mark
Protection class	III
Installation category	I
Pollution degree	1
Enclosure	For indoor use only

a. Note: This is based on sample temperature of 65 °C + 25 °C (safety factor). Special considerations apply to samples with fire points below 90 °C.

Table 1 Sensor configurations

Model	Range (µm)	Particles/mL[a]	Flow rates (mL/min)	Calibrated flow rate (mL/min)
HRLD 100	4 µm(c)–70 µm(c)[b]	10,000	20–100	60
HRLD 100HC	4 µm(c)–70 µm(c)[b]	18,000	10–50	20
HRLD 150	1.3–150	18,000	10–50	10 or 25
HRLD 150JA	1.3–150	18,000	10–50	10 or 25
HRLD 400	2.0–400	10,000	20–100	60
HRLD 400HC	2.0–400	18,000	10–50	20
HRLD 600JS	2.0–600	6000	30–200	100

a. Concentration at <10% coincidence error

b. As per ISO 11171 and ISO 4406 (NIST Suspension SRM 2806). Note that 4 µm(c) is less than 1 µm per ISO 4402 (ACFTD in oil).

Table 2 Fluid compatibility

Group 1: HRLD 100/100HC/150/400/400HC	Group 2: HRLD 150JA/600JS (Group 1, plus)
Hydraulic fluids, Stoddard solvent, jet fuel (JP4, JP5), phosphate esters[a], kerosene, purified water, water for injection, diesel fuel, alcohols	Aldehydes, aromatics, hydrochloric acid, ammonium hydroxide, ketones, sulfuric acid, sodium hydroxide, esters, phosphoric acid, hydrogen peroxide

a. Viton™ seals standard. Kalrez™ seals available for phosphate esters.

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTE The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTE Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A

on the instrument is referenced in the manual with a precautionary statement.


	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the need for protective eye wear.
	This symbol indicates a laser device is used in the equipment.
	This symbol indicates the presence of devices sensitive to Electro-static Discharge (ESD) and indicates that care must be taken to prevent damage with the equipment.
UK CA	A “UKCA” mark indicates that a product has been assessed before being placed in the UK market, and has been found to meet UK safety, health, and/or environmental protection requirements.

Product overview

The HIAC HRLD+ sensor uses light obscuration technology to sense particles in liquids and can be used with a wide range of HIAC counters. The primary use of the sensor is for liquid contamination and particle size distribution analysis applications. Refer to Figure 1.The sensors are supported by Models 8000A, 2000, 9703 and 9703+, which are necessary to supply particle count and size data for the user.

The sensors supply an accurate measurement of particulate contamination in virtually any liquid, except certain solvents and acids. Samples with a very high level of particle contaminants must be diluted.

Figure 1 Product overview


1. Shipping plugs
2. Outlet port
3. Microcell
4. Power and signal connector
5. Inlet port

System configuration

Refer to Figure 2 for the sensor installation overview with external equipment.

Figure 2 System configuration


1. Counter
2. Sensor control and signal cable
3. Sensor
4. Drain line
5. Sample flow
6. Waste container
7. Sample
8. Sampler
9. Sampler control cable

Product components


Make sure that all components have been received. Refer to Figure 3 and Figure 4. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 Product components (150 and 150-JA models)

1. Sensor
2. Sensor clearing shim

Figure 4 Product components (all other models)

1. Sensor **2.** Cleaning floss

Installation

DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

WARNING

Personal injury hazard. Never remove covers from the instrument. This is a laser-based instrument and the user risks injury if exposed to the laser.

NOTE The laser diode is extremely sensitive to static charges and out of tolerance voltage variations. Never make or remove any connections with the line power applied to the counter. The counter provides power ±15 VDC to the sensor.

NOTE Do not put the sensor in an area where it is susceptible to electronic noise and mechanical vibration because this will damage the sensor.

Make sure that the sensor is adjacent to the particle counter at a distance less than the length of the interconnecting communication/electrical cable. Make sure that the sensor is adjacent to each of the associated liquid sampling and/or handling equipment so that the tubing length is kept to a minimum. Sufficient access for operation, maintenance and testing must be provided.

Make sure to use the external equipment documentation for the correct mechanical and electrical procedures for those devices. If necessary, contact the manufacturer or sales representative to resolve any questions of compatibility or suitability of this product for a specific application.

Install the sensor

The sensor has a flow direction indicator. Make sure that the flow direction on the sensor is monitored during installation (Figure 2).

1. Remove the two shipping plugs from the microcell.

 NOTE To prevent contamination, install the shipping caps when the sensor is not connected to the system.

2. Connect the sensor adapter fitting from the sample source to the sensor inlet port (Figure 1).

3. Connect the sensor outlet port (Figure 1) to the drain or to the flow meter tubing inlet (Figure 2).

4. Make sure to flush the cell and sample lines completely with clean, particle-free water immediately after making the cell connections. Make sure that there are no leaks in the sampling system to obtain accurate samples.

5. Make sure that all fittings are properly connected and leak-free. Do not over-tighten any fittings. Finger-tighten, then an additional ¼-turn.

Connect power to the sensor

NOTE Damage to the sensor can occur if connected or disconnected while powered. Always remove power from the instrument before making electrical connections.

1. Set the counter power to off.

2. Connect the sensor cable connector to the sensor connector on the particle counter. The sensor is ready to use.

3. Set the counter power to on.

Operation

DANGER

Fire/explosion hazard. Do not use the instrument with samples or liquids that have an ignition point below 90 °C (194 °F).

WARNING

Personal injury hazard. Never remove covers from the instrument. This is a laser-based instrument and the user risks injury if exposed to the laser.

DANGER

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Make sure to flush the cell and sample lines completely with clean, particle-free water immediately before use. Make sure that there are no leaks in the sampling system to get accurate samples. Make sure that the sample is compatible with water before use.

NOTE Cover removal will void the warranty.

The sensor operates automatically with the connected particle counter and the sampler systems. To optimize the sensor performance:

- Use standard particle counters with the sensor such as the HIAC Model 8000A or 9064. When a sensor is used with one of these particle counters, the applicable calibration algorithm is interpolation (curve fit). Do not use the equation algorithm.
- Make sure to flush the tubing and the microcell with clean miscible liquids between different sample analyses.
- Powder sample: Make sure to select an applicable liquid to prepare a powder suspension for the analyses. For a stable dispersion, use a surfactant. Ultra-sonification just prior to sample analysis is recommended.
- Keep the sample at ambient temperature. Cold liquids can cause water condensation on the microcell window.
- Make sure that the sensor microcell is kept clean when the sensor is not in use. If the microcell is not clean and is allowed to dry, particle residue forms on the microcell windows and causes bad interference with the laser beam. This becomes difficult to clean.
- Make sure that all samples are compatible with the components used in the sample path.

- Do not let particles go in the microcell if the dimensions of the particles are close to the maximum cross section of the microcell. The microcell dimension (width) is very small.

NOTE Clean the microcell if number counts are too high in the first channel for a "clean" liquid, or if an over concentration is found. Refer to *Clean the sensor and the microcell*.

Condition the sensor

When oil applications are used, the sensor flow cell (microcell) must be conditioned for use with oils. When the application is water soluble or cleaning functions are used, the sensor flow cell must be conditioned for water.

Contact technical support for assistance about the miscibility of a liquid and its compatibility with the sensor.

NOTE All solutions used to condition the sensor for oil to water and water to oil activities must be particle free liquids.

To change the application from water to oil, complete this procedure. Otherwise, erroneous counts will result.

Items to collect:

- Clean water
- Alcohol (isopropyl)
- Filtered mineral spirits
- Clean particle free oil

1. Flush the sensor with clean water.
2. Flush the sensor with alcohol (typically isopropyl alcohol).
3. Flush the sensor with filtered mineral spirits several times.
4. Flush the sensor with clean (particle free) oil before the actual sample analysis with oil. The sensor is ready for oil application.
5. To change from an oil to an aqueous application, complete the above procedure in the reverse order (i.e., flush the sensor with filtered mineral spirits, then flush with alcohol, and last flush with water). The basis for this procedure is that one fluid is used after the other must be miscible. Immiscibility will make droplets and cause false particle counts, which may collect layers on the windows of the cell.

Maintenance

DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Clean the sensor and the microcell

DANGER

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

DANGER

Personal injury hazard. Never remove covers from the instrument. This is a laser-based instrument and the user risks injury if exposed to the laser.

NOTE Damage to the sensor can occur if connected or disconnected while powered. Always remove power from the instrument before making electrical connections.

NOTE Potential Instrument Damage. Delicate internal electronic components can be damaged by static electricity, resulting in degraded performance or eventual failure.

Sensor body: Examine the sensor periodically for debris and deposits. Clean the sensor when there is a build-up of deposits. Use a clean, damp, lint-free soft cloth to remove loose debris from the sensor body.

Microcell: Clean the microcell monthly and also when:

- a clean sample goes through the microcell and returns very high counts are noted in the lower channels
- an over concentration indications are shown on the counter
- a Baseline Hi sensor message is found

Clean the HRLD+ 150/150-JA Sensor cell

NOTE If a facility air source is used to clean the sensor, make sure that the facility air system does not contain suspended oil or produce water droplets that will interfere with sensor cleaning.

1. Set the counter power to off.

2. Purge the sensor of any sample and flush the sensor with a suitable miscible liquid. If necessary, prepare the sensor microcell for water miscible liquids.

3. Remove the sensor from the sampler.

4. Blow either clean dry filtered air or "canned air" through the outlet port of the sensor for about 1 second to dry the microcell. Refer to Figure 5.

5. Make sure that the sensor microcell is not blocked. Refer to Figure 6. If the microcell slit is clogged, refer to Clean a blocked microcell.

6. Install the sensor and make sure that the orientation of the inlet and outlet ports is reversed. The flush flow should be in the opposite direction of the sample flow through the sensor.

7. Flush with clean (filtered and deionized) water repeatedly.

 NOTE Use a 5% solution of Contrad 70 for periodic cleaning and use a 50% solution of Contrad 70 for heavily contaminated sensors. If this does not solve the problem, flush the sensor with alcohol.

8. After the sensor microcell is flushed with the cleaning liquid, rinse the microcell with clean particle free water.

9. Install the sensor on the sampler configured in the original orientation (inlet/outlet ports). Make sure that the system has no leaks and, as necessary, prepare the sensor for oil applications.

10. If necessary, repeat this procedure.

 NOTE If the cleaning process has been completed three times without satisfactory results, contact technical support for a sensor microcell overhaul.

At the completion of this procedure, refer to the counter and sampler operation manuals for the proper start-up procedure.

Figure 5 Air jet to clean the microcell



Figure 6 Look at the rectangular slit

Clean the HRLD+ 100/400/600 Sensor cell

NOTE If a facility air source is used to clean the sensor, make sure that the facility air system does not contain suspended oil or produce water droplets that will cause problems with sensor cleaning.

For HRLD+ 100, 100HC, 400, 400HC and 600JS Sensors

1 Set the counter power to off.

2 Purge the sensor of any sample and flush the sensor with an applicable miscible liquid. If necessary, prepare the sensor microcell for water miscible liquids. The manufacturer recommends filtered mineral spirits, petroleum ether, Contrad 70 or Stoddard Solvent. Contact technical support for the applicability of other solvents. Flush the microcell with alcohol.

3 Remove the sensor from the sampler.

4 Blow either clean dry filtered air or "canned air" through the outlet port of the sensor for about 15 seconds to dry the microcell. Refer to Figure 5.

5 Make sure that the sensor microcell is not blocked. Refer to Figure 6. If the microcell slit is clogged, refer to Clean a blocked microcell.

6 Mix clean water and the Micro™ cleaner in a washer bottle. Refer to Figure 7.

7 Push the white end of the cleaning floss leader through the DRY outlet side of the sensor microcell until the leader is seen at the inlet. It may be necessary to turn the floss to get it past the microcell orifice and out the inlet. Refer to Figure 8.

8 Pour the mixed Contrad 70 into the sensor port cavity of the flow cell until the flow cell is full of the solution. Never use the floss dry. A diluted micro solution or deionized water must always be on the floss when it is used.

NOTE Never pull the waxed cleaning floss end through the microcell. The wax could leave a deposit on the microcell windows which is difficult to remove.

9 Rinse the microcell with clean particle free water.

10 Install the sensor on the sampler configured in the original orientation (inlet/outlet ports). Make sure that the system has no leaks. Prepare the sensor for oil applications, if necessary.

11 Repeat this procedure, if necessary.

NOTE If the cleaning process has been completed three times without satisfactory results, contact technical support for a sensor microcell overhaul.

Refer to the counter and sampler operation manuals for the proper start-up procedure.

Figure 7 Prepare cleaning solution



Figure 8 Clean the sensor with floss

Clean a blocked microcell


Do the steps that follow to clean a blocked microcell.

1. Remove the sensor.
2. Install the sensor in the opposite direction of the flow path (Figure 9).
3. Flush the microcell with 20 mL of clean, particle-free water that is miscible with the sample solution in the microcell. If the sampler stops while the microcell is flushed, proceed to step 5.
4. Use a clean dry filtered air jet to blow air in the opposite direction of the sample flow through the microcell (Figure 5). This helps to remove loose particles that may be stuck in the microcell.
5. Carefully insert a clearing shim in the outlet port of the sensor and move the clearing shim toward the inlet port, through the microcell. This action should push the blockage out of the microcell.
6. Install the sensor in the correct flow path.

7 Repeat steps 3 and 4.

If the problem still occurs, contact technical support.

Figure 9 Reverse flow installation

Troubleshooting

Error/Warning	Description	Solution
Sampler stalls	Viscosity is too high.	Dilute the sample with a miscible solvent (e.g., kerosene). Warm the sample, but not more than the sample flash point and flammability specifications.
	The sensor microcell is clogged.	Refer to Clean a blocked microcell.
	Bubbles in the sample or a leak in the fluid path.	Tighten fittings or replaced damaged tubing or failed fittings.
Baseline Hi Error message	Sample is too dark. Sample has a high particle concentration.	Dilute the sample.
	Sensor microcell is clogged or dirty.	Refer to Clean a blocked microcell.
High noise level, with clean fluids, but no sample flow	Sensor microcell is clogged or dirty.	Refer to Clean a blocked microcell.
High noise level, with flow only	The “Clean sample” is dirty.	Make sure to do a filtered blank sample.
	Turbulence/bubbles occur from debris clogged in the sensor microcell flow channel.	Refer to Clean a blocked microcell.

Error/Warning	Description	Solution
Bad distribution is exhibited with polystyrene spheres.	Turbulence/bubbles occur from debris clogged in the sensor microcell flow channel.	Refer to Clean a blocked microcell.
	The sensor microcell is clogged or dirty.	Refer to Clean a blocked microcell.

Operational errors and solutions

Error/Warning	Description	Solution
Fluids are left in the microcell for extended periods of time	Fluids leave a residue that will coat the windows of the microcell.	• Do not use samples more than one time. • Make sure that the particles in a sample are smaller than the cell width (500 µm x 1000 µm). • Refer to Clean the sensor and the microcell. • Use only approved cleaning equipment to clean the sensor flow channel. Refer to Replacement parts and accessories.
Particles that are too large or are too irregularly shaped to go through the sensor	The microcell clogs and can create particle build up behind the clog, which disrupts or blocks the flow thorugh the microcell. If particles clog at the center of the microcell, then the laser light may not move to the photo-diode. This can cause turbulence and bubble formation.	

Replacement parts and accessories

⚠ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

NOTE Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
HRLD+ 100 sensor	2089829-100
HRLD+ 100HC sensor	2089829-100HC
HRLD+ 150 sensor	2089829-150

Description	Item no.
HRLD+ 150JA sensor	2089829-150JA
HRLD+ 400 sensor	2089829-400
HRLD+ 400HC sensor	2089829-400HC
HRLD+ 600JS sensor	2089829-600JS

Accessories

Description	Item no.
Accessory kit, In-line adapter kit fittings	690-0001
Cleaning kit or Sensor clearing shim	80003108
Cleaning floss	660-850-0001

Beckman Coulter, Inc.
250 S. Kraemer Blvd.
Brea, CA 92821 U.S.A.

EC REP
Beckman Coulter Eurocenter S.A.
22, rue Juste-Olivier
Case Postale 1044
CH - 1260 Nyon 1, Switzerland
Tel: +41 (0) 22 365 36 11


Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other trademarks are the property of their respective owners.

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