



CELL LAB Rat Anti-Mouse CD5/Lyt-1

<u>Cat. No.</u>	<u>Form</u>	<u>Quantity</u>
732007	Purified (UNLB) Antibody	0.5 mg
732008	Fluorescein (FITC) Conjugate	0.5 mg
732009	Fluorescein (FITC) Conjugate	0.1 mg
732010	Biotin (BIOT) Conjugate	0.5 mg
732011	Phycoerythrin (PE) Conjugate	0.1 mg
733262	Phycoerythrin (PE) Conjugate	0.2 mg
732012	Allophycocyanin (APC) Conjugate	0.1 mg
732013	Spectral Red™ (SPRD) Conjugate	0.1 mg

For Laboratory Use Only

DESCRIPTION

Clone: B19.1
Isotype: Rat IgG2a κ
Specificity: CD5/Lyt-1, Mr 67 kDa

CD5/Lyt-1 antigen is a monomeric type I transmembrane glycoprotein expressed on thymocytes, T lymphocytes, and a subset of B lymphocytes, but not on natural killer (NK) cells.¹⁻³ It has been identified as the major ligand of the B-cell antigen CD72.^{4,5} The frequency of CD5⁺ B cells exhibits strain-dependent variation, and the phenotypic, anatomical, functional, developmental, and pathological characteristics of the CD5⁺ B cells suggest that they may represent a distinct lineage, known as B-1 cells.⁶ Binding of CD5 on the T cell surface can augment alloantigen-or mitogen-induced lymphocyte proliferation and induces increased cytosolic free calcium, IL-2 secretion, and IL-2R expression.⁷⁻¹² It has been proposed that CD5 negatively regulates signal transduction mediated by the T-cell and B-cell receptors.^{13,14}

APPLICATIONS

- Identification and enumeration of CD5⁺ cells by flow cytometry
- Identification of CD4⁺ cells in acetone-fixed frozen sections
- Immunoprecipitation
- Activation of T cells *in vitro*
- Western blotting

CHARACTERIZATION

To ensure lot-to-lot consistency, each batch of product is tested to conform with characteristics of a standard reference reagent using immunofluorescence staining and flow cytometry.

WORKING DILUTIONS

Flow Cytometry:

FITC conjugate	≤1 μg/10 ⁶ cells
BIOT conjugate	≤1 μg/10 ⁶ cells
PE conjugate	≤0.2 μg/10 ⁶ cells
APC conjugate	≤0.2 μg/10 ⁶ cells
SPRD conjugate	≤0.2 μg/10 ⁶ cells

Other Applications: Since applications vary, determine the optimum working dilution of the product that is appropriate for your specific needs.

HANDLING AND STORAGE

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of 100 mM borate buffered saline, pH 8.0. No preservatives or amine-containing buffer salts added.
- The fluorescein (FITC) conjugates are supplied as 0.5 mg or 0.1 mg in 1.0 mL of PBS/NaN₃.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN₃.
- The phycoerythrin (PE) conjugates are supplied as 0.1 mg in 1.0 mL or 0.2 mg in 2.0 mL of PBS/NaN₃ and a stabilizing agent.
- The allophycocyanin (APC) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent.
- The Spectral Red (SPRD) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN₃ and a stabilizing agent.
- Protect fluorochrome-conjugated forms from light. Do not freeze.
- Reagent is stable until the expiration date on the vial when stored at 2-8°C.

STATEMENT OF WARNINGS

1. Specimens, samples and all material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
2. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
3. Do not use reagent beyond the expiration date on the vial label.
4. Minimize exposure of reagent to light during storage or incubation.
5. Avoid microbial contamination of reagent or erroneous results may occur.
6. Use Good Laboratory Practice (GLP) when handling this reagent.
7. Harmful if swallowed.
8. After contact with skin, wash immediately with plenty of water.
9. Contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, immediately wash excessively with water.

TRADEMARKS

The Beckman Coulter logo is a trademark of Beckman Coulter, Inc.
Spectral Red is a trademark of Southern Biotechnology Associates, Inc.
Cy5 is a trademark of GE Healthcare, Inc.

For additional information or if damaged product is received, contact your local Beckman Coulter Representative.

Spectral Red is a PE/Cy[™]5 tandem conjugate. Cy5 is for non-commercial research use only, not for therapeutic or in vivo applications. Other use needs license from Amersham Biosciences Corp., under U.S. Patent Nos. 4,981,977 and 5,268,486 and other patents pending. This material (or portions of this material) is subject to proprietary rights of Amersham Biosciences Corp. and Carnegie Mellon University and made and sold under license from Amersham Biosciences Corp. This product is licensed for sale only for research. It is not licensed for any other use. There is no implied license hereunder for any commercial use. Commercial use shall include: 1) sale, lease, license or other transfer of the material or any material derived or produced from it 2) sale, lease, license or other grant of rights to use this material or any material derived or produced from it 3) use of this material to perform services for a fee for third parties. If you require a commercial license to use this material and do not have one, return this material, unopened to Beckman Coulter, Inc. 11800 SW 147 Ave. Miami, FL 33196, USA and any money paid for the material will be refunded.

REFERENCES

1. Ledbetter JA and Herzenberg LA. 1979. Xenogeneic monoclonal antibodies to mouse lymphoid differentiation antigens. *Immunol Rev*, 47:63-90.
2. Seaman WE, Wofsy D, Greenspan JS and Ledbetter JA. 1983. Treatment of autoimmune MRL/lpr mice with monoclonal antibody to Thy-1.2: a single injection has sustained effects on lymphoproliferation and renal disease. *J Immunol*, 130:1713-1718.
3. Ledbetter JA, Rouse RV, Micklem HS and Herzenberg LA. 1980. T cell subsets defined by expression of Lyt-1,2,3 and Thy-1 antigens. Two-parameter immunofluorescence and cytotoxicity analysis with monoclonal antibodies modifies current views. *J Exp Med*, 152:280-295.
4. Van de Velde H, von Hoegen, I, Luo W, Parnes JR and Thielemans K. 1991. The B-cell surface protein CD72/Lyb-2 is the ligand for CD5. *Nature*, 351:662-665.
5. Luo W, Van de Velde H, von Hoegen I, Parnes JR and Thielemans K. 1992. Ly-1 (CD5), a membrane glycoprotein of mouse T lymphocytes and a subset of B cells, is a natural ligand of the B cell surface protein Lyb-2 (CD72). *J Immunol*, 148:1630-1634.
6. Kantor AB and Herzenberg LA. 1993. Origin of murine B cell lineages. *Annu Rev Immunol*, 11:501-538.
7. Hollander N. 1982. Effects of anti-Lyt antibodies on T-cell functions. *Immunol Rev*, 68:43-66.
8. Stanton T, Stevens TL, Ledbetter JA and Wofsy D. 1986. Anti-Ly-1 antibody induces interleukin 2 release from T cells. *J Immunol*, 136:1734-1737.
9. Verwilghen J, Vandesande R, Vandenberghe P and Ceuppens JL. 1990. Crosslinking of the CD5 antigen on human T cells induces functional IL2 receptors. *Cell Immunol*, 131:109-119.
10. Logdberg L and Shevach EM. 1985. Role of the Ly 1 antigen in interleukin 1-induced thymocyte activation. *Eur J Immunol*, 15:1007-1013.
11. June CH, Rabinovitch PS and Ledbetter JA. 1987. CD5 antibodies increase intracellular ionized calcium concentration in T cells. *J Immunol*, 138:2782-2792.
12. Hollander N, Pillemer E and Weissman IL. 1981. Effects of Lyt antibodies on T-cell functions: augmentation by anti-Lyt-1 as opposed to inhibition by anti-Lyt-2. *Proc Natl Acad Sci USA*, 78:1148-1151.
13. Tarakhovskiy A, Kanner SB, Hombach J, Ledbetter JA, Muller W, Killeen N and Rajewsky K. 1995. A role for CD5 in TCR-mediated signal transduction and thymocyte selection. *Science*, 269:535-537.
14. Bikah G, Carey J, Ciallella JR, Tarakhovskiy A and Bondada S. 1996. CD5-mediated negative regulation of antigen receptor-induced growth signals in B-1 B cells. *Science*, 274:1906-1909.



Manufactured for:
Beckman Coulter, Inc.
4300 N. Harbor Blvd.
Fullerton, CA 92835
www.beckmancoulter.com

Printed in USA
Made in USA

© 2005 Beckman Coulter, Inc.
All Rights Reserved.

PN 733913-A