



## CELL LAB Rat Anti-Mouse CD86/B7-2

Cat. No.	Form	Quantity
732514	Purified (UNLB) Antibody	0.5 mg
732515	Fluorescein (FITC) Conjugate	0.5 mg
732516	Biotin (BIOT) Conjugate	0.5 mg
732517	Phycoerythrin (PE) Conjugate	0.1 mg
734742	Allophycocyanin (APC) Conjugate	0.1 mg

### For Laboratory Use Only

#### DESCRIPTION

**Clone:** 2D10  
**Isotype:** Rat IgG2b $\kappa$   
**Specificity:** Mouse CD86, the B7-2 co-stimulatory molecule and a ligand for CD28 and CD152

CD86/B7-2 is an 80 kDa type I integral membrane protein and a member of the immunoglobulin superfamily of cell surface receptors. Along with CD80/B7-1, CD86 participates in T cell co-stimulation via interaction with CD28 and CD152/CTLA-4. It is expressed at very low density in normal tissues but is rapidly upregulated on stimulated B cells, monocytes and dendritic cells. Since CD86/B7-2 has rapid kinetics of induction, it is believed to be the major CD28 ligand expressed early in the immune response.<sup>1-9</sup>

#### APPLICATIONS

- Identification and enumeration of CD86<sup>+</sup> cells by flow cytometry<sup>1,3,4,9</sup>
- Identification of CD86<sup>+</sup> cells in acetone-fixed, frozen sections<sup>8,9</sup>
- Immunoprecipitation<sup>3</sup>
- *In vitro* blocking assays<sup>1,2,8,9</sup>

#### CHARACTERIZATION

To ensure lot-to-lot consistency, each batch of product is tested to conform with characteristics of a standard reference reagent using immunofluorescence staining and flow cytometry.

#### WORKING DILUTIONS

**Flow Cytometry:**

FITC conjugate	$\leq 1 \mu\text{g}/10^6$ cells
BIOT conjugate	$\leq 1 \mu\text{g}/10^6$ cells
PE conjugate	$\leq 0.2 \mu\text{g}/10^6$ cells
APC conjugate	$\leq 0.2 \mu\text{g}/10^6$ cells

**Other Applications:** Since applications vary, determine the optimum working dilution of the product that is appropriate for your specific needs.

#### HANDLING AND STORAGE

- The purified (UNLB) antibody is supplied as 0.5 mg of purified immunoglobulin in 1.0 mL of 100 mM borate buffered saline, pH 8.0. No preservatives or amine-containing buffer salts added.
- The fluorescein (FITC) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>.
- The biotin (BIOT) conjugate is supplied as 0.5 mg in 1.0 mL of PBS/NaN<sub>3</sub>.

- The phycoerythrin (PE) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN<sub>3</sub> and a stabilizing agent.
- The allophycocyanin (APC) conjugate is supplied as 0.1 mg in 1.0 mL of PBS/NaN<sub>3</sub> and a stabilizing agent.
- Protect fluorochrome-conjugated forms from light. Do not freeze.
- Reagent is stable until the expiration date on the vial when stored at 2-8°C.

## STATEMENT OF WARNINGS

1. Specimens, samples and all material coming in contact with them should be handled as if capable of transmitting infection and disposed of with proper precautions.
2. Never pipet by mouth and avoid contact of samples with skin and mucous membranes.
3. Do not use reagent beyond the expiration date on the vial label.
4. Minimize exposure of reagent to light during storage or incubation.
5. Avoid microbial contamination of reagent or erroneous results may occur.
6. Use Good Laboratory Practice (GLP) when handling this reagent.
7. Harmful if swallowed.
8. After contact with skin, wash immediately with plenty of water.
9. Contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be flushed with running water while being discarded. These precautions are recommended to avoid deposits in metal piping in which explosive conditions can develop. If skin or eye contact occurs, immediately wash excessively with water.

## TRADEMARKS

The Beckman Coulter logo is a trademark of Beckman Coulter, Inc.

For additional information or if damaged product is received, contact your local Beckman Coulter Representative.

## REFERENCES

1. Hathcock KS, Laszlo G, Pucillo C, Linsley P and Hodes RJ. 1994. Comparative analysis of B7-1 and B7-2 co-stimulatory ligands: expression and function. *J Exp Med*, 180:631-640.
2. Freeman GJ, Borriello F, Hodes RJ, Reiser H, Hathcock KS, Laszlo G, McKnight AJ, Kim J, Du L, and Lombard DB. 1993. Uncovering of functional alternative CTLA-4 counter-receptor in B7-deficient mice. *Science*, 262:907-909.
3. Laszlo G, Hathcock KS, Dickler HB and Hodes RJ. 1993. Characterization of a novel cell-surface molecule expressed on subpopulations of activated T and B cells. *J Immunol*, 150:5252-5262.
4. Hathcock KS, Laszlo G, Dickler HB, Bradshaw J, Linsley P and Hodes RJ. 1993. Identification of an alternative CTLA-4 ligand co-stimulatory for T cell activation. *Science*, 262:905-907.
5. Thompson CB. 1995. Distinct roles for the co-stimulatory ligands B7-1 and B7-2 in T helper cell differentiation? *Cell*, 81:979-982.
6. June CH, Bluestone JA, Nadler LM and Thompson CB. 1994. The B7 and CD28 receptor families. *Immunol Today*, 15:321-331.
7. Han S, Hathcock K, Zheng B, Kepler TB, Hodes R and Kelsoe G. 1995. Cellular interaction in germinal centers. Roles of CD40 ligand and B7-2 in established germinal centers. *J Immunol*, 155:556-567.
8. Inaba K, Witmer-Pack M, Inaba M, Hathcock KS, Sakuta H, Azuma M, Yagita H, Okumura K, Linsley PS, Ikehara S, Muramatsu S, Hodes RJ and Steinman RM. 1994. The tissue distribution of the B7-2 co-stimulator in mice: abundant expression on dendritic cells in situ and during maturation in vitro. *J Exp Med*, 180:1849-1860.
9. Larsen CP, Ritchie SC, Hendrix R, Linsley PS, Hathcock KS, Hodes RJ, Lowry RP and Pearson TC. 1994. Regulation of immunostimulatory function and co-stimulatory molecule (B7-1 and B7-2) expression on murine dendritic cells. *J Immunol*, 152:5208-5219.



Manufactured for:  
 Beckman Coulter, Inc.  
 4300 N. Harbor Blvd.  
 Fullerton, CA 92835  
[www.beckmancoulter.com](http://www.beckmancoulter.com)

Printed in USA  
 Made in USA

© 2005 Beckman Coulter, Inc.  
 All Rights Reserved.

PN 733991-A