

COULTER® GEN•S™ HEMATOLOGY ANALYZER

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
IDENTIFY POTENTIAL HAZARDS AND SAFETY FEATURES								
Review MSDS	www.beckmancoulter.com MSDS manual in lab	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Recognize Hazard Symbols	On-line Help	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
IDENTIFY AND BRIEFLY EXPLAIN FUNCTION OF ANALYZER COMPONENTS	Training Manual On-line Help	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
OVERVIEW OF SCREENS AND SOFTWARE	Training Manual On-line Help							
Identify Workstation Application Buttons		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Use HELP (Contents, Index, Search, Help mode button)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Buttons on Common Toolbar of Results & Graphics window		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Buttons on Specific Toolbar of Results & Graphics window		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Buttons on Common Toolbar of QA window		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Buttons on Specific Toolbar of QC window		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Identify Buttons on Specific Toolbar of Daily Checks window		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

COULTER® GEN•S™ HEMATOLOGY ANALYZER

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
MANAGE REAGENTS AND CONSUMABLES	Training Manual On-line Help							
Setup new lot numbers of Latron, 5C and Retic C controls		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Replace reagent containers		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
STARTUP AND SHUTDOWN PROCEDURES	Training Manual On-line Help							
Perform and review Startup		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform Shutdown		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
QA/QC PROCEDURES								
Demonstrate proper handling of Latron, 5C, Retic C	Control package inserts	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run Latron, 5C, Retic C	Training Manual On-line Help Control package inserts	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Review Latron, 5C, Retic C results	Laboratory SOP	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Download QC files for IQAP	Training Manual On-line Help	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Perform Calibration	Training Manual On-line Help S-Cal package inserts	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

COULTER® GEN•S™ HEMATOLOGY ANALYZER

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
PROCESS SAMPLES AND QUALITY CONTROL MATERIALS	Training Manual On-line Help							
Load cassette with sample/control		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run Latron Primer and Control in Manual/Secondary		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run 5C and Retic C in the Automatic/Primary Mode		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run patient sample in the Automatic/Primary Mode		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run patient sample in the Manual/Secondary Mode		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Run patient sample using pre-dilute function		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
AS NECESSARY CLEANING PROCEDURES	Training Manual On-line Help							
Clean and Replace Air Filters		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Clean BSV (Blood Sampling Valve)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Clean Drip Plate		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Clean Probe Wipe		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Clean Needle Assembly		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

COULTER® GEN•S™ HEMATOLOGY ANALYZER

Task	Resource	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Clean Stripper Plate		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
AS NECESSARY REPLACEMENT PROCEDURES	Training Manual On-line Help							
Replace needle		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Replace tubing		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			