

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name iQ Calibrator
Part number 800-3103
Series name iQ 200 and DxU Microscopy Series Analyzers

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For In Vitro Diagnostic Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare
Ireland

Supplier

Beckman Coulter, Inc.
250 S. Kraemer Blvd
Brea, CA 92821, U.S.A.
Tel: 800-854-3633

CANADA
Beckman Coulter Canada LP
7075 Financial Drive
Mississauga, ON L5N 6V8
Canada
1-800-463-7828

UNITED KINGDOM
Beckman Coulter (UK) Ltd.
Amersham Place
Little Chalfont
Buckinghamshire
United Kingdom, HP7 9NA
01494 441181

AUSTRALIA
Beckman Coulter Australia Pty Ltd
23-27 Chaplin Drive
Lane Cove NSW 2066
Australia
ABN 81 002 011 672
24 Hour emergency contact phone
number:
1800 060 881

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

SWITZERLAND
Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale
1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11
Monday through Friday, 9:00 am to
7:00pm)

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637

SAFETY DATA SHEET

Document ID: 300-7651 Version AN
Revision Date (year/month/day) 2025/06/08
Last Revision Date (year/month/day) 2024/10/30

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

Telephone: +356 21343270
Hours available: 08:30 – 17:00

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare
Ireland
Tel: 353 (0)65 6831100

SDSNT@beckman.com

e-mail address

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description

Mixture
Yellow-orange; Liquid; Odorless

Classification according to EC 1272/2008 (CLP/GHS)

Not classified as hazardous per EC 1272/2008 (CLP/GHS)

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS

Not classified as hazardous per US-OSHA HCS 2012 and UN GHS

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS Hazardous ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)

Pictogram

None

SAFETY DATA SHEET

Document ID: 300-7651 Version AN
Revision Date (year/month/day) 2025/06/08
Last Revision Date (year/month/day) 2024/10/30

Section 2 Hazards identification (Continued)

Signal word

None

Hazard statements

EUH208 May produce an allergic reaction.

Precautionary statements

Prevention

None

Response

None

Storage

None

Disposal

None

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

This product contains material of human origin and should be considered as potentially capable of transmitting infectious diseases.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note
reaction mass of: 5-chloro-2-methyl-4-isothiazolin -3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) CAS # 55965-84-9 EINECS # Not available Index # 613-167-00-5	< 0.05	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 M-factor Acute=100 M-factor Chronic=100 Skin Corr. 1C, H314 Skin Sens. 1A, H317 EUH071 Specific Concentration Limit (SCL) Skin Irrit. 2 H315 >= 0.06% - < 0.6% Eye Dam. 1 H318 >= 0.6% Skin Corr. 1C H314 >= 0.6% Eye Irrit. 2 H319 >= 0.06% - < 0.6% Skin Sens. 1A H317 >= 0.0015% Acute Toxicity Estimates (ATE) ATE Dermal = 87.12 mg/kg ATE Inhalation - Vapors = 0.5 mg/L ATE Oral = 53 mg/kg	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 Skin Corr. 1C, H314 Skin Sens. 1A, H317	10, 15, 8

Section 3 Composition and information on ingredients (Continued)

- 10 - Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6] (3:1) is the active ingredient of ProClin 150.
15 - May produce an allergic reaction.
8 - Present at concentration below the cut-off limits.

See section 8 for available Occupational exposure limits
See Section 15 for additional regulatory information
See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

- Inhalation** If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.
- Eye contact** If product enters eyes, rinse eyes gently with water as a precaution.
- Skin contact** In case of skin contact, rinse with water as a precaution.
- Ingestion** If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

4.2 Most important symptoms and effects, both acute and delayed

This product contains a sensitizing substance below concentration limit, may produce an allergic reaction in some people. Refer Section 3.

4.3 Indication of any immediate medical attention and special treatment needed

No specific medical attention or treatment required.

Section 5 Firefighting measures

- 5.1 Extinguishing media** In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam. For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

No special hazards determined.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product (an aqueous solution).

5.3 Advice for firefighters

Protective equipment Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

SAFETY DATA SHEET

Document ID: 300-7651 Version AN
 Revision Date (year/month/day) 2025/06/08
 Last Revision Date (year/month/day) 2024/10/30

Section 6 Accidental release measures (Continued)

- | | |
|---|--|
| Personal precautions | This product contains material of human origin and should be handled as though capable of transmitting infectious diseases. Observe general safety guidelines for protection during clean up procedures.
Wear protective gloves, protective clothing and eye/face protection. |
| 6.2 Environmental precautions | Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.
Dispose of contents/container in accordance with local regulations |
| 6.3 Methods and material for containment and cleaning up | |
| Spill and leak procedures | As a precautionary measure, treat spilled material with a 1:10 bleach/water solution. Absorb liquid and place in container suitable for disposal. Avoid generation of aerosols during clean up. Comply with applicable waste disposal regulations. |
| 6.4 Reference to other sections | Refer sections 8 and 13. |

Section 7 Handling and storage

- | | |
|---|--|
| 7.1 Precautions for safe handling | This product should be handled as though capable of transmitting infectious diseases. Universal precautions should be followed when using this product. |
| 7.2 Conditions for safe storage, including any incompatibilities | Store at 2 to 8°C, as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10). |
| 7.3 Specific end uses | No further relevant information available. |

Section 8 Exposure controls and personal protection

- | | |
|--|------------------|
| 8.1 Control parameters | |
| Exposure limits | |
| US OSHA | None established |
| ACGIH | None established |
| ACGIH Biological Exposure Indices (BEI) | None established |
| DFG MAK | None established |
| Ireland | None established |
| IOELVs | None established |
| NIOSH | None established |
| China | None established |
| Croatia | None established |

SAFETY DATA SHEET

Document ID: 300-7651 Version AN
Revision Date (year/month/day) 2025/06/08
Last Revision Date (year/month/day) 2024/10/30

Section 8 Exposure controls and personal protection (Continued)

Japan	None established
Sweden (AFS 2015:7 and amendments)	None established
Turkey	None established
8.2 Exposure controls	
Engineering controls	No special engineering controls are required. Use with good general ventilation.
Eye protection	Safety glasses or chemical goggles should be worn to prevent eye contact. Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.
Skin protection	Wear protective clothing and impervious gloves, as appropriate.
Respiratory protection	Under normal conditions, the use of this product should not require respiratory protection.

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Liquid	Density and/or relative density	≈ 1.007
Color	Yellow-orange	Solubility	
Odor	Odorless	Water	Miscible
pH	≈ 7.0	Organic	Not determined
Freezing point	Not determined	Partition coefficient n-octanol/water (log value)	Not determined
Boiling point or initial boiling point and boiling range	Not determined	Auto-ignition temp.	Not applicable
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	Not determined
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not applicable		
Relative vapor density	Not determined		
Particle characteristics	Not applicable		

9.2 Other information

Information with regard to physical hazard classes

Not applicable

Section 9 Physical and chemical properties (Continued)

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

- 10.1 Reactivity** No further relevant information available.
- 10.2 Chemical stability** The product is stable in accordance with recommended storage conditions.
- 10.3 Possibility of hazardous reactions**
No further relevant information available.
- 10.4 Conditions to avoid** To maintain product performance keep away from strong acids, strong bases, strong oxidizers.
Avoid exposure to heat and direct sunlight.
- 10.5 Incompatible materials** No further relevant information available.
- 10.6 Hazardous decomposition products**
No decomposition products posing significant hazards would be expected from this product (an aqueous solution).

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)
CAS # 55965-84-9

Dermal LD50 Rabbit 87.12 mg/kg (ECHA_API); Oral LD50 Rat 53 mg/kg (NLM_CIP)

- Primary routes of exposure** Common routes of entry include inhalation, ingestion and eye/skin contact. Specific paths of concern for potentially infectious materials are skin puncture, contact with broken skin, contact with mucous membranes and inhalation of aerosolized material.
- Acute toxicity** Not classified based on available data.
- Skin corrosion/irritation** Not classified based on available data.
- Serious eye damage/irritation** Not classified based on available data.
- Respiratory or skin sensitisation** This product contains a sensitizing substance below concentration limit, may produce an allergic reaction in some people. Refer Section 3.
- Germ cell mutagenicity** Not classified based on available data.
- Carcinogenicity** No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.
- Reproductive toxicity** Not classified based on available data.

Section 11 Toxicological information (Continued)

Specific target organ toxicity (STOT) – single exposure

Not classified based on available data.

Specific target organ toxicity (STOT) – repeated exposure

Not classified based on available data.

Aspiration hazard

Not classified based on available data.

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information

This product contains material of human origin and should be considered as potentially capable of transmitting infectious diseases.

Section 12 Ecological information

12.1 Toxicity

Fresh water species

No information available.

Microtox/organisms

No information available.

Water flea

No information available.

Fresh water algae

No information available.

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulative potential Not determined for the product.

12.4 Mobility in soil Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects

This product contains environmentally hazardous substance below the cutoff level. Refer section 3 for ingredient information. Do not allow undiluted product to enter sewer/surface or ground water.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Dispose of as potentially biohazardous waste and in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we

SAFETY DATA SHEET

Document ID: 300-7651 Version AN
Revision Date (year/month/day) 2025/06/08
Last Revision Date (year/month/day) 2024/10/30

	recommend that you contact the relevant (local) authorities and/or and approved waste-disposal company for information.
Package disposal	Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.
Additional information	Suggested European waste catalogue 18 01 03* - wastes whose collection and disposal is subject to special requirements in order to prevent infection. Dispose in accordance with national, state and local waste regulations

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

- 14.1 **UN/ID number:** Not regulated for transportation
- 14.2 **UN proper shipping name:** Not regulated for transportation
- 14.3 **Transport hazard class(es):** Not regulated for transportation
- 14.4 **Packing group:** Not regulated for transportation
- 14.5 **Environmental hazards:** Not regulated for transportation
- 14.6 **Special precautions for user:** None
- 14.7 **Maritime transport in bulk according to IMO instruments:** Not applicable

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

Section 15 Regulatory information (Continued)

No ingredients listed.

Massachusetts Right To Know (RTK) List

CAS # 10377-60-3 Magnesium Nitrate

New Jersey Dept. of Health Right To Know (RTK) List

CAS # 10377-60-3 Magnesium Nitrate

Pennsylvania Right To Know (RTK) List

CAS # 10377-60-3 Magnesium Nitrate

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany)

WGK 1, low water endangering

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

Not applicable.

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product is exempt from WHMIS label and SDS requirements.

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

CAS # 10377-60-3 Magnesium Nitrate

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Section 15 Regulatory information (Continued)

Türkiye-REACH - KKDIK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Safety Ratings	Flammability: 0 Health: 1 Reactivity with water: 0 Physical contact: 1	Code 0=None 1=Slight 2=Caution 3=Severe
-----------------------	---	---

Revision changes Update supplier addresses in Section 1.3

Document version and issue/revision date

Revision Date (year/month/day) 2025/06/08
Last Revision Date (year/month/day) 2024/10/30
Document ID: 300-7651
Version: AN

Hazard Classification Procedure

This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.

Description of hazard class and hazard statements from Section 3

Aquatic Acute 1 - Aquatic Hazard Acute, Category 1
Acute Tox. Dermal 2 - Acute Toxicity Dermal, Category 2
Acute Tox. Inhal. 2 - Acute Toxicity Inhalation, Category 2
Acute Tox. Oral 3 - Acute Toxicity Oral, Category 3
Eye Dam. 1 - Eye Damage Category 1
Aquatic Longterm 1 - Aquatic Hazard Long term, Category 1
Skin Corr. 1C - Skin Corrosion Category 1C
Skin Sens. 1A - Skin Sensitization Category 1A
EUH071 - Corrosive to the respiratory tract.
H301 - Toxic if swallowed.
H310 - Fatal in contact with skin.
H314 - Causes severe skin burns and eye damage.
H317 - May cause an allergic skin reaction.
H318 - Causes serious eye damage.
H330 - Fatal if inhaled.
H400 - Very toxic to aquatic life.
H410 - Very toxic to aquatic life with long lasting effects.

Section 16 Other information (Continued)

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
US OSHA - United States Occupational Safety and Health Administration
vPvB - very Persistent and very Bioaccumulative substances
WHMIS - Workplace Hazardous Material Information System

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.
