

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name MicroScan 0.8% Sulfanilic Acid Reagent

Part number 6668-30/W, B1010-44A, B1015-44

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For In Vitro Diagnostic Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter, Inc.
250 S. Kraemer Blvd
Brea, CA 92821, U.S.A.
Tel: 800-854-3633

Supplier

CANADA
Beckman Coulter Canada LP
7075 Financial Drive
Mississauga, ON L5N 6V8
Canada
1-800-463-7828

UNITED KINGDOM
Beckman Coulter (UK) Ltd.
Amersham Place
Little Chalfont
Buckinghamshire
United Kingdom, HP7 9NA
01494 441181

AUSTRALIA
Beckman Coulter Australia Pty Ltd
23-27 Chaplin Drive
Lane Cove NSW 2066
Australia
ABN 81 002 011 672
24 Hour emergency contact phone
number:
1800 060 881

SWITZERLAND
Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale
1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11
Monday through Friday, 9:00 am to
7:00pm)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare
Ireland
Tel: 353 (0)65 6831100

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

e-mail address SDSNT@beckman.com

1.4 Emergency telephone number

Telephone number (24H) Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description Mixture
Colorless; Liquid; Pungent

Classification according to EC 1272/2008 (CLP/GHS)
Skin Corrosion Category 1, H314
Eye Damage Category 1, H318

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS
Skin Corrosion Category 1
Eye Damage Category 1
Skin Sensitization Category 1

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS
Hazardous ingredients

Sulphanilic Acid
Acetic Acid

Pictogram



SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 2 Hazards identification (Continued)

Signal word

DANGER

Hazard statements

H314 Causes severe skin burns and eye damage.

H317 May cause an allergic skin reaction.

Precautionary statements

Prevention

P261 Avoid breathing vapours.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves, protective clothing and eye/face protection.

Response

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P303+P361+P353 IF ON SKIN (or hair): Rinse skin with water.

P304+P340 IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before use.

P363 Wash contaminated clothing before reuse.

Storage

P405 Store locked up.

Disposal

P501 Dispose of contents/container in accordance with local/national regulations

Product label will display most significant precautionary statements.

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note

SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 3 Composition and information on ingredients (Continued)

<p>Acetic Acid CAS # 64-19-7 EINECS # 200-580-7 Index # 607-002-00-6</p>	<p>25 - 35</p>	<p>Flam. Liq. 3, H226 Skin Corr. 1A, H314 Specific Concentration Limit (SCL) Skin Irrit. 2 H315 >= 10% - < 25% Skin Corr. 1A H314 >= 90% Skin Corr. 1B H314 >= 25% - < 90% Eye Irrit. 2 H319 >= 10% - < 25%</p>	<p>Flam. Liq. 3, H226 Skin Corr. 1A, H314</p>
<p>Sulphanilic Acid CAS # 121-57-3 EINECS # 204-482-5 Index # 612-014-00-X</p>	<p>0.1 - 1</p>	<p>Eye Irrit. 2, H319 Skin Irrit. 2, H315 Skin Sens. 1, H317</p>	<p>Eye Irrit. 2, H319 Skin Irrit. 2, H315 Skin Sens. 1, H317</p>

See section 8 for available Occupational exposure limits
See Section 15 for additional regulatory information
See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

- Inhalation** If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.
- Eye contact** If product enters eyes, rinse eyes gently with water for 15 minutes or longer, making sure that the eyelid is held open. Obtain medical advice/attention.
- Skin contact** In case of skin contact, rinse with plenty of water for at least 15 minutes. Remove contaminated clothing and shoes. Obtain medical advice/attention.
- Ingestion** If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

4.2 Most important symptoms and effects, both acute and delayed

Corrosive. May cause burns to eyes, skin and respiratory tract.
May cause an allergic skin reaction.
See Section 11 Toxicological Information for more detailed health information.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available. Refer to Section 4.1.

Section 5 Firefighting measures

- 5.1 Extinguishing media** In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam.
For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

No special hazards determined.

Section 5 Firefighting measures (Continued)

Hazardous combustion products

No combustion products posing significant hazards are expected from this product (an aqueous solution).

5.3 Advice for firefighters

Protective equipment

Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information

No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions

Observe general safety guidelines for protection; avoid eye and skin contact. Wear protective gloves, protective clothing and tightly sealed eye/face protection.

6.2 Environmental precautions

Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.

6.3 Methods and material for containment and cleaning up

Spill and leak procedures

Absorb spilled material with an appropriate inert, non-flammable absorbent and dispose according to local regulations.

6.4 Reference to other sections

Refer sections 8 and 13.

Section 7 Handling and storage

7.1 Precautions for safe handling

Use good laboratory procedures; avoid eye and skin contact.

7.2 Conditions for safe storage, including any incompatibilities

Store at 15 to 30°C, as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).

7.3 Specific end uses

No further relevant information available.

Section 8 Exposure controls and personal protection

8.1 Control parameters

Exposure limits

US OSHA

Acetic Acid
CAS # 64-19-7 10 ppm TWA; 25 mg/m³ TWA

ACGIH

Acetic Acid
CAS # 64-19-7 15 ppm STEL; 10 ppm TWA

ACGIH Biological Exposure Indices (BEI)

None established

Section 8 Exposure controls and personal protection (Continued)

DFG MAK

Acetic Acid
CAS # 64-19-7 20 ppm Peak; 50 mg/m³ Peak; 10 ppm TWA MAK; 25 mg/m³ TWA MAK

Ireland

Acetic Acid
CAS # 64-19-7 20 ppm TWA; 50 mg/m³ TWA; 20 ppm STEL; 50 mg/m³ STEL

IOELVs

Acetic Acid
CAS # 64-19-7 25 mg/m³ TWA; 10 ppm TWA; 50 mg/m³ STEL; 20 ppm STEL

NIOSH

Acetic Acid
CAS # 64-19-7 50 ppm IDLH; 15 ppm STEL; 37 mg/m³ STEL; 10 ppm TWA; 25 mg/m³ TWA

China

Acetic Acid
CAS # 64-19-7 20 mg/m³ STEL; 10 mg/m³ TWA

Croatia

Acetic Acid
CAS # 64-19-7 10 ppm TWA [GVI]; 25 mg/m³ TWA [GVI]; 20 ppm STEL [KGI]; 50 mg/m³ STEL [KGI]

Japan

Acetic Acid
CAS # 64-19-7 10 ppm OEL; 25 mg/m³ OEL

Sweden (AFS 2015:7 and amendments)

Acetic Acid
CAS # 64-19-7 5 ppm TLV NGV; 13 mg/m³ TLV NGV; 10 ppm Binding STEL Bindande KGV; 25 mg/m³ Binding STEL Bindande KGV

Turkey

Acetic Acid
CAS # 64-19-7 10 ppm TWA; 25 mg/m³ TWA

8.2 Exposure controls

Engineering controls

No special engineering controls are required. Use with good general ventilation.

Eye protection

Safety glasses or chemical goggles should be worn to prevent eye contact.

Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.

Skin protection

Wear impervious gloves such as Nitrile or equivalent and protective clothing. Refer to U.S. OSHA 29 CFR 1910.138, European Standard EN 374, EN 14605:2005+A1:2009 or appropriate government standards.

Respiratory protection

Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Liquid	Density and/or relative density	Not determined
Color	Colorless	Solubility	
Odor	Pungent	Water	Miscible
pH	<2	Organic	Not determined
Freezing point	Not determined	Partition coefficient n-octanol/water (log value)	Not determined
Boiling point or initial boiling point and boiling range	105.47°C (221.846°F)	Auto-ignition temp.	Not applicable
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	2.13 kPa@20°C
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not applicable		
Relative vapor density	2.07 (air=1)		
Particle characteristics	Not applicable		

9.2 Other information

Information with regard to physical hazard classes

No further relevant information available.

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

10.1 Reactivity	No further relevant information available.
10.2 Chemical stability	The product is stable in accordance with recommended storage conditions.
10.3 Possibility of hazardous reactions	No further relevant information available.
10.4 Conditions to avoid	To maintain product performance keep away from strong acids, strong bases, strong oxidizers. Avoid exposure to heat and direct sunlight.
10.5 Incompatible materials	No further relevant information available.

Section 10 Stability and reactivity (Continued)

10.6 Hazardous decomposition products

No decomposition products posing significant hazards would be expected from this product (an aqueous solution).

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

Acetic Acid CAS # 64-19-7	Dermal LD50 Rabbit 1060 mg/kg (JAPAN_GHS); Inhalation LC50 Rat 11.4 mg/L 4 h (NLM_CIP); Oral LD50 Rat 3310 mg/kg (JAPAN_GHS)
Sulphanilic Acid CAS # 121-57-3	Dermal LD50 Rat >2000 mg/kg (ECHA_API); Oral LD50 Rat 12300 mg/kg (NLM_CIP)

Primary routes of exposure Eye contact, ingestion, inhalation, and skin contact.

Acute toxicity Not classified based on available data.

Skin corrosion/irritation Contact may cause severe skin burns.

Serious eye damage/irritation Contact may cause serious eye damage.

Respiratory or skin sensitisation May cause an allergic skin reaction.

Germ cell mutagenicity Not classified based on available data.

Carcinogenicity No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.

Reproductive toxicity Not classified based on available data.

Specific target organ toxicity (STOT) – single exposure

Not classified based on available data.

Specific target organ toxicity (STOT) – repeated exposure

Not classified based on available data.

Aspiration hazard Not classified based on available data.

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information No further relevant information available.

SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 12 Ecological information

12.1 Toxicity

Fresh water species

Acetic Acid
CAS # 64-19-7

LC50 96 h Pimephales promelas: 79 mg/L [static] (EPA); LC50 96 h Lepomis macrochirus: 75 mg/L [static] (EPA)

Sulphanilic Acid
CAS # 121-57-3

LC50 96 h Pimephales promelas: 77.8 - 129.6 mg/L [static] (EPA)

Microtox/organisms

No information available.

Water flea

Acetic Acid
CAS # 64-19-7

EC50 48 h Daphnia magna: 65 mg/L [Static] (EPA)

Sulphanilic Acid
CAS # 121-57-3

EC50 48 h Daphnia magna: 85.66 mg/L (IUCLID)

Fresh water algae

Sulphanilic Acid
CAS # 121-57-3

EC50 72 h Desmodesmus subspicatus: 91 mg/L (IUCLID)

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulative potential Not determined for the product.

12.4 Mobility in soil Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects No further relevant information available.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Package disposal

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Additional information

Suggested European waste catalogue 18 01 06* - chemicals consisting of or containing dangerous substances. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport information

Shipping information	IATA	IMDG	US DOT	European ADR	Canadian TDG
14.1 UN/ID number	2790	2790	2790	2790	PIN - 2790

SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 14 Transport information (Continued)

14.2 UN proper shipping name	Acetic acid solution				
14.3 Transport hazard class(es)	8 Corrosives	8 Corrosive	8 Corrosives	8 Corrosive	8 Corrosives
Subsidiary risk	None	None	None	None	None
Classification code	Not applicable	Not applicable	Not applicable	C3	Not applicable
14.4 Packing group	III	III	III	III	III
Special provisions	A803	None	148	597, 647	None
Additional information					
IATA ERG code	8L	Not applicable	Not applicable	Not applicable	Not applicable
EmS	Not applicable	F-A, S-B	Not applicable	Not applicable	Not applicable
NAERG code	Not applicable	Not applicable	153	Not applicable	153
14.5 Environmental hazards					
Marine pollutant	Not applicable	No	Not applicable	Not applicable	Not applicable
14.6 Special precautions for user	Warning: Corrosive substances.				
14.7 Maritime transport in bulk according to IMO instruments	Not applicable				

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

CAS # 64-19-7 Acetic Acid

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Section 15 Regulatory information (Continued)

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

CAS # 64-19-7 Acetic Acid

New Jersey Dept. of Health Right To Know (RTK) List

CAS # 64-19-7 Acetic Acid

Pennsylvania Right To Know (RTK) List

CAS # 64-19-7 Acetic Acid

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany)

WGK 1, low water endangering

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

CAS # 64-19-7 Acetic Acid Entry No.: 75 (B)

CAS # 121-57-3 Sulphanilic Acid Entry No.: 75

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product is exempt from WHMIS label and SDS requirements.

China

Section 15 Regulatory information (Continued)

Catalog of Hazardous Chemicals - Hazardous Chemicals

CAS # 64-19-7 Acetic Acid
CAS # 121-57-3 Sulphanilic Acid

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Türkiye-REACH - KKDİK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 0 Health: 3 Reactivity with water: 0 Physical contact: 3	Code 0=None 1=Slight 2=Caution 3=Severe
--------------------------------------	---	---

Revision changes Update supplier addresses in Section 1.3

Document version and issue/revision date

Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20
Document ID: B71520
Version: AH

Hazard Classification Procedure

This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.

Description of hazard class and hazard statements from Section 3

Eye Irrit. 2 - Eye Irritation Category 2
Flam. Liq. 3 - Flammable Liquids, Category 3
Skin Corr. 1A - Skin Corrosion Category 1A
Skin Irrit. 2 - Skin Irritation Category 2
Skin Sens. 1 - Skin Sensitization Category 1
H226 - Flammable liquid and vapour.
H314 - Causes severe skin burns and eye damage.
H315 - Causes skin irritation.

Section 16 Other information (Continued)

Abbreviations and acronyms

H317 - May cause an allergic skin reaction.
H319 - Causes serious eye irritation.
ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
US OSHA - United States Occupational Safety and Health Administration
vPvB - very Persistent and very Bioaccumulative substances
WHMIS - Workplace Hazardous Material Information System

SAFETY DATA SHEET

Document ID: B71520 Version AH
Revision Date (year/month/day) 2025/07/13
Last Revision Date (year/month/day) 2024/09/20

Section 16 Other information (Continued)

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.