

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name	MicroScan Conventional Panels Negative
Part number	B1017-406, B1017-407, B1017-408, B1017-409, B1017-411, B1017-412, B1017-414, B1017-419, B1017-420, B1017-421, B1017-422, B1017-423, B1017-424, B1017-425, B1017-426, B1017-427, B1017-428, B1017-429, B1017-430, B1017-431, B1017-432, B1017-433, B1017-434, B1017-435, B1017-436, B1017-437, B1017-438, B1020-102B, C42464, C48639, C48640, C89643, C89644, C89645
Synonym	Panel Negative Combo 50 Panel Negative Urine Combo 51 Panel Negative BP Combo 41 Panel Negative Urine Combo 55 Panel Negative BP Combo 44 Panel Negative MIC 38 Panel Negative Urine Combo 61 Panel Negative MIC 43 Panel Negative Combo 67 Panel Negative Combo 68 Panel Negative/Urine Combo 73 Negative MIC Panel Type 45 Negative MIC Panel Type 46 Negative Urine Combo Panel Type 74 Negative Urine Combo Panel Type 75 Negative Urine Combo Panel Type 76 Negative Urine Combo Panel Type 77 Negative Urine Combo Panel Type 78 Negative Urine Combo Panel Type 79 PNL, DETECT NEG MIC 2 (RUO) Panel, Neg MIC 53 Panel, Neg MIC 54 Panel, Neg Urine Combo 84 Panel, Neg Urine Combo 85 Panel, Neg Urine Combo 86 Panel, Neg Urine Combo 87 Panel, Neg Urine Combo 88 PNL, NEG MIC 56 PNL, NEG URINE COMBO 89 PNL, NEG URINE COMBO 90 Neg MIC Panel Type 42 PNL, NEG URINE COMBO 101 PNL, NEG URINE COMBO 102 PNL, NEG URINE COMBO 103

1.2 Relevant identified uses of the substance or mixture and uses advised against

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

Product use For In Vitro Diagnostic Use. See product literature for details.
This safety data sheet is provided for information only, not required by regulations

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter, Inc.
250 S. Kraemer Blvd
Brea, CA 92821, U.S.A.
Tel: 800-854-3633

Supplier

CANADA
Beckman Coulter Canada LP
7075 Financial Drive
Mississauga, ON L5N 6V8
Canada
1-800-463-7828

UNITED KINGDOM
Beckman Coulter (UK) Ltd.
Amersham Place
Little Chalfont
Buckinghamshire
United Kingdom, HP7 9NA
01494 441181

AUSTRALIA
Beckman Coulter Australia Pty Ltd
23-27 Chaplin Drive
Lane Cove NSW 2066
Australia
ABN 81 002 011 672
24 Hour emergency contact phone
number:
1800 060 881

SWITZERLAND
Beckman Coulter Eurocenter SA
22, rue Juste-Olivier, Case Postale
1044,
CH-1260 Nyon 1, Switzerland.
Telephone: +41 (0)22 365 36 11
Monday through Friday, 9:00 am to
7:00pm)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

Beckman Coulter Ireland Inc.
Lismeehan
O'Callaghan's Mills
Co. Clare
Ireland
Tel: 353 (0)65 6831100

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00

e-mail address

SDSNT@beckman.com

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001)
703-527-3887

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description

Panels
Colorless; Solid; Broth odor

Classification according to EC 1272/2008 (CLP/GHS)

Not classified as hazardous per EC 1272/2008 (CLP/GHS)

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS

Not classified as hazardous per US-OSHA HCS 2012 and UN GHS

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS

Not classified as hazardous per EC 1272/2008 (CLP/GHS), US-OSHA and GHS

2.3 Other hazards

This product contains material(s) of animal origin. Observe general safety guidelines for protection when handling this product.

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:

None

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 4 First aid measures

4.1 Description of first aid measures

Inhalation	If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.
Eye contact	If product enters eyes, rinse eyes gently with water as a precaution.
Skin contact	In case of skin contact, rinse with water as a precaution.
Ingestion	If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

4.2 Most important symptoms and effects, both acute and delayed

No adverse symptoms or effects have been identified.

4.3 Indication of any immediate medical attention and special treatment needed

No specific medical attention or treatment required.

Section 5 Firefighting measures

5.1 Extinguishing media In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam. For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

No special hazards determined.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product.

5.3 Advice for firefighters

Protective equipment Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions This product contains a material of animal origin. Observe general safety guidelines for protection during clean up procedures.
Wear protective gloves, protective clothing and eye/face protection.

6.2 Environmental precautions Avoid release to the environment.

6.3 Methods and material for containment and cleaning up

Spill and leak procedures Sweep up spill and place in container suitable for disposal.
Dispose of all waste material in accordance with local guidelines.

6.4 Reference to other sections Refer sections 8 and 13.

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 7 Handling and storage

- 7.1 Precautions for safe handling** This product should be handled as though capable of transmitting infectious diseases. Universal precautions should be followed when using this product.
- 7.2 Conditions for safe storage, including any incompatibilities**
Store at 2 to 25°C , as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).
- 7.3 Specific end uses** No further relevant information available.

Section 8 Exposure controls and personal protection

- 8.1 Control parameters**
- Exposure limits**
- | | |
|--|------------------|
| US OSHA | None established |
| ACGIH | None established |
| ACGIH Biological Exposure Indices (BEI) | None established |
| DFG MAK | None established |
| Ireland | None established |
| IOELVs | None established |
| NIOSH | None established |
| China | None established |
| Croatia | None established |
| Japan | None established |
| Sweden (AFS 2015:7 and amendments) | None established |
| Turkey | None established |
- 8.2 Exposure controls**
- Engineering controls** No special engineering controls are required. Use with good general ventilation.
- Eye protection** Safety glasses or chemical goggles should be worn to prevent eye contact.
Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.
- Skin protection** Wear protective clothing and impervious gloves, as appropriate.
- Respiratory protection** Under normal conditions, the use of this product should not require respiratory protection.

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Solid	Density and/or relative density	Not determined
Color	Colorless	Solubility	
Odor	Broth odor	Water	Insoluble in water.
pH	Not applicable	Organic	Not determined
Melting Point	Not applicable	Partition coefficient n-octanol/water (log value)	Not applicable
Boiling point or initial boiling point and boiling range	Not applicable	Auto-ignition temp.	Not applicable
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	Not applicable
		Kinematic viscosity	Not applicable
Lower and upper explosion limit	Not applicable		
Relative vapor density	Not applicable		
Particle characteristics	Not determined		

9.2 Other information

Information with regard to physical hazard classes

Not applicable
No further relevant information available.

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

10.1 Reactivity	No further relevant information available.
10.2 Chemical stability	The product is stable in accordance with recommended storage conditions.
10.3 Possibility of hazardous reactions	No further relevant information available.
10.4 Conditions to avoid	No further relevant information available.
10.5 Incompatible materials	No further relevant information available.
10.6 Hazardous decomposition products	When stored as labeled, no known hazardous decomposition products are formed during the shelf-life of this product.

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

Not applicable

Primary routes of exposure Common routes of entry include inhalation, ingestion and eye/skin contact. Specific paths of concern for potentially infectious materials are skin puncture, contact with broken skin, contact with mucous membranes and inhalation of aerosolized material.

Acute toxicity Not classified based on available data.

Skin corrosion/irritation Not classified based on available data.

Serious eye damage/irritation Not classified based on available data.

Respiratory or skin sensitisation Not classified based on available data.

Germ cell mutagenicity Not classified based on available data.

Carcinogenicity No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.

Reproductive toxicity Not classified based on available data.

Specific target organ toxicity (STOT) – single exposure

Not classified based on available data.

Specific target organ toxicity (STOT) – repeated exposure

Not classified based on available data.

Aspiration hazard Not classified based on available data.

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information This product contains material(s) of animal origin. Observe general safety guidelines for protection when handling this product.

Section 12 Ecological information

12.1 Toxicity

Fresh water species No information available.

Microtox/organisms No information available.

Water flea No information available.

Fresh water algae No information available.

12.2 Persistence and degradability Not determined for the product.

12.3 Bioaccumulative potential Not determined for the product.

Section 12 Ecological information (Continued)

- 12.4 Mobility in soil** Not determined for the product.
- 12.5 Results of PBT and vPvB assessment**
Not determined for the product. PBT: Not applicable, vPvB: Not applicable.
- 12.6 Endocrine disrupting properties**
This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).
- 12.7 Other adverse effects** No further relevant information available.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Dispose of as potentially biohazardous waste and in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or approved waste-disposal company for information.

Package disposal

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Additional information

No further relevant information available.

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

- 14.1 UN/ID number:** Not regulated for transportation
- 14.2 UN proper shipping name:** Not regulated for transportation
- 14.3 Transport hazard class(es):** Not regulated for transportation
- 14.4 Packing group:** Not regulated for transportation
- 14.5 Environmental hazards:** Not regulated for transportation
- 14.6 Special precautions for user:** None
- 14.7 Maritime transport in bulk according to IMO instruments:** Not applicable

Section 15 Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
US Federal and State Regulations

Section 15 Regulatory information (Continued)

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

No ingredients listed.

New Jersey Dept. of Health Right To Know (RTK) List

No ingredients listed.

Pennsylvania Right To Know (RTK) List

No ingredients listed.

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

No ingredients listed.

Section 15 Regulatory information (Continued)

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product is exempt from WHMIS label and SDS requirements.

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

No ingredients listed.

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Türkiye-REACH - KKDIK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 0 Health: 1 Reactivity with water: 0 Physical contact: 1	Code 0=None 1=Slight 2=Caution 3=Severe
--------------------------------------	---	---

Revision changes Update supplier addresses in Section 1.3

Document version and issue/revision date

Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15
Document ID: B72544
Version: AU

Section 16 Other information (Continued)

Hazard Classification Procedure

This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
US OSHA - United States Occupational Safety and Health Administration
vPvB - very Persistent and very Bioaccumulative substances
WHMIS - Workplace Hazardous Material Information System

SAFETY DATA SHEET

Document ID: B72544 Version AU
Revision Date (year/month/day) 2025/07/18
Last Revision Date (year/month/day) 2024/08/15

Section 16 Other information (Continued)

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.