

**SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**

**1.1. Product identifier**

**Product Description:** Infinity Ammonia Reagent  
**Cat No. :** OSR61154

**1.2. Relevant identified uses of the substance or mixture and uses advised against**

**Recommended Use** In vitro diagnostic.  
**Uses advised against** No Information available

**1.3. Details of the supplier of the safety data sheet**

**Company** Fisher Diagnostics  
A Division of Fisher Scientific Company, LLC  
A Part of Thermo Fisher Scientific, Inc.  
8365 Valley Pike  
Middletown, VA 22645-1905  
Tel: (800) 528-0494

**E-mail address** begel.sdsdesk@thermofisher.com

**1.4. Emergency telephone number**

Chemtrec US: (800) 424-9300  
Chemtrec Outside US: +1 (703) 741-5970

**SECTION 2: HAZARDS IDENTIFICATION**

**2.1. Classification of the substance or mixture**

**CLP Classification - Regulation (EC) No 1272/2008**

**Physical hazards**

Based on available data, the classification criteria are not met

**Health hazards**

Based on available data, the classification criteria are not met

**Environmental hazards**

Based on available data, the classification criteria are not met

*Full text of Hazard Statements: see section 16*

**2.2. Label elements**

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

None required

EUH210 - Safety data sheet available on request

## 2.3. Other hazards

No information available

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

### 3.2. Mixtures

Component	CAS No	EC No	Weight %	CLP Classification - Regulation (EC) No 1272/2008
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-	77-86-1	201-064-4	0.1<1.0	-
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	EEC No. 206-167-8	0.1<1.0	-
Sodium azide	26628-22-8	247-852-1	<0.1	Acute Tox. 2 (H300) Aquatic Acute 1 (H400) Aquatic Chronic 1 (H410) (EUH032)
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	EEC No. 214-684-5	0.1<1.0	-
Water	7732-18-5	231-791-2	90<100	-
Adenosine 5-(trihydrogen diphosphate), P->5-ester with 1,4-dihydro-1-.beta.-D-ribofuranosyl-3-pyridinecarboxamide, disodium salt	606-68-8	EEC No. 210-123-3	<0.1	-

Component	Specific concentration limits (SCL's)	M-Factor	Component notes
Sodium azide	-	1	-

Full text of Hazard Statements: see section 16

## SECTION 4: FIRST AID MEASURES

### 4.1. Description of first aid measures

<b>General Advice</b>	Show this safety data sheet to the doctor in attendance.
<b>Eye Contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
<b>Skin Contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
<b>Ingestion</b>	Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. Get medical attention immediately if symptoms occur.

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

**Self-Protection of the First Aider** No special precautions required.

## **4.2. Most important symptoms and effects, both acute and delayed**

None reasonably foreseeable.

## **4.3. Indication of any immediate medical attention and special treatment needed**

**Notes to Physician** Treat symptomatically.

## **SECTION 5: FIREFIGHTING MEASURES**

### **5.1. Extinguishing media**

#### **Suitable Extinguishing Media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### **Extinguishing media which must not be used for safety reasons**

No information available.

### **5.2. Special hazards arising from the substance or mixture**

Thermal decomposition can lead to release of irritating gases and vapors.

#### **Hazardous Combustion Products**

Carbon oxides.

### **5.3. Advice for firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## **SECTION 6: ACCIDENTAL RELEASE MEASURES**

### **6.1. Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation. Use personal protective equipment as required.

### **6.2. Environmental precautions**

Should not be released into the environment.

### **6.3. Methods and material for containment and cleaning up**

Sweep up and shovel into suitable containers for disposal.  
Prevent further leakage or spillage if safe to do so

### **6.4. Reference to other sections**

Refer to protective measures listed in Sections 8 and 13.

## **SECTION 7: HANDLING AND STORAGE**

### **7.1. Precautions for safe handling**

Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid contact with skin, eyes or clothing. Avoid

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

ingestion and inhalation.

## Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Wash hands before breaks and after work.

## 7.2. Conditions for safe storage, including any incompatibilities

Keep container tightly closed.

**Technical Rules for Hazardous Substances (TRGS) 510**  
**Storage Class (LGK) (Germany)**

Storage Class/LGK 12

**Switzerland - Storage of hazardous substances**

Storage class - SC 10/12

<https://www.kvu.ch/de/themen/stoffe-und-produkte>

<https://www.kvu.ch/fr/themes/substances-et-produits>

<https://www.kvu.ch/it/temi/sostanze-e-prodotti>

## 7.3. Specific end use(s)

Use in laboratories

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1. Control parameters

#### Exposure limits

List source(s): **CH** - The Government of Switzerland has set a directive on limit values for working materials (Grenzwerte am arbeitsplatz) which is based on the Swiss Federal Regulation "Verordnung über die Verhütung von Unfällen und Berufskrankheiten". This directive is administered, periodically revised and enforced by SUVA (Swiss National Accident Insurance Fund). **EU** - Commission Directive (EU) 2019/1831 of 24 October 2019 establishing a fifth list of indicative occupational exposure limit values pursuant to Council Directive 98/24/EC and amending Commission Directive 2000/39/EC **UK** - EH40/2005 Work Exposure Limits, Forth edition. Published 2020. **IRE** - 2010 Code of Practice for the Safety, Health and Welfare at Work (Chemical Agents) Regulations 2001. Published by the Health and Safety Authority.

Component	European Union	The United Kingdom	France	Belgium	Spain
Sodium azide	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	TWA / VME: 0.1 mg/m <sup>3</sup> (8 heures). restrictive limit STEL / VLCT: 0.3 mg/m <sup>3</sup> . restrictive limit Peau	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	STEL / VLA-EC: 0.3 mg/m <sup>3</sup> (15 minutos). TWA / VLA-ED: 0.1 mg/m <sup>3</sup> (8 horas) Piel

Component	Italy	Germany	Portugal	The Netherlands	Finland
Sodium azide	TWA: 0.1 mg/m <sup>3</sup> 8 ore. Media Ponderata nel Tempo STEL: 0.3 mg/m <sup>3</sup> 15 minuti. Breve termine Pelle	MAK 0.2 mg/m <sup>3</sup> (inhalable)	STEL: 0.3 mg/m <sup>3</sup> 15 minutos Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm TWA: 0.1 mg/m <sup>3</sup> 8 horas Pele	huid STEL: 0.3 mg/m <sup>3</sup> 15 minuten TWA: 0.1 mg/m <sup>3</sup> 8 uren	TWA: 0.1 mg/m <sup>3</sup> 8 tunteina STEL: 0.3 mg/m <sup>3</sup> 15 minuutteina Iho

Component	Austria	Denmark	Switzerland	Poland	Norway
Sodium azide	Haut MAK-KZGW: 0.3 mg/m <sup>3</sup> 15 Minuten MAK-TMW: 0.1 mg/m <sup>3</sup> 8 Stunden	TWA: 0.1 mg/m <sup>3</sup> 8 timer Hud	STEL: 0.4 mg/m <sup>3</sup> 15 Minuten TWA: 0.2 mg/m <sup>3</sup> 8 Stunden	STEL: 0.3 mg/m <sup>3</sup> 15 minutach TWA: 0.1 mg/m <sup>3</sup> 8 godzinach	TWA: 0.1 mg/m <sup>3</sup> 8 timer STEL: 0.3 mg/m <sup>3</sup> 15 minutter. value from the regulation

Component	Bulgaria	Croatia	Ireland	Cyprus	Czech Republic
Sodium azide	TWA: 0.1 mg/m <sup>3</sup> STEL : 0.3 mg/m <sup>3</sup> Skin notation	kože TWA-GVI: 0.1 mg/m <sup>3</sup> 8 satima.	TWA: 0.1 mg/m <sup>3</sup> 8 hr. STEL: 0.3 mg/m <sup>3</sup> 15 min Skin	Skin-potential for cutaneous absorption STEL: 0.3 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> 8 hodinách. Potential for cutaneous

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

		STEL-KGVI: 0.3 mg/m <sup>3</sup> 15 minutama.		TWA: 0.1 mg/m <sup>3</sup>	absorption Ceiling: 0.3 mg/m <sup>3</sup>
--	--	--	--	----------------------------	--

Component	Estonia	Gibraltar	Greece	Hungary	Iceland
Sodium azide	Nahk TWA: 0.1 mg/m <sup>3</sup> 8 tundides. STEL: 0.3 mg/m <sup>3</sup> 15 minutites.	Skin notation TWA: 0.1 mg/m <sup>3</sup> 8 hr STEL: 0.3 mg/m <sup>3</sup> 15 min	STEL: 0.1 ppm STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 ppm TWA: 0.3 mg/m <sup>3</sup>	STEL: 0.3 mg/m <sup>3</sup> 15 percekben. CK TWA: 0.1 mg/m <sup>3</sup> 8 óraban. AK	STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup> 8 klukkustundum. Skin notation

Component	Latvia	Lithuania	Luxembourg	Malta	Romania
Sodium azide	skin - potential for cutaneous exposure STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> IPRD Oda STEL: 0.3 mg/m <sup>3</sup>	Possibility of significant uptake through the skin TWA: 0.1 mg/m <sup>3</sup> 8 Stunden STEL: 0.3 mg/m <sup>3</sup> 15 Minuten	possibility of significant uptake through the skin TWA: 0.1 mg/m <sup>3</sup> STEL: 0.3 mg/m <sup>3</sup> 15 minuti	Skin notation TWA: 0.1 mg/m <sup>3</sup> 8 ore STEL: 0.3 mg/m <sup>3</sup> 15 minute

Component	Russia	Slovak Republic	Slovenia	Sweden	Turkey
Sodium azide		Ceiling: 0.3 mg/m <sup>3</sup> Potential for cutaneous absorption TWA: 0.1 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> 8 urah Koža STEL: 0.3 mg/m <sup>3</sup> 15 minutah	Binding STEL: 0.3 mg/m <sup>3</sup> 15 minuter TLV: 0.1 mg/m <sup>3</sup> 8 timmar. NGV	Deri TWA: 0.1 mg/m <sup>3</sup> 8 saat STEL: 0.3 mg/m <sup>3</sup> 15 dakika

## Biological limit values

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

## Monitoring methods

BS EN 14042:2003 Title Identifier: Workplace atmospheres. Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents.

## Derived No Effect Level (DNEL) / Derived Minimum Effect Level (DMEL)

See table for values

Component	Acute effects local (Dermal)	Acute effects systemic (Dermal)	Chronic effects local (Dermal)	Chronic effects systemic (Dermal)
1,3-Propanediol, 2-amino-2-(hydroxymethyl)- 77-86-1 ( 0.1<1.0 )				DNEL = 166.7mg/kg bw/day
Sodium azide 26628-22-8 ( <0.1 )				DNEL = 46.7µg/kg bw/day
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride 1185-53-1 ( 0.1<1.0 )				DNEL = 216.6mg/kg bw/day

Component	Acute effects local (Inhalation)	Acute effects systemic (Inhalation)	Chronic effects local (Inhalation)	Chronic effects systemic (Inhalation)
1,3-Propanediol, 2-amino-2-(hydroxymethyl)- 77-86-1 ( 0.1<1.0 )				DNEL = 117.5mg/m <sup>3</sup>
Sodium azide 26628-22-8 ( <0.1 )				DNEL = 0.164mg/m <sup>3</sup>
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride				DNEL = 152.8mg/m <sup>3</sup>

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

1185-53-1 ( 0.1<1.0 )				
-----------------------	--	--	--	--

## Predicted No Effect Concentration (PNEC)

See values below.

Component	Fresh water	Fresh water sediment	Water Intermittent	Microorganisms in sewage treatment	Soil (Agriculture)
1,3-Propanediol, 2-amino-2-(hydroxymethyl) - 77-86-1 ( 0.1<1.0 )				PNEC = 300mg/L	
Sodium azide 26628-22-8 ( <0.1 )	PNEC = 0.35µg/L	PNEC = 16.7µg/kg sediment dw	PNEC = 3.5µg/L	PNEC = 30µg/L	

Component	Marine water	Marine water sediment	Marine water Intermittent	Food chain	Air
Sodium azide 26628-22-8 ( <0.1 )	PNEC = 15ng/L	PNEC = 0.72µg/kg sediment dw	PNEC = 150ng/L		

## 8.2. Exposure controls

### Engineering Measures

None under normal use conditions.

### Personal protective equipment

#### Eye Protection

Wear safety glasses with side shields (or goggles) (European standard - EN 166)

#### Hand Protection

Protective gloves

Glove material	Breakthrough time	Glove thickness	EU standard	Glove comments
Natural rubber	See manufacturers	-	EN 374	(minimum requirement)
Nitrile rubber	recommendations			
Neoprene				
PVC				

#### Skin and body protection

Long sleeved clothing.

Inspect gloves before use. observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information) gloves are suitable for the task: Chemical compatability, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion. gloves with care avoiding skin contamination.

#### Respiratory Protection

No protective equipment is needed under normal use conditions.

#### Large scale/emergency use

Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced  
**Recommended Filter type:** Particle filter

#### Small scale/Laboratory use

Maintain adequate ventilation

#### Environmental exposure controls

No information available.

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

### 9.1. Information on basic physical and chemical properties

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

<b>Physical State</b>	Liquid	
<b>Appearance</b>	Clear, colorless	
<b>Odor</b>	Mild	
<b>Odor Threshold</b>	No data available	
<b>Melting Point/Range</b>	No data available	
<b>Softening Point</b>	No data available	
<b>Boiling Point/Range</b>	100 °C / 212 °F	
<b>Flammability (liquid)</b>	No data available	
<b>Flammability (solid,gas)</b>	Not applicable	Liquid
<b>Explosion Limits</b>	No data available	
<b>Flash Point</b>	No data available	<b>Method -</b> No information available
<b>Autoignition Temperature</b>	No data available	
<b>Decomposition Temperature</b>	No data available	
<b>pH</b>	8.7	
<b>Viscosity</b>	No data available	
<b>Water Solubility</b>	Soluble in water	
<b>Solubility in other solvents</b>	No information available	
<b>Partition Coefficient (n-octanol/water)</b>		
<b>Component</b>	<b>log Pow</b>	
1,3-Propanediol,	-3.6	
2-amino-2-(hydroxymethyl)-, hydrochloride		
<b>Vapor Pressure</b>	No data available	
<b>Density / Specific Gravity</b>	1.006 kg/L @ 20C	
<b>Bulk Density</b>	Not applicable	Liquid
<b>Vapor Density</b>	No data available	(Air = 1.0)
<b>Particle characteristics</b>	Not applicable (liquid)	

## 9.2. Other information

**Explosive Properties** Not explosive

## SECTION 10: STABILITY AND REACTIVITY

### 10.1. Reactivity

None known, based on information available

### 10.2. Chemical stability

Stable under recommended storage conditions.

### 10.3. Possibility of hazardous reactions

**Hazardous Polymerization** Hazardous polymerization does not occur.  
**Hazardous Reactions** None under normal processing.

### 10.4. Conditions to avoid

Incompatible products. Excess heat.

### 10.5. Incompatible materials

None known.

### 10.6. Hazardous decomposition products

Carbon oxides.

## SECTION 11: TOXICOLOGICAL INFORMATION

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

## 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

### Product Information

#### (a) acute toxicity;

Oral

Based on available data, the classification criteria are not met

Dermal

Based on available data, the classification criteria are not met

Inhalation

Based on available data, the classification criteria are not met

### Toxicology data for the components

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
1,3-Propanediol, 2-amino-2-(hydroxymethyl)- Sodium azide	LD50 = 5900 mg/kg ( Rat )	LD50 > 5000 mg/kg ( Rat )	-
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	OECD 425 (Rat) LD50 > 5000 mg/kg bw	OECD 402 (Rat) LD50 > 5000 mg/kg bw	LC50 0.054 - 0.52 mg/L ( Rat ) 4 h
Water	-	-	-

(b) skin corrosion/irritation; No data available

(c) serious eye damage/irritation; No data available

#### (d) respiratory or skin sensitization;

Respiratory

No data available

Skin

No data available

Component	Test method	Test species	Study result
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride 1185-53-1 ( 0.1<1.0 )	OECD Test Guideline 406	guinea pig	non-sensitising

(e) germ cell mutagenicity; No data available

Component	Test method	Test species	Study result
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride 1185-53-1 ( 0.1<1.0 )	OECD Test Guideline 471 Bacterial Reverse Mutation Test	Mammalian in vitro	negative

(f) carcinogenicity; No data available

There are no known carcinogenic chemicals in this product

(g) reproductive toxicity; No data available

(h) STOT-single exposure; No data available

(i) STOT-repeated exposure; No data available

Target Organs

None known.

(j) aspiration hazard; No data available



# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

**Symptoms / effects, both acute and delayed** No information available.

## 11.2. Information on other hazards

**Endocrine Disrupting Properties** Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

## SECTION 12: ECOLOGICAL INFORMATION

### 12.1. Toxicity

#### Ecotoxicity effects

The environmental impact of this product has not been fully investigated.

Component	Freshwater Fish	Water Flea	Freshwater Algae
Sodium azide	LC50: = 0.7 mg/L, 96h (Lepomis macrochirus) LC50: = 0.8 mg/L, 96h (Oncorhynchus mykiss) LC50: = 5.46 mg/L, 96h flow-through (Pimephales promelas)		
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride		Daphnia Magna EC50 >100 mg/L (48h)	

Component	Microtox	M-Factor
Sodium azide		1
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	OECD 209 EC50 > 1000 mg/L (3h)	

### 12.2. Persistence and degradability

#### Persistence

Soluble in water, Persistence is unlikely, based on information available.

### 12.3. Bioaccumulative potential

Bioaccumulation is unlikely

Component	log Pow	Bioconcentration factor (BCF)
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	-3.6	No data available

### 12.4. Mobility in soil

The product is water soluble, and may spread in water systems. Will likely be mobile in the environment due to its water solubility. Highly mobile in soils

### 12.5. Results of PBT and vPvB assessment

No data available for assessment.

### 12.6. Endocrine disrupting properties

#### Endocrine Disruptor Information

This product does not contain any known or suspected endocrine disruptors

### 12.7. Other adverse effects

#### Persistent Organic Pollutant Ozone Depletion Potential

This product does not contain any known or suspected substance

This product does not contain any known or suspected substance

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

## SECTION 13: DISPOSAL CONSIDERATIONS

### 13.1. Waste treatment methods

<b>Waste from Residues/Unused Products</b>	Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. Consult local, regional, and national hazardous waste regulations to ensure complete and accurate classification.
<b>Contaminated Packaging</b>	Dispose of contents/containers in accordance with local regulations.
<b>European Waste Catalogue (EWC)</b>	According to the European Waste Catalog, Waste Codes are not product specific, but application specific.
<b>Other Information</b>	Waste codes should be assigned by the user based on the application for which the product was used.
<b>Switzerland - Waste Ordinance</b>	Disposal should be in accordance with applicable regional, national and local laws and regulations. Ordinance on the Avoidance and the Disposal of Waste (Waste Ordinance, ADWO) SR 814.600 <a href="https://www.fedlex.admin.ch/eli/cc/2015/891/en">https://www.fedlex.admin.ch/eli/cc/2015/891/en</a>

## SECTION 14: TRANSPORT INFORMATION

**IMDG/IMO** Not regulated

**14.1. UN number**  
**14.2. UN proper shipping name**  
**14.3. Transport hazard class(es)**  
**14.4. Packing group**

N/A

**ADR** Not regulated

**14.1. UN number**  
**14.2. UN proper shipping name**  
**14.3. Transport hazard class(es)**  
**14.4. Packing group**

**IATA** Not regulated

**14.1. UN number**  
**14.2. UN proper shipping name**  
**14.3. Transport hazard class(es)**  
**14.4. Packing group**

N/A

**14.5. Environmental hazards** No hazards identified

**14.6. Special precautions for user** No special precautions required

**14.7. Maritime transport in bulk according to IMO instruments** Not applicable, packaged goods

## SECTION 15: REGULATORY INFORMATION

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

**International Inventories**

OSR61154 AMMONIA REAGENT

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

X = listed, U.S.A. (TSCA), Canada (DSL/NDL), Europe (EINECS/ELINCS/NLP), Australia (AICS), Korea (KECL), China (IECSC), Japan (ENCS), Philippines (PICCS), Taiwan (TCSI), Japan (ISHL), New Zealand (NZIoC), Japan (ISHL). US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

Component	CAS No	EINECS	ELINCS	NLP	IECSC	TCSI	KECL	ENCS	ISHL
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-	77-86-1	201-064-4	-	-	X	X	KE-01403	X	X
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	206-167-8	-	-	-	X	-	-	-
Sodium azide	26628-22-8	247-852-1	-	-	X	X	KE-31357	X	X
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	214-684-5	-	-	X	X	KE-34819	X	-
Water	7732-18-5	231-791-2	-	-	X	X	KE-35400	X	-
Adenosine 5-(trihydrogen diphosphate), P->5-ester with 1,4-dihydro-1-.beta.-D-ribofuranosyl-3-pyridinecarboxamide, disodium salt	606-68-8	210-123-3	-	-	-	X	-	-	-

Component	CAS No	TSCA	TSCA Inventory notification - Active-Inactive	DSL	NDL	AICS	NZIoC	PICCS
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-	77-86-1	X	ACTIVE	X	-	X	X	X
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	-	-	X	-	X	X	-
Sodium azide	26628-22-8	X	ACTIVE	X	-	X	X	X
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	X	ACTIVE	X	-	X	X	X
Water	7732-18-5	X	ACTIVE	X	-	X	X	X
Adenosine 5-(trihydrogen diphosphate), P->5-ester with 1,4-dihydro-1-.beta.-D-ribofuranosyl-3-pyridinecarboxamide, disodium salt	606-68-8	X	ACTIVE	X	-	X	X	-

Legend: X - Listed '-' - Not Listed

KECL - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

## Authorisation/Restrictions according to EU REACH

Component	CAS No	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-	77-86-1	Not applicable	Not applicable
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	Not applicable	Not applicable
Sodium azide	26628-22-8	Not applicable	Not applicable
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	Not applicable	Not applicable
Water	7732-18-5	Not applicable	Not applicable
Adenosine 5-(trihydrogen diphosphate), P->5-ester with 1,4-dihydro-1-.beta.-D-ribofurano syl-3-pyridinecarboxamide, disodium salt	606-68-8	Not applicable	Not applicable

Regulation (EC) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of dangerous chemicals

Not applicable

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work .

Take note of Directive 2000/39/EC establishing a first list of indicative occupational exposure limit values

## National Regulations

**UK** - Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment

**WGK Classification** Water endangering class = 3 (self classification)

Component	Germany - Water Classification (VwVwS)	Germany - TA-Luft Class
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-	WGK1	
Sodium azide	WGK2	
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	WGK1	

## Swiss Regulations

Article 4 para. 4 of the Ordinance on the protection of young people in the workplace (SR 822.115) and Article 1 lit. f of the EAER regulation on hazardous work and young people (SR 822.115.2).

Take note on Article 13 Maternity Ordinance (SR 822.111.52) with regards expectant and nursing mothers.

## 15.2. Chemical safety assessment

Chemical Safety Assessment/Reports (CSA/CSR) are not required for mixtures

## SECTION 16: OTHER INFORMATION

### Full text of H-Statements referred to under sections 2 and 3

H300 - Fatal if swallowed

H400 - Very toxic to aquatic life

H410 - Very toxic to aquatic life with long lasting effects

EUH032 - Contact with acids liberates very toxic gas

### Legend

**CAS** - Chemical Abstracts Service

**EINECS/ELINCS** - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**IECSC** - Chinese Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**ENCS** - Japanese Existing and New Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**NZIoC** - New Zealand Inventory of Chemicals

**WEL** - Workplace Exposure Limit

**ACGIH** - American Conference of Governmental Industrial Hygienists

**DNEL** - Derived No Effect Level

**RPE** - Respiratory Protective Equipment

**LC50** - Lethal Concentration 50%

**NOEC** - No Observed Effect Concentration

**PBT** - Persistent, Bioaccumulative, Toxic

**TWA** - Time Weighted Average

**IARC** - International Agency for Research on Cancer Predicted No Effect Concentration (PNEC)

**LD50** - Lethal Dose 50%

**EC50** - Effective Concentration 50%

**POW** - Partition coefficient Octanol:Water

**vPvB** - very Persistent, very Bioaccumulative

**ADR** - European Agreement Concerning the International Carriage of Dangerous Goods by Road

**IMO/IMDG** - International Maritime Organization/International Maritime Dangerous Goods Code

**ICAO/IATA** - International Civil Aviation Organization/International Air Transport Association

**MARPOL** - International Convention for the Prevention of Pollution from Ships

# SAFETY DATA SHEET

Infinity Ammonia Reagent

Revision Date 28-Jan-2022

**OECD** - Organisation for Economic Co-operation and Development

**BCF** - Bioconcentration factor

**ATE** - Acute Toxicity Estimate

**VOC** - (volatile organic compound)

## Key literature references and sources for data

<https://echa.europa.eu/information-on-chemicals>

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

## Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

**Physical hazards** On basis of test data

**Health Hazards** Calculation method

**Environmental hazards** Calculation method

## Training Advice

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

**Revision Date** 28-Jan-2022

**Revision Summary** SDS sections updated.

**This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.  
COMMISSION REGULATION (EU) 2020/878 amending Annex II to Regulation (EC) No  
1907/2006 .**

**For Switzerland - Compiled in accordance with the technical provisions referred to in Annex 2,  
Number 3, ChemO (SR 813.11 - Ordinance on Protection against Dangerous Substances and  
Preparations).**

## Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet**

**SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING**

**1.1. Product identifier**

**Product Description:** Infinity Ammonia Standard  
**Cat No. :** OSR61154

**1.2. Relevant identified uses of the substance or mixture and uses advised against**

**Recommended Use** In vitro diagnostic.  
**Uses advised against** No Information available

**1.3. Details of the supplier of the safety data sheet**

**Company** Fisher Diagnostics  
A Division of Fisher Scientific Company, LLC  
A Part of Thermo Fisher Scientific, Inc.  
8365 Valley Pike  
Middletown, VA 22645-1905  
Tel: (800) 528-0494

**E-mail address** begel.sdsdesk@thermofisher.com

**1.4. Emergency telephone number**

Chemtrec US: (800) 424-9300  
Chemtrec Outside US: +1 (703) 741-5970

**SECTION 2: HAZARDS IDENTIFICATION**

**2.1. Classification of the substance or mixture**

**CLP Classification - Regulation (EC) No 1272/2008**

**Physical hazards**

Based on available data, the classification criteria are not met

**Health hazards**

Based on available data, the classification criteria are not met

**Environmental hazards**

Based on available data, the classification criteria are not met

*Full text of Hazard Statements: see section 16*

**2.2. Label elements**

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

None required

EUH210 - Safety data sheet available on request

## 2.3. Other hazards

No information available

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

### 3.2. Mixtures

Component	CAS No	EC No	Weight %	CLP Classification - Regulation (EC) No 1272/2008
Sodium azide	26628-22-8	247-852-1	<0.1	Acute Tox. 2 (H300) Aquatic Acute 1 (H400) Aquatic Chronic 1 (H410) (EUH032)
Water	7732-18-5	231-791-2	90<100	-
Ammonium chloride	12125-02-9	235-186-4	<0.001	Acute Tox. 4 (H302) Eye Irrit. 2 (H319)

Component	Specific concentration limits (SCL's)	M-Factor	Component notes
Sodium azide	-	1	-

Full text of Hazard Statements: see section 16

## SECTION 4: FIRST AID MEASURES

### 4.1. Description of first aid measures

<b>General Advice</b>	Show this safety data sheet to the doctor in attendance.
<b>Eye Contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
<b>Skin Contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
<b>Ingestion</b>	Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. Get medical attention immediately if symptoms occur.
<b>Self-Protection of the First Aider</b>	No special precautions required.

### 4.2. Most important symptoms and effects, both acute and delayed

None reasonably foreseeable.

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

## 4.3. Indication of any immediate medical attention and special treatment needed

**Notes to Physician** Treat symptomatically.

## SECTION 5: FIREFIGHTING MEASURES

### 5.1. Extinguishing media

#### **Suitable Extinguishing Media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### **Extinguishing media which must not be used for safety reasons**

No information available.

### 5.2. Special hazards arising from the substance or mixture

None reasonably foreseeable.

#### **Hazardous Combustion Products**

Carbon oxides.

### 5.3. Advice for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

## SECTION 6: ACCIDENTAL RELEASE MEASURES

### 6.1. Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Use personal protective equipment as required.

### 6.2. Environmental precautions

Should not be released into the environment.

### 6.3. Methods and material for containment and cleaning up

Sweep up and shovel into suitable containers for disposal.  
Prevent further leakage or spillage if safe to do so

### 6.4. Reference to other sections

Refer to protective measures listed in Sections 8 and 13.

## SECTION 7: HANDLING AND STORAGE

### 7.1. Precautions for safe handling

Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid contact with skin, eyes or clothing. Avoid ingestion and inhalation.

#### **Hygiene Measures**

Handle in accordance with good industrial hygiene and safety practice. Keep away from food, drink and animal feeding stuffs. Do not eat, drink or smoke when using this product. Remove and wash contaminated clothing and gloves, including the inside, before re-use. Wash hands before breaks and after work.



# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

## 7.2. Conditions for safe storage, including any incompatibilities

Keep container tightly closed.

**Technical Rules for Hazardous Substances (TRGS) 510**  
**Storage Class (LGK) (Germany)**

Storage Class/LGK 12

**Switzerland - Storage of hazardous substances**

Storage class - SC 10/12  
<https://www.kvu.ch/de/themen/stoffe-und-produkte>  
<https://www.kvu.ch/fr/themes/substances-et-produits>  
<https://www.kvu.ch/it/temi/sostanze-e-prodotti>

## 7.3. Specific end use(s)

Use in laboratories

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1. Control parameters

#### Exposure limits

List source(s): **CH** - The Government of Switzerland has set a directive on limit values for working materials (Grenzwerte am arbeitsplatz) which is based on the Swiss Federal Regulation "Verordnung über die Verhütung von Unfällen und Berufskrankheiten". This directive is administered, periodically revised and enforced by SUVA (Swiss National Accident Insurance Fund). **EU** - Commission Directive (EU) 2019/1831 of 24 October 2019 establishing a fifth list of indicative occupational exposure limit values pursuant to Council Directive 98/24/EC and amending Commission Directive 2000/39/EC **UK** - EH40/2005 Work Exposure Limits, Forth edition. Published 2020. **IRE** - 2010 Code of Practice for the Safety, Health and Welfare at Work (Chemical Agents) Regulations 2001. Published by the Health and Safety Authority.

Component	European Union	The United Kingdom	France	Belgium	Spain
Sodium azide	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	TWA / VME: 0.1 mg/m <sup>3</sup> (8 heures). restrictive limit STEL / VLCT: 0.3 mg/m <sup>3</sup> . restrictive limit Peau	Skin TWA 0.1 mg/m <sup>3</sup> STEL 0.3 mg/m <sup>3</sup>	STEL / VLA-EC: 0.3 mg/m <sup>3</sup> (15 minutos). TWA / VLA-ED: 0.1 mg/m <sup>3</sup> (8 horas) Piel
Ammonium chloride		STEL: 20 mg/m <sup>3</sup> 15 min TWA: 10 mg/m <sup>3</sup> 8 hr	TWA / VME: 10 mg/m <sup>3</sup> (8 heures).	TWA: 10 mg/m <sup>3</sup> 8 uren STEL: 20 mg/m <sup>3</sup> 15 minuten	STEL / VLA-EC: 20 mg/m <sup>3</sup> (15 minutos). TWA / VLA-ED: 10 mg/m <sup>3</sup> (8 horas)

Component	Italy	Germany	Portugal	The Netherlands	Finland
Sodium azide	TWA: 0.1 mg/m <sup>3</sup> 8 ore. Media Ponderata nel Tempo STEL: 0.3 mg/m <sup>3</sup> 15 minuti. Breve termine Pelle	MAK 0.2 mg/m <sup>3</sup> (inhalable)	STEL: 0.3 mg/m <sup>3</sup> 15 minutos Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm TWA: 0.1 mg/m <sup>3</sup> 8 horas Pele	huid STEL: 0.3 mg/m <sup>3</sup> 15 minuten TWA: 0.1 mg/m <sup>3</sup> 8 uren	TWA: 0.1 mg/m <sup>3</sup> 8 tunteina STEL: 0.3 mg/m <sup>3</sup> 15 minuutteina Iho
Ammonium chloride			STEL: 20 mg/m <sup>3</sup> 15 minutos TWA: 10 mg/m <sup>3</sup> 8 horas		

Component	Austria	Denmark	Switzerland	Poland	Norway
Sodium azide	Haut MAK-KZGW: 0.3 mg/m <sup>3</sup> 15 Minuten MAK-TMW: 0.1 mg/m <sup>3</sup> 8 Stunden	TWA: 0.1 mg/m <sup>3</sup> 8 timer Hud	STEL: 0.4 mg/m <sup>3</sup> 15 Minuten TWA: 0.2 mg/m <sup>3</sup> 8 Stunden	STEL: 0.3 mg/m <sup>3</sup> 15 minutach TWA: 0.1 mg/m <sup>3</sup> 8 godzinach	TWA: 0.1 mg/m <sup>3</sup> 8 timer STEL: 0.3 mg/m <sup>3</sup> 15 minutter. value from the regulation
Ammonium chloride		TWA: 10 mg/m <sup>3</sup> 8 timer	TWA: 3 mg/m <sup>3</sup> 8 Stunden	STEL: 20 mg/m <sup>3</sup> 15 minutach TWA: 10 mg/m <sup>3</sup> 8 godzinach	TWA: 10 mg/m <sup>3</sup> 8 timer STEL: 20 mg/m <sup>3</sup> 15 minutter. set equal to the limit value for Nuisance dust

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

Component	Bulgaria	Croatia	Ireland	Cyprus	Czech Republic
Sodium azide	TWA: 0.1 mg/m <sup>3</sup> STEL : 0.3 mg/m <sup>3</sup> Skin notation	kože TWA-GVI: 0.1 mg/m <sup>3</sup> 8 satima. STEL-KGVI: 0.3 mg/m <sup>3</sup> 15 minutama.	TWA: 0.1 mg/m <sup>3</sup> 8 hr. STEL: 0.3 mg/m <sup>3</sup> 15 min Skin	Skin-potential for cutaneous absorption STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> 8 hodinách. Potential for cutaneous absorption Ceiling: 0.3 mg/m <sup>3</sup>
Ammonium chloride	TWA: 10.0 mg/m <sup>3</sup>	TWA-GVI: 10 mg/m <sup>3</sup> 8 satima. STEL-KGVI: 20 mg/m <sup>3</sup> 15 minutama.	TWA: 10 mg/m <sup>3</sup> 8 hr. fume STEL: 20 mg/m <sup>3</sup> 15 min		TWA: 5 mg/m <sup>3</sup> 8 hodinách. fume Ceiling: 10 mg/m <sup>3</sup> fume

Component	Estonia	Gibraltar	Greece	Hungary	Iceland
Sodium azide	Nahk TWA: 0.1 mg/m <sup>3</sup> 8 tundides. STEL: 0.3 mg/m <sup>3</sup> 15 minutites.	Skin notation TWA: 0.1 mg/m <sup>3</sup> 8 hr STEL: 0.3 mg/m <sup>3</sup> 15 min	STEL: 0.1 ppm STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 ppm TWA: 0.3 mg/m <sup>3</sup>	STEL: 0.3 mg/m <sup>3</sup> 15 percekben. CK TWA: 0.1 mg/m <sup>3</sup> 8 órában. AK	STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup> 8 klukkustundum. Skin notation
Ammonium chloride			STEL: 20 mg/m <sup>3</sup> TWA: 10 mg/m <sup>3</sup>		TWA: 10 mg/m <sup>3</sup> 8 klukkustundum. fume Ceiling: 20 mg/m <sup>3</sup> fume

Component	Latvia	Lithuania	Luxembourg	Malta	Romania
Sodium azide	skin - potential for cutaneous exposure STEL: 0.3 mg/m <sup>3</sup> TWA: 0.1 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> IPRD Oda STEL: 0.3 mg/m <sup>3</sup>	Possibility of significant uptake through the skin TWA: 0.1 mg/m <sup>3</sup> 8 Stunden STEL: 0.3 mg/m <sup>3</sup> 15 Minuten	possibility of significant uptake through the skin TWA: 0.1 mg/m <sup>3</sup> STEL: 0.3 mg/m <sup>3</sup> 15 minuti	Skin notation TWA: 0.1 mg/m <sup>3</sup> 8 ore STEL: 0.3 mg/m <sup>3</sup> 15 minute
Ammonium chloride	TWA: 10 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup> IPRD			TWA: 5 mg/m <sup>3</sup> 8 ore STEL: 10 mg/m <sup>3</sup> 15 minute

Component	Russia	Slovak Republic	Slovenia	Sweden	Turkey
Sodium azide		Ceiling: 0.3 mg/m <sup>3</sup> Potential for cutaneous absorption TWA: 0.1 mg/m <sup>3</sup>	TWA: 0.1 mg/m <sup>3</sup> 8 urah Koža STEL: 0.3 mg/m <sup>3</sup> 15 minutah	Binding STEL: 0.3 mg/m <sup>3</sup> 15 minuter TLV: 0.1 mg/m <sup>3</sup> 8 timmar. NGV	Deri TWA: 0.1 mg/m <sup>3</sup> 8 saat STEL: 0.3 mg/m <sup>3</sup> 15 dakika
Ammonium chloride	MAC: 10 mg/m <sup>3</sup>				

## Biological limit values

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

## Monitoring methods

BS EN 14042:2003 Title Identifier: Workplace atmospheres. Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents.

## Derived No Effect Level (DNEL) / Derived Minimum Effect Level (DMEL)

See table for values

Component	Acute effects local (Dermal)	Acute effects systemic (Dermal)	Chronic effects local (Dermal)	Chronic effects systemic (Dermal)
Sodium azide 26628-22-8 ( <0.1 )				DNEL = 46.7µg/kg bw/day
Ammonium chloride 12125-02-9 ( <0.001 )				DNEL = 128.9mg/kg bw/day DNEL = 190mg/kg bw/day

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

Component	Acute effects local (Inhalation)	Acute effects systemic (Inhalation)	Chronic effects local (Inhalation)	Chronic effects systemic (Inhalation)
Sodium azide 26628-22-8 ( <0.1 )				DNEL = 0.164mg/m <sup>3</sup>
Ammonium chloride 12125-02-9 ( <0.001 )				DNEL = 33.5mg/m <sup>3</sup> DNEL = 43.97mg/m <sup>3</sup>

## Predicted No Effect Concentration (PNEC)

See values below.

Component	Fresh water	Fresh water sediment	Water Intermittent	Microorganisms in sewage treatment	Soil (Agriculture)
Sodium azide 26628-22-8 ( <0.1 )	PNEC = 0.35µg/L	PNEC = 16.7µg/kg sediment dw	PNEC = 3.5µg/L	PNEC = 30µg/L	
Ammonium chloride 12125-02-9 ( <0.001 )	PNEC = 0.25mg/L PNEC = 1.2mg/L		PNEC = 0.43mg/L PNEC = 1.2mg/L	PNEC = 16.2mg/L	PNEC = 0.163mg/kg soil dw PNEC = 50.7mg/kg soil dw

Component	Marine water	Marine water sediment	Marine water Intermittent	Food chain	Air
Sodium azide 26628-22-8 ( <0.1 )	PNEC = 15ng/L	PNEC = 0.72µg/kg sediment dw	PNEC = 150ng/L		
Ammonium chloride 12125-02-9 ( <0.001 )	PNEC = 0.025mg/L PNEC = 11.2mg/L				

## 8.2. Exposure controls

### Engineering Measures

None under normal use conditions.

### Personal protective equipment

#### Eye Protection

Wear safety glasses with side shields (or goggles) (European standard - EN 166)

#### Hand Protection

Protective gloves

Glove material	Breakthrough time	Glove thickness	EU standard	Glove comments
Natural rubber	See manufacturers	-	EN 374	(minimum requirement)
Nitrile rubber	recommendations			
Neoprene				
PVC				

#### Skin and body protection

Long sleeved clothing.

Inspect gloves before use. observe the instructions regarding permeability and breakthrough time which are provided by the supplier of the gloves. (Refer to manufacturer/supplier for information) gloves are suitable for the task: Chemical compatability, Dexterity, Operational conditions, User susceptibility, e.g. sensitisation effects, also take into consideration the specific local conditions under which the product is used, such as the danger of cuts, abrasion. gloves with care avoiding skin contamination.

#### Respiratory Protection

No protective equipment is needed under normal use conditions.

### Large scale/emergency use

Use a NIOSH/MSHA or European Standard EN 136 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced

**Recommended Filter type:** Particle filter

### Small scale/Laboratory use

Maintain adequate ventilation

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

Environmental exposure controls No information available.

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

### 9.1. Information on basic physical and chemical properties

Physical State	Liquid	
Appearance	Clear, colorless	
Odor	Mild	
Odor Threshold	No data available	
Melting Point/Range	No data available	
Softening Point	No data available	
Boiling Point/Range	100 °C / 212 °F	
Flammability (liquid)	No data available	
Flammability (solid,gas)	Not applicable	Liquid
Explosion Limits	No data available	
Flash Point	No data available	Method - No information available
Autoignition Temperature	No data available	
Decomposition Temperature	No data available	
pH	7.4	
Viscosity	No data available	
Water Solubility	Soluble in water	
Solubility in other solvents	No information available	
Partition Coefficient (n-octanol/water)		
Component	log Pow	
Ammonium chloride	-4.38	
Vapor Pressure	No data available	
Density / Specific Gravity	0.999 kg/L @ 20C	
Bulk Density	Not applicable	Liquid
Vapor Density	No data available	(Air = 1.0)
Particle characteristics	Not applicable (liquid)	

### 9.2. Other information

Explosive Properties Not explosive

## SECTION 10: STABILITY AND REACTIVITY

### 10.1. Reactivity

None known, based on information available

### 10.2. Chemical stability

Stable under recommended storage conditions.

### 10.3. Possibility of hazardous reactions

Hazardous Polymerization Hazardous polymerization does not occur.  
Hazardous Reactions None under normal processing.

### 10.4. Conditions to avoid

Incompatible products. Excess heat.

### 10.5. Incompatible materials

None known.

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

## 10.6. Hazardous decomposition products

Carbon oxides.

## SECTION 11: TOXICOLOGICAL INFORMATION

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Product Information

##### (a) acute toxicity;

Oral

Based on available data, the classification criteria are not met

Dermal

Based on available data, the classification criteria are not met

Inhalation

Based on available data, the classification criteria are not met

#### Toxicology data for the components

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Sodium azide	LD50 = 27 mg/kg ( Rat )	-	LC50 0.054 - 0.52 mg/L ( Rat ) 4 h
Water	-	-	-
Ammonium chloride	1650 mg/kg ( Rat )	> 2000 mg/kg	-

(b) skin corrosion/irritation; No data available

(c) serious eye damage/irritation; No data available

##### (d) respiratory or skin sensitization;

Respiratory

No data available

Skin

No data available

(e) germ cell mutagenicity; No data available

(f) carcinogenicity; No data available

There are no known carcinogenic chemicals in this product

(g) reproductive toxicity; No data available

(h) STOT-single exposure; No data available

(i) STOT-repeated exposure; No data available

Target Organs

None known.

(j) aspiration hazard; No data available

Symptoms / effects, both acute and delayed No information available.

### 11.2. Information on other hazards

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

**Endocrine Disrupting Properties** Assess endocrine disrupting properties for human health. This product does not contain any known or suspected endocrine disruptors.

## SECTION 12: ECOLOGICAL INFORMATION

### 12.1. Toxicity

#### Ecotoxicity effects

The environmental impact of this product has not been fully investigated.

Component	Freshwater Fish	Water Flea	Freshwater Algae
Sodium azide	LC50: = 0.7 mg/L, 96h (Lepomis macrochirus) LC50: = 0.8 mg/L, 96h (Oncorhynchus mykiss) LC50: = 5.46 mg/L, 96h flow-through (Pimephales promelas)		
Ammonium chloride	Cyprinus carpio: LC50 = 209 mg/L	EC50 = 202 mg/L/24h	-

Component	Microtox	M-Factor
Sodium azide		1
Ammonium chloride	-	

### 12.2. Persistence and degradability

#### Persistence

Soluble in water, Persistence is unlikely, based on information available.

### 12.3. Bioaccumulative potential

Bioaccumulation is unlikely

Component	log Pow	Bioconcentration factor (BCF)
Ammonium chloride	-4.38	No data available

### 12.4. Mobility in soil

The product is water soluble, and may spread in water systems. Will likely be mobile in the environment due to its water solubility. Highly mobile in soils

### 12.5. Results of PBT and vPvB assessment

No data available for assessment.

### 12.6. Endocrine disrupting properties

#### Endocrine Disruptor Information

This product does not contain any known or suspected endocrine disruptors

### 12.7. Other adverse effects

#### Persistent Organic Pollutant Ozone Depletion Potential

This product does not contain any known or suspected substance  
This product does not contain any known or suspected substance

## SECTION 13: DISPOSAL CONSIDERATIONS

### 13.1. Waste treatment methods

#### Waste from Residues/Unused Products

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. Consult local, regional, and national hazardous waste regulations to ensure complete and accurate classification.

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

<b>Contaminated Packaging</b>	Dispose of contents/containers in accordance with local regulations.
<b>European Waste Catalogue (EWC)</b>	According to the European Waste Catalog, Waste Codes are not product specific, but application specific.
<b>Other Information</b>	Waste codes should be assigned by the user based on the application for which the product was used.
<b>Switzerland - Waste Ordinance</b>	Disposal should be in accordance with applicable regional, national and local laws and regulations. Ordinance on the Avoidance and the Disposal of Waste (Waste Ordinance, ADWO) SR 814.600 <a href="https://www.fedlex.admin.ch/eli/cc/2015/891/en">https://www.fedlex.admin.ch/eli/cc/2015/891/en</a>

## SECTION 14: TRANSPORT INFORMATION

**IMDG/IMO** Not regulated

**14.1. UN number**

**14.2. UN proper shipping name**

**14.3. Transport hazard class(es)** N/A

**14.4. Packing group**

**ADR** Not regulated

**14.1. UN number**

**14.2. UN proper shipping name**

**14.3. Transport hazard class(es)**

**14.4. Packing group**

**IATA** Not regulated

**14.1. UN number**

**14.2. UN proper shipping name**

**14.3. Transport hazard class(es)** N/A

**14.4. Packing group**

**14.5. Environmental hazards** No hazards identified

**14.6. Special precautions for user** No special precautions required

**14.7. Maritime transport in bulk according to IMO instruments** Not applicable, packaged goods

## SECTION 15: REGULATORY INFORMATION

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

### International Inventories

China, X = listed, Australia, U.S.A. (TSCA), Canada (DSL/NDL), Europe (EINECS/ELINCS/NLP), Australia (AICS), Korea (KECL), China (IECSC), Japan (ENCS), Philippines (PICCS), Taiwan (TCSI), Japan (ISHL), New Zealand (NZIoC), Japan (ISHL). US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

Component	CAS No	EINECS	ELINCS	NLP	IECSC	TCSI	KECL	ENCS	ISHL
Sodium azide	26628-22-8	247-852-1	-	-	X	X	KE-31357	X	X
Water	7732-18-5	231-791-2	-	-	X	X	KE-35400	X	-
Ammonium chloride	12125-02-9	235-186-4	-	-	X	X	KE-01645	X	X

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

Component	CAS No	TSCA	TSCA Inventory notification - Active-Inactive	DSL	NDSL	AICS	NZIoC	PICCS
Sodium azide	26628-22-8	X	ACTIVE	X	-	X	X	X
Water	7732-18-5	X	ACTIVE	X	-	X	X	X
Ammonium chloride	12125-02-9	X	ACTIVE	X	-	X	X	X

Legend: X - Listed '-' - Not Listed

KECL - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

## Authorisation/Restrictions according to EU REACH

Component	REACH (1907/2006) - Annex XIV - Substances Subject to Authorization	REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances	REACH Regulation (EC 1907/2006) article 59 - Candidate List of Substances of Very High Concern (SVHC)
Ammonium chloride	-	Use restricted. See item 75. (see link for restriction details) Use restricted. See item 65. (see link for restriction details)	-

<https://echa.europa.eu/substances-restricted-under-reach>

Component	CAS No	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements
Sodium azide	26628-22-8	Not applicable	Not applicable
Water	7732-18-5	Not applicable	Not applicable
Ammonium chloride	12125-02-9	Not applicable	Not applicable

## Regulation (EC) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of dangerous chemicals

Not applicable

Take note of Directive 98/24/EC on the protection of the health and safety of workers from the risks related to chemical agents at work .

Take note of Directive 2000/39/EC establishing a first list of indicative occupational exposure limit values

## National Regulations

**UK** - Take note of Control of Substances Hazardous to Health Regulations (COSHH) 2002 and 2005 Amendment

## WGK Classification

Water endangering class = 3 (self classification)

Component	Germany - Water Classification (VwVwS)	Germany - TA-Luft Class
Sodium azide	WGK2	
Ammonium chloride	WGK1	

## Swiss Regulations

Article 4 para. 4 of the Ordinance on the protection of young people in the workplace (SR 822.115) and Article 1 lit. f of the EAER regulation on hazardous work and young people (SR 822.115.2).

Take note on Article 13 Maternity Ordinance (SR 822.111.52) with regards expectant and nursing mothers.

Component	Switzerland - Ordinance on the Reduction of Risk from handling of hazardous substances preparation (SR 814.81)	Switzerland - Ordinance on Incentive Taxes on Volatile Organic Compounds (OVOC)	Switzerland - Ordinance of the Rotterdam Convention on the Prior Informed Consent Procedure
Ammonium chloride 12125-02-9 ( <0.001 )	Prohibited and Restricted Substances		



# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

## 15.2. Chemical safety assessment

Chemical Safety Assessment/Reports (CSA/CSR) are not required for mixtures

## SECTION 16: OTHER INFORMATION

### Full text of H-Statements referred to under sections 2 and 3

H300 - Fatal if swallowed  
H302 - Harmful if swallowed  
H319 - Causes serious eye irritation  
H400 - Very toxic to aquatic life  
H410 - Very toxic to aquatic life with long lasting effects  
EUH032 - Contact with acids liberates very toxic gas

### Legend

**CAS** - Chemical Abstracts Service

**EINECS/ELINCS** - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

**IECSC** - Chinese Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

**DSL/NDL** - Canadian Domestic Substances List/Non-Domestic Substances List

**ENCS** - Japanese Existing and New Chemical Substances

**AICS** - Australian Inventory of Chemical Substances

**NZIoC** - New Zealand Inventory of Chemicals

**WEL** - Workplace Exposure Limit

**ACGIH** - American Conference of Governmental Industrial Hygienists

**DNEL** - Derived No Effect Level

**RPE** - Respiratory Protective Equipment

**LC50** - Lethal Concentration 50%

**NOEC** - No Observed Effect Concentration

**PBT** - Persistent, Bioaccumulative, Toxic

**TWA** - Time Weighted Average

**IARC** - International Agency for Research on Cancer  
Predicted No Effect Concentration (PNEC)

**LD50** - Lethal Dose 50%

**EC50** - Effective Concentration 50%

**POW** - Partition coefficient Octanol:Water

**vPvB** - very Persistent, very Bioaccumulative

**ADR** - European Agreement Concerning the International Carriage of Dangerous Goods by Road

**IMO/IMDG** - International Maritime Organization/International Maritime Dangerous Goods Code

**OECD** - Organisation for Economic Co-operation and Development

**BCF** - Bioconcentration factor

### Key literature references and sources for data

<https://echa.europa.eu/information-on-chemicals>

Suppliers safety data sheet, Chemadvisor - LOLI, Merck index, RTECS

**ICAO/IATA** - International Civil Aviation Organization/International Air Transport Association

**MARPOL** - International Convention for the Prevention of Pollution from Ships

**ATE** - Acute Toxicity Estimate

**VOC** - (volatile organic compound)

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

**Physical hazards** On basis of test data

**Health Hazards** Calculation method

**Environmental hazards** Calculation method

### Training Advice

Chemical hazard awareness training, incorporating labelling, Safety Data Sheets (SDS), Personal Protective Equipment (PPE) and hygiene.

Revision Date 28-Jan-2022

Revision Summary SDS sections updated.

**This safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006. COMMISSION REGULATION (EU) 2020/878 amending Annex II to Regulation (EC) No 1907/2006 .**

**For Switzerland - Compiled in accordance with the technical provisions referred to in Annex 2, Number 3, ChemO (SR 813.11 - Ordinance on Protection against Dangerous Substances and**

# SAFETY DATA SHEET

Infinity Ammonia Standard

Revision Date 28-Jan-2022

---

Preparations).

## Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of Safety Data Sheet**

**SAFETY DATA SHEET**

Revision Date 28-Jan-2022

Revision Number 1

**1. Identification**

**Product Name** Infinity Ammonia Reagent

**Cat No. :** OSR61154

**Synonyms** No information available

**Recommended Use** In vitro diagnostic.

**Uses advised against** Food, drug, pesticide or biocidal product use.

**Details of the supplier of the safety data sheet**

Fisher Diagnostics  
A Division of Fisher Scientific Company, LLC  
A Part of Thermo Fisher Scientific, Inc.  
8365 Valley Pike  
Middletown, VA 22645-1905  
Tel: (800) 528-0494

**Emergency Telephone Number** Chemtrec US: (800) 424-9300  
Chemtrec Outside US: +1 (703) 741-5970

**2. Hazard(s) identification****Classification**

Classification under 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

**Label Elements**

None required

**Hazards not otherwise classified (HNOC)**

None identified

**3. Composition/Information on Ingredients**

**OSHA Hazard Classification** The product contains no substances which at their given concentration, are considered to be hazardous to health

Component	CAS No	Weight %
Tris (hydroxymethyl) aminomethane	77-86-1	0.1<1.0
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	0.1<1.0
Sodium azide	26628-22-8	<0.1
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	0.1<1.0
Water	7732-18-5	90<100

#### 4. First-aid measures

<b>General Advice</b>	Show this safety data sheet to the doctor in attendance.
<b>Eye Contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
<b>Skin Contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. Get medical attention immediately if symptoms occur.
<b>Ingestion</b>	Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
<b>Most important symptoms and effects</b>	None reasonably foreseeable.
<b>Notes to Physician</b>	Treat symptomatically

#### 5. Fire-fighting measures

<b>Suitable Extinguishing Media</b>	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
<b>Unsuitable Extinguishing Media</b>	Water may be ineffective
<b>Flash Point</b>	No information available
<b>Method -</b>	No information available
<b>Autoignition Temperature</b>	No information available
<b>Explosion Limits</b>	
<b>Upper</b>	No data available
<b>Lower</b>	No data available
<b>Sensitivity to Mechanical Impact</b>	No
<b>Sensitivity to Static Discharge</b>	No

#### Specific Hazards Arising from the Chemical

Keep product and empty container away from heat and sources of ignition.

#### Hazardous Combustion Products

Carbon oxides.

Carbon oxides.

#### Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

#### NFPA

**Health**  
0

**Flammability**  
0

**Instability**  
0

**Physical hazards**  
N/A

## 6. Accidental release measures

**Personal Precautions** Ensure adequate ventilation. Use personal protective equipment as required.  
**Environmental Precautions** Should not be released into the environment.

**Methods for Containment and Clean Up** Sweep up and shovel into suitable containers for disposal.

## 7. Handling and storage

**Handling** Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid contact with skin, eyes or clothing. Avoid ingestion and inhalation.

**Storage.** Keep container tightly closed.

## 8. Exposure controls / personal protection

### Exposure Guidelines

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH	Mexico OEL (TWA)
Sodium azide	Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm	Skin (Vacated) Ceiling: 0.1 ppm (Vacated) Ceiling: 0.3 mg/m <sup>3</sup>	Ceiling: 0.1 ppm Ceiling: 0.3 mg/m <sup>3</sup>	Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm

### Legend

ACGIH - American Conference of Governmental Industrial Hygienists

OSHA - Occupational Safety and Health Administration

NIOSH IDLH: NIOSH - National Institute for Occupational Safety and Health

**Other Exposure Guidelines** Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962 (11th Cir., 1992)

**Engineering Measures** None under normal use conditions.

### Personal Protective Equipment

**Eye/face Protection** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

**Skin and body protection** Wear appropriate protective gloves and clothing to prevent skin exposure.

**Respiratory Protection** No protective equipment is needed under normal use conditions.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

## 9. Physical and chemical properties

<b>Physical State</b>	Liquid
<b>Appearance</b>	Clear, colorless
<b>Odor</b>	Mild
<b>Odor Threshold</b>	No information available
<b>pH</b>	8.7
<b>Melting Point/Range</b>	No data available
<b>Boiling Point/Range</b>	100 °C / 212 °F
<b>Flash Point</b>	No information available
<b>Evaporation Rate</b>	No information available
<b>Flammability (solid,gas)</b>	Not applicable
<b>Flammability or explosive limits</b>	
<b>Upper</b>	No data available

<b>Lower</b>	No data available
<b>Vapor Pressure</b>	No information available
<b>Vapor Density</b>	No information available
<b>Density</b>	1.006 kg/L @ 20C
<b>Specific Gravity</b>	No information available
<b>Solubility</b>	Soluble in water
<b>Partition coefficient; n-octanol/water</b>	No data available
<b>Autoignition Temperature</b>	No information available
<b>Decomposition Temperature</b>	No information available
<b>Viscosity</b>	No information available

## 10. Stability and reactivity

<b>Reactive Hazard</b>	None known, based on information available
<b>Stability</b>	Stable under recommended storage conditions.
<b>Conditions to Avoid</b>	Incompatible products.
<b>Incompatible Materials</b>	Strong oxidizing agents
<b>Hazardous Decomposition Products</b>	Carbon oxides
<b>Hazardous Polymerization</b>	Hazardous polymerization does not occur.
<b>Hazardous Reactions</b>	None under normal processing.

## 11. Toxicological information

### Acute Toxicity

#### Product Information

**Oral LD50** Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.

**Dermal LD50** Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.

**Vapor LC50** Based on ATE data, the classification criteria are not met. ATE > 20 mg/l.

#### Component Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Tris (hydroxymethyl) aminomethane	LD50 = 5900 mg/kg ( Rat )	LD50 > 5000 mg/kg ( Rat )	Not listed
Sodium azide	LD50 = 27 mg/kg ( Rat )	-	LC50 0.054 - 0.52 mg/L ( Rat ) 4 h
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	OECD 425 (Rat) LD50 > 5000 mg/kg bw	OECD 402 (Rat) LD50 > 5000 mg/kg bw	Not listed
Water	-	-	-

**Toxicologically Synergistic** No information available

#### Products

#### Delayed and immediate effects as well as chronic effects from short and long-term exposure

**Irritation** No information available

**Sensitization** No information available

**Carcinogenicity** The table below indicates whether each agency has listed any ingredient as a carcinogen.

Component	CAS No	IARC	NTP	ACGIH	OSHA	Mexico
Tris (hydroxymethyl) aminomethane	77-86-1	Not listed	Not listed	Not listed	Not listed	Not listed
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	Not listed	Not listed	Not listed	Not listed	Not listed
Sodium azide	26628-22-8	Not listed	Not listed	Not listed	Not listed	Not listed
1,3-Propanediol,	1185-53-1	Not listed	Not listed	Not listed	Not listed	Not listed

2-amino-2-(hydroxymethyl)-, hydrochloride						
Water	7732-18-5	Not listed	Not listed	Not listed	Not listed	Not listed

**Mutagenic Effects** No information available

**Reproductive Effects** No information available.

**Developmental Effects** No information available.

**Teratogenicity** No information available.

**STOT - single exposure** None known

**STOT - repeated exposure** None known

**Aspiration hazard** No information available

**Symptoms / effects, both acute and delayed** No information available

**Endocrine Disruptor Information** No information available

**Other Adverse Effects** The toxicological properties have not been fully investigated.

## 12. Ecological information

### Ecotoxicity

The environmental impact of this product has not been fully investigated.

Component	Freshwater Algae	Freshwater Fish	Microtox	Water Flea
Sodium azide	Not listed	LC50: = 0.7 mg/L, 96h (Lepomis macrochirus) LC50: = 0.8 mg/L, 96h (Oncorhynchus mykiss) LC50: = 5.46 mg/L, 96h flow-through (Pimephales promelas)	Not listed	Not listed
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	Not listed	Not listed	OECD 209 EC50 > 1000 mg/L (3h)	Daphnia Magna EC50 > 100 mg/L (48h)

**Persistence and Degradability** Soluble in water Persistence is unlikely based on information available.

**Bioaccumulation/ Accumulation** No information available.

**Mobility** Will likely be mobile in the environment due to its water solubility.

Component	log Pow
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	-3.6

## 13. Disposal considerations

**Waste Disposal Methods** This material, as supplied, is not a hazardous waste according to Federal regulations (40 CFR 261). This material could become a hazardous waste if it is mixed with or otherwise comes in contact with a hazardous waste, if chemical additions are made to this material, or if the material is processed or otherwise altered. Consult 40 CFR 261 to determine whether the altered material is a hazardous waste. Disposal should be in accordance with applicable regional, national and local laws and regulations.

## 14. Transport information

**DOT** NOT REGULATED  
**Hazard Class** N/A  
**TDG** Not regulated

<b>IATA</b>	Not regulated
<b>Hazard Class</b>	N/A
<b>IMDG/IMO</b>	Not regulated
<b>Hazard Class</b>	N/A

## 15. Regulatory information

### United States of America Inventory

Component	CAS No	TSCA	TSCA Inventory notification - Active-Inactive	TSCA - EPA Regulatory Flags
Tris (hydroxymethyl) aminomethane	77-86-1	X	ACTIVE	-
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	-	-	-
Sodium azide	26628-22-8	X	ACTIVE	-
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	X	ACTIVE	-
Water	7732-18-5	X	ACTIVE	-

#### Legend:

**TSCA** US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

X - Listed

- - Not Listed

**TSCA 12(b)** - Notices of Export Not applicable

### International Inventories

X = listed, U.S.A. (TSCA), Canada (DSL/NDL), Europe (EINECS/ELINCS/NLP), Australia (AICS), Korea (KECL), China (IECSC), Japan (ENCS), Philippines (PICCS), Taiwan (TCSI), Japan (ISHL), New Zealand (NZIoC), Japan (ISHL).

Component	CAS No	DSL	NDL	EINECS	PICCS	ENCS	ISHL	AICS	IECSC	KECL
Tris (hydroxymethyl) aminomethane	77-86-1	X	-	201-064-4	X	X	X	X	X	KE-01403
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	X	-	206-167-8	-	-		X	-	-
Sodium azide	26628-22-8	X	-	247-852-1	X	X	X	X	X	KE-31357
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	X	-	214-684-5	X	X		X	X	KE-34819
Water	7732-18-5	X	-	231-791-2	X	X		X	X	KE-35400

**KECL** - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)

### U.S. Federal Regulations

#### SARA 313

Component	CAS No	Weight %	SARA 313 - Threshold Values %
Sodium azide	26628-22-8	<0.1	1.0

**SARA 311/312 Hazard Categories** See section 2 for more information

**CWA (Clean Water Act)** Not applicable

**Clean Air Act** Not applicable

**OSHA** - Occupational Safety and Health Administration Not applicable

**CERCLA** Not applicable



Component	Hazardous Substances RQs	CERCLA EHS RQs
Sodium azide	1000 lb	1000 lb

**U.S. State Regulations**

**California Proposition 65**

This product does not contain any Proposition 65 chemicals.

**U.S. State Right-to-Know Regulations**

Component	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Sodium azide	X	X	X	-	X
Water	-	-	X	-	-

**U.S. Department of Transportation**

Reportable Quantity (RQ): Y  
DOT Marine Pollutant N  
DOT Severe Marine Pollutant N

**U.S. Department of Homeland Security**

This product contains the following DHS chemicals:  
**Legend** - STQs = Screening Threshold Quantities, APA = A placarded amount

Component	DHS Chemical Facility Anti-Terrorism Standard
Sodium azide	Theft STQs - 400lb

**Other International Regulations**

**Mexico - Grade**

No information available

**Authorisation/Restrictions according to EU REACH**

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

Component	CAS No	OECD HPV	Persistent Organic Pollutant	Ozone Depletion Potential	Restriction of Hazardous Substances (RoHS)
Tris (hydroxymethyl) aminomethane	77-86-1	Listed	Not applicable	Not applicable	Not applicable
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	Not applicable	Not applicable	Not applicable	Not applicable
Sodium azide	26628-22-8	Not applicable	Not applicable	Not applicable	Not applicable
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	Not applicable	Not applicable	Not applicable	Not applicable
Water	7732-18-5	Listed	Not applicable	Not applicable	Not applicable

Component	CAS No	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements	Rotterdam Convention (PIC)	Basel Convention (Hazardous Waste)
Tris (hydroxymethyl) aminomethane	77-86-1	Not applicable	Not applicable	Not applicable	Not applicable
Pentanedioic acid, 2-oxo-, disodium salt	305-72-6	Not applicable	Not applicable	Not applicable	Not applicable
Sodium azide	26628-22-8	Not applicable	Not applicable	Not applicable	Not applicable
1,3-Propanediol, 2-amino-2-(hydroxymethyl)-, hydrochloride	1185-53-1	Not applicable	Not applicable	Not applicable	Not applicable
Water	7732-18-5	Not applicable	Not applicable	Not applicable	Not applicable

**16. Other information**

**Prepared By** Regulatory Affairs  
Regulatory Affairs Fisher Diagnostics A Division of Fisher Scientific Co. LLC A Part of  
Thermo Fisher Scientific, Inc. Tel: (800) 528-0494

**Revision Date** 28-Jan-2022

**Print Date** 28-Jan-2022

**Revision Summary** This document has been updated to comply with the US OSHA HazCom 2012 Standard replacing the current legislation under 29 CFR 1910.1200 to align with the Globally Harmonized System of Classification and Labeling of Chemicals (GHS).

**Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of SDS**

**SAFETY DATA SHEET**

Revision Date 28-Jan-2022

Revision Number 1

**1. Identification**

**Product Name** Infinity Ammonia Standard

**Cat No. :** OSR61154

**Synonyms** No information available

**Recommended Use** In vitro diagnostic.

**Uses advised against** Food, drug, pesticide or biocidal product use.

**Details of the supplier of the safety data sheet**

Fisher Diagnostics  
A Division of Fisher Scientific Company, LLC  
A Part of Thermo Fisher Scientific, Inc.  
8365 Valley Pike  
Middletown, VA 22645-1905  
Tel: (800) 528-0494

**Emergency Telephone Number** Chemtrec US: (800) 424-9300  
Chemtrec Outside US: +1 (703) 741-5970

**2. Hazard(s) identification****Classification**

Classification under 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

**Label Elements**

None required

**Hazards not otherwise classified (HNOC)**

None identified

**3. Composition/Information on Ingredients**

**OSHA Hazard Classification** The product contains no substances which at their given concentration, are considered to be hazardous to health

Component	CAS No	Weight %
Sodium azide	26628-22-8	<0.1
Water	7732-18-5	90<100
Ammonium chloride	12125-02-9	<0.001

#### 4. First-aid measures

<b>General Advice</b>	Show this safety data sheet to the doctor in attendance.
<b>Eye Contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
<b>Skin Contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
<b>Inhalation</b>	Remove to fresh air. Get medical attention immediately if symptoms occur.
<b>Ingestion</b>	Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
<b>Most important symptoms and effects</b>	None reasonably foreseeable.
<b>Notes to Physician</b>	Treat symptomatically

#### 5. Fire-fighting measures

<b>Suitable Extinguishing Media</b>	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
<b>Unsuitable Extinguishing Media</b>	Water may be ineffective
<b>Flash Point</b>	No information available
<b>Method -</b>	No information available
<b>Autoignition Temperature</b>	No information available
<b>Explosion Limits</b>	
<b>Upper</b>	No data available
<b>Lower</b>	No data available
<b>Sensitivity to Mechanical Impact</b>	No
<b>Sensitivity to Static Discharge</b>	No

#### Specific Hazards Arising from the Chemical

None reasonably foreseeable.

#### Hazardous Combustion Products

Carbon oxides.

Carbon oxides.

#### Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

#### NFPA

**Health**  
0

**Flammability**  
0

**Instability**  
0

**Physical hazards**  
N/A

#### 6. Accidental release measures

<b>Personal Precautions</b>	Ensure adequate ventilation. Use personal protective equipment as required.
<b>Environmental Precautions</b>	Should not be released into the environment.

**Methods for Containment and Clean Up** Sweep up and shovel into suitable containers for disposal.

## 7. Handling and storage

**Handling** Wear personal protective equipment/face protection. Ensure adequate ventilation. Avoid contact with skin, eyes or clothing. Avoid ingestion and inhalation.

**Storage.** Keep container tightly closed.

## 8. Exposure controls / personal protection

### Exposure Guidelines

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH	Mexico OEL (TWA)
Sodium azide	Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm	Skin (Vacated) Ceiling: 0.1 ppm (Vacated) Ceiling: 0.3 mg/m <sup>3</sup>	Ceiling: 0.1 ppm Ceiling: 0.3 mg/m <sup>3</sup>	Ceiling: 0.29 mg/m <sup>3</sup> Ceiling: 0.11 ppm
Ammonium chloride	TWA: 10 mg/m <sup>3</sup> STEL: 20 mg/m <sup>3</sup>	(Vacated) TWA: 10 mg/m <sup>3</sup> (Vacated) STEL: 20 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup> STEL: 20 mg/m <sup>3</sup>	TWA: 10 mg/m <sup>3</sup> STEL: 20 mg/m <sup>3</sup>

### Legend

ACGIH - American Conference of Governmental Industrial Hygienists

OSHA - Occupational Safety and Health Administration

NIOSH IDLH: NIOSH - National Institute for Occupational Safety and Health

**Other Exposure Guidelines** Vacated limits revoked by the Court of Appeals decision in AFL-CIO v. OSHA, 965 F.2d 962 (11th Cir., 1992)

**Engineering Measures** None under normal use conditions.

### Personal Protective Equipment

**Eye/face Protection** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

**Skin and body protection** Wear appropriate protective gloves and clothing to prevent skin exposure.

**Respiratory Protection** No protective equipment is needed under normal use conditions.

**Hygiene Measures** Handle in accordance with good industrial hygiene and safety practice.

## 9. Physical and chemical properties

Physical State	Liquid
Appearance	Clear, colorless
Odor	Mild
Odor Threshold	No information available
pH	7.4
Melting Point/Range	No data available
Boiling Point/Range	100 °C / 212 °F
Flash Point	No information available
Evaporation Rate	No information available
Flammability (solid,gas)	Not applicable
Flammability or explosive limits	
Upper	No data available
Lower	No data available
Vapor Pressure	No information available

Vapor Density	No information available
Density	0.999 kg/L @ 20C
Specific Gravity	No information available
Solubility	Soluble in water
Partition coefficient; n-octanol/water	No data available
Autoignition Temperature	No information available
Decomposition Temperature	No information available
Viscosity	No information available

## 10. Stability and reactivity

Reactive Hazard	None known, based on information available
Stability	Stable under recommended storage conditions.
Conditions to Avoid	Incompatible products.
Incompatible Materials	Strong oxidizing agents
Hazardous Decomposition Products	Carbon oxides
Hazardous Polymerization	Hazardous polymerization does not occur.
Hazardous Reactions	None under normal processing.

## 11. Toxicological information

### Acute Toxicity

#### Product Information

##### Oral LD50

Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.

##### Dermal LD50

Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.

##### Vapor LC50

Based on ATE data, the classification criteria are not met. ATE > 20 mg/l.

#### Component Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Sodium azide	LD50 = 27 mg/kg ( Rat )	-	LC50 0.054 - 0.52 mg/L ( Rat ) 4 h
Water	-	-	-
Ammonium chloride	1650 mg/kg ( Rat )	> 2000 mg/kg	Not listed

**Toxicologically Synergistic** No information available

#### Products

#### Delayed and immediate effects as well as chronic effects from short and long-term exposure

**Irritation** No information available

**Sensitization** No information available

**Carcinogenicity** The table below indicates whether each agency has listed any ingredient as a carcinogen.

Component	CAS No	IARC	NTP	ACGIH	OSHA	Mexico
Sodium azide	26628-22-8	Not listed	Not listed	Not listed	Not listed	Not listed
Water	7732-18-5	Not listed	Not listed	Not listed	Not listed	Not listed
Ammonium chloride	12125-02-9	Not listed	Not listed	Not listed	Not listed	Not listed

**Mutagenic Effects** No information available

**Reproductive Effects** No information available.

**Developmental Effects** No information available.

**Teratogenicity** No information available.

**STOT - single exposure** None known  
**STOT - repeated exposure** None known

**Aspiration hazard** No information available

**Symptoms / effects, both acute and delayed** No information available

**Endocrine Disruptor Information** No information available

**Other Adverse Effects** The toxicological properties have not been fully investigated.

## 12. Ecological information

### Ecotoxicity

The environmental impact of this product has not been fully investigated.

Component	Freshwater Algae	Freshwater Fish	Microtox	Water Flea
Sodium azide	Not listed	LC50: = 0.7 mg/L, 96h (Lepomis macrochirus) LC50: = 0.8 mg/L, 96h (Oncorhynchus mykiss) LC50: = 5.46 mg/L, 96h flow-through (Pimephales promelas)	Not listed	Not listed
Ammonium chloride	-	Cyprinus carpio: LC50 = 209 mg/L	-	EC50 = 202 mg/L/24h

**Persistence and Degradability** Soluble in water Persistence is unlikely based on information available.

**Bioaccumulation/ Accumulation** No information available.

**Mobility** Will likely be mobile in the environment due to its water solubility.

Component	log Pow
Ammonium chloride	-4.38

## 13. Disposal considerations

**Waste Disposal Methods** This material, as supplied, is not a hazardous waste according to Federal regulations (40 CFR 261). This material could become a hazardous waste if it is mixed with or otherwise comes in contact with a hazardous waste, if chemical additions are made to this material, or if the material is processed or otherwise altered. Consult 40 CFR 261 to determine whether the altered material is a hazardous waste. Disposal should be in accordance with applicable regional, national and local laws and regulations.

## 14. Transport information

**DOT** NOT REGULATED  
**Hazard Class** N/A  
**TDG** Not regulated  
**IATA** Not regulated  
**Hazard Class** N/A  
**IMDG/IMO** Not regulated  
**Hazard Class** N/A

## 15. Regulatory information

### United States of America Inventory

Component	CAS No	TSCA	TSCA Inventory notification - Active-Inactive	TSCA - EPA Regulatory Flags
Sodium azide	26628-22-8	X	ACTIVE	-

Water	7732-18-5	X	ACTIVE	-
Ammonium chloride	12125-02-9	X	ACTIVE	-

**Legend:**

TSCA US EPA (TSCA) - Toxic Substances Control Act, (40 CFR Part 710)

X - Listed

- - Not Listed

TSCA 12(b) - Notices of Export Not applicable

**International Inventories**

China, X = listed, Australia, U.S.A. (TSCA), Canada (DSL/NDSL), Europe (EINECS/ELINCS/NLP), Australia (AICS), Korea (KECL), China (IECSC), Japan (ENCS), Philippines (PICCS), Taiwan (TCSI), Japan (ISHL), New Zealand (NZIoC), Japan (ISHL).

Component	CAS No	DSL	NDSL	EINECS	PICCS	ENCS	ISHL	AICS	IECSC	KECL
Sodium azide	26628-22-8	X	-	247-852-1	X	X	X	X	X	KE-31357
Water	7732-18-5	X	-	231-791-2	X	X		X	X	KE-35400
Ammonium chloride	12125-02-9	X	-	235-186-4	X	X	X	X	X	KE-01645

KECL - NIER number or KE number (<http://ncis.nier.go.kr/en/main.do>)**U.S. Federal Regulations****SARA 313**

Component	CAS No	Weight %	SARA 313 - Threshold Values %
Sodium azide	26628-22-8	<0.1	1.0
Ammonium chloride	12125-02-9	<0.001	1.0

SARA 311/312 Hazard Categories See section 2 for more information

**CWA (Clean Water Act)**

Component	CWA - Hazardous Substances	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants
Ammonium chloride	X	5000 lb	-	-

Clean Air Act Not applicable

OSHA - Occupational Safety and Health Administration Not applicable

CERCLA Not applicable

Component	Hazardous Substances RQs	CERCLA EHS RQs
Sodium azide	1000 lb	1000 lb
Ammonium chloride	5000 lb	-

**U.S. State Regulations**

California Proposition 65 This product does not contain any Proposition 65 chemicals.

**U.S. State Right-to-Know Regulations**

Component	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Sodium azide	X	X	X	-	X
Water	-	-	X	-	-
Ammonium chloride	X	X	X	-	X

U.S. Department of Transportation



Reportable Quantity (RQ): Y  
 DOT Marine Pollutant N  
 DOT Severe Marine Pollutant N

### U.S. Department of Homeland Security

This product contains the following DHS chemicals:

**Legend** - STQs = Screening Threshold Quantities, APA = A placarded amount

Component	DHS Chemical Facility Anti-Terrorism Standard
Sodium azide	Theft STQs - 400lb

### Other International Regulations

**Mexico - Grade** No information available

### Authorisation/Restrictions according to EU REACH

Component	REACH (1907/2006) - Annex XIV - Substances Subject to Authorization	REACH (1907/2006) - Annex XVII - Restrictions on Certain Dangerous Substances	REACH Regulation (EC 1907/2006) article 59 - Candidate List of Substances of Very High Concern (SVHC)
Ammonium chloride	-	Use restricted. See item 75. (see link for restriction details) Use restricted. See item 65. (see link for restriction details)	-

<https://echa.europa.eu/substances-restricted-under-reach>

### Safety, health and environmental regulations/legislation specific for the substance or mixture

Component	CAS No	OECD HPV	Persistent Organic Pollutant	Ozone Depletion Potential	Restriction of Hazardous Substances (RoHS)
Sodium azide	26628-22-8	Not applicable	Not applicable	Not applicable	Not applicable
Water	7732-18-5	Listed	Not applicable	Not applicable	Not applicable
Ammonium chloride	12125-02-9	Listed	Not applicable	Not applicable	Not applicable

Component	CAS No	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Major Accident Notification	Seveso III Directive (2012/18/EC) - Qualifying Quantities for Safety Report Requirements	Rotterdam Convention (PIC)	Basel Convention (Hazardous Waste)
Sodium azide	26628-22-8	Not applicable	Not applicable	Not applicable	Not applicable
Water	7732-18-5	Not applicable	Not applicable	Not applicable	Not applicable
Ammonium chloride	12125-02-9	Not applicable	Not applicable	Not applicable	Not applicable

## 16. Other information

### Prepared By

Regulatory Affairs

Regulatory Affairs Fisher Diagnostics A Division of Fisher Scientific Co. LLC A Part of Thermo Fisher Scientific, Inc. Tel: (800) 528-0494

### Revision Date

28-Jan-2022

### Print Date

28-Jan-2022

### Revision Summary

This document has been updated to comply with the US OSHA HazCom 2012 Standard replacing the current legislation under 29 CFR 1910.1200 to align with the Globally Harmonized System of Classification and Labeling of Chemicals (GHS).

### Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of SDS**