

Kit SDS Cover Sheet

Document ID: OSR6119-75: Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Product information

Product name	GGT
Part number	OSR6119, OSR6219

Components

Description	GGT R1 GGT R2
--------------------	------------------

Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.



SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
 Revision Date (year/month/day) 2023/12/23
 Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name GGT R1
Part number Component of P/N OSR6119, OSR6219

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For In Vitro Diagnostic Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter, Inc.
 250 S. Kraemer Blvd
 Brea, CA 92821, U.S.A.
 Tel: 800-854-3633

Supplier

CANADA
 Beckman Coulter Canada LP
 7075 Financial Drive
 Mississauga, ON L5N 6V8
 Canada
 1-800-463-7828

UNITED KINGDOM
 Beckman Coulter (UK) Ltd.
 Oakley Court
 Kingsmead Business Park, London
 Road
 High Wycombe
 United Kingdom HP11 1JU
 01494 441181

AUSTRALIA
 Beckman Coulter Australia Pty Ltd
 23-27 Chaplin Drive
 Lane Cove NSW 2066
 Australia
 ABN 81 002 011 672
 24 Hour emergency contact phone
 number:
 1800 060 881

SWITZERLAND
 Beckman Coulter Eurocenter SA
 22, rue Juste-Olivier, Case Postale
 1044,
 CH-1260 Nyon 1, Switzerland.
 Telephone: +41 (0)22 365 36 11
 Monday through Friday, 9:00 am to
 7:00pm)

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00

e-mail address

SDSNT@beckman.com
Further information Contact:
Customer support Unit, Beckman Coulter Ireland Inc.
Technical Service Department Tel. +001-800-854-3633 (PST)
E-mail: Techsupportuk@beckman.com

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887
Tel +353 (0)65 683 1170; 08:00 - 16:30 hrs Mon-Thur, 08:00 - 15:30 hrs Fri (GMT) Tel +001-800-223-0130 (PST)

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description

In vitro diagnostic reagent.

Colorless; Liquid; Odorless

Classification according to EC 1272/2008 (CLP/GHS)

Skin Sensitization Category 1, H317

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS

Skin Irritation Category 3

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS

Hazardous ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)

Tris(hydroxymethyl)- aminomethane

Pictogram



Signal word

WARNING

Hazard statements

H316 Causes mild skin irritation.

H317 May cause an allergic skin reaction.

Precautionary statements

Prevention

P261 Avoid breathing vapours.

P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves, protective clothing and eye/face protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P332+P313 If skin irritation occurs: Get medical advice/attention.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 Take off contaminated clothing and wash it before use.

Storage

None

Disposal

P501 Dispose of contents/container in accordance with local/national regulations

Product label will display most significant precautionary statements.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 2 Hazards identification (Continued)

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note
N-Glycylglycine CAS # 556-50-3 EINECS # 209-127-8 Index # Not available	1 - 5	Eye Irrit. 2, H319	Eye Irrit. 2, H319	3, 8
Tris(hydroxymethyl)- aminomethane CAS # 77-86-1 EINECS # 201-064-4 Index # Not available	1 - 5	Eye Irrit. 2, H319 STOT SE 3, H335 Skin Irrit. 2, H315	Eye Irrit. 2, H319 STOT SE 3, H335 Skin Irrit. 2, H315	
reaction mass of: 5-chloro-2-methyl-4-isothiazolin -3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) CAS # 55965-84-9 EINECS # Not available Index # 613-167-00-5	< 0.05	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 M-factor Acute=100 M-factor Chronic=100 Skin Corr. 1C, H314 Skin Sens. 1A, H317 EUH071 Specific Concentration Limit (SCL) Skin Irrit. 2 H315 >= 0.06% - < 0.6% Eye Dam. 1 H318 >= 0.6% Skin Corr. 1C H314 >= 0.6% Eye Irrit. 2 H319 >= 0.06% - < 0.6% Skin Sens. 1A H317 >= 0.0015% Acute Toxicity Estimates (ATE) ATE Dermal = 87.12 mg/kg ATE Inhalation - Vapors = 0.5 mg/L ATE Oral = 53 mg/kg	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 Skin Corr. 1C, H314 Skin Sens. 1A, H317	9

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 3 Composition and information on ingredients (Continued)

3 - Health hazard

8 - Present at concentration below the cut-off limits.

9 - Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6] (3:1) is the active ingredient of ProClin 300.

See section 8 for available Occupational exposure limits

See Section 15 for additional regulatory information

See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

Inhalation

If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.

Eye contact

If product enters eyes, rinse eyes gently with water as a precaution.

Skin contact

In case of skin contact, rinse with plenty of water. Remove contaminated clothing and shoes. If pain or irritation occurs, obtain medical advice/attention.

Ingestion

If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

4.2 Most important symptoms and effects, both acute and delayed

May cause an allergic skin reaction.

May cause mild skin irritation.

See Section 11 Toxicological Information for more detailed health information.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available. Refer to Section 4.1.

Section 5 Firefighting measures

5.1 Extinguishing media

In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam.

For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture

Special fire and explosion hazards

No special hazards determined.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product (an aqueous solution).

5.3 Advice for firefighters

Protective equipment

Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information

No further relevant information available.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions Observe general safety guidelines for protection; avoid eye and skin contact.
Wear protective gloves, protective clothing and eye/face protection.

6.2 Environmental precautions Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.
Dispose of contents/container in accordance with local regulations

6.3 Methods and material for containment and cleaning up

Spill and leak procedures Absorb spilled material with an appropriate inert, non-flammable absorbent and dispose according to local regulations.

6.4 Reference to other sections Refer sections 8 and 13.

Section 7 Handling and storage

7.1 Precautions for safe handling Use good laboratory procedures; avoid eye and skin contact.

7.2 Conditions for safe storage, including any incompatibilities

Store at 2 to 8°C, as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).

7.3 Specific end uses No further relevant information available.

Section 8 Exposure controls and personal protection

8.1 Control parameters

Exposure limits

US OSHA None established

ACGIH None established

ACGIH Biological Exposure Indices (BEI)
None established

DFG MAK None established

Ireland None established

IOELVs
None established

NIOSH None established

China None established

Croatia None established

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 8 Exposure controls and personal protection (Continued)

Japan	None established
Sweden (AFS 2015:7 and amendments)	None established
Turkey	None established
8.2 Exposure controls	
Engineering controls	No special engineering controls are required. Use with good general ventilation.
Eye protection	Safety glasses or chemical goggles should be worn to prevent eye contact. Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.
Skin protection	Wear impervious gloves such as Nitrile or equivalent and protective clothing. Refer to U.S. OSHA 29 CFR 1910.138, European Standard EN 374, EN 14605:2005+A1:2009 or appropriate government standards.
Respiratory protection	Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties			
Physical state	Liquid	Density and/or relative density	1.02
Color	Colorless	Solubility	
Odor	Odorless	Water	Miscible
pH	8.4	Organic	Not determined
Freezing point	Similar to water, approximately 0°C	Partition coefficient n-octanol/water (log value)	Not determined
Boiling point or initial boiling point and boiling range	Similar to water, approximately 100°C	Auto-ignition temp.	Product is not selfigniting
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	Similar to water, approximately 23 hPa
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not applicable		
Relative vapor density	Not determined		

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 9 Physical and chemical properties (Continued)

Particle characteristics Not applicable

9.2 Other information

Information with regard to physical hazard classes

Not applicable

No further relevant information available.

Other safety characteristics

No further relevant information available.

Section 10 Stability and reactivity

10.1 Reactivity No further relevant information available.

10.2 Chemical stability The product is stable in accordance with recommended storage conditions.

10.3 Possibility of hazardous reactions

No further relevant information available.

10.4 Conditions to avoid

To maintain product performance keep away from strong acids, strong bases, strong oxidizers.

Avoid exposure to heat and direct sunlight.

10.5 Incompatible materials

No further relevant information available.

10.6 Hazardous decomposition products

No decomposition products posing significant hazards would be expected from this product (an aqueous solution).

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

Tris(hydroxymethyl)-
aminomethane
CAS # 77-86-1

Dermal LD50 Rat >5000 mg/kg (ECHA_API); Oral LD50 Rat 5900 mg/kg (NLM_CIP)

reaction mass of: 5-chloro-
2-methyl-4-isothiazolin
-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC#
220-239-6](3:1)
CAS # 55965-84-9

Dermal LD50 Rabbit 87.12 mg/kg (ECHA_API); Oral LD50 Rat 53 mg/kg (NLM_CIP)

Primary routes of exposure

Eye contact, ingestion, inhalation, and skin contact.

Acute toxicity

Not classified based on available data.

Skin corrosion/irritation

May cause mild skin irritation.

Serious eye damage/irritation

Not classified based on available data.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 11 Toxicological information (Continued)

Respiratory or skin sensitisation	May cause sensitization by skin contact.
Germ cell mutagenicity	Not classified based on available data.
Carcinogenicity	No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.
Reproductive toxicity	Not classified based on available data.
Specific target organ toxicity (STOT) – single exposure	Not classified based on available data.
Specific target organ toxicity (STOT) – repeated exposure	Not classified based on available data.
Aspiration hazard	Not classified based on available data.
11.2 Information on other hazards	
Endocrine disrupting properties	This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).
Other information	No further relevant information available.

Section 12 Ecological information

12.1 Toxicity	
Fresh water species	No information available.
Microtox/organisms	No information available.
Water flea	No information available.
Fresh water algae	No information available.
12.2 Persistence and degradability	Not determined for the product.
12.3 Bioaccumulative potential	Not determined for the product.
12.4 Mobility in soil	Not determined for the product.
12.5 Results of PBT and vPvB assessment	Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 12 Ecological information (Continued)

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects

This product contains environmentally hazardous substance below the cutoff level. Refer section 3 for ingredient information. Do not allow undiluted product to enter sewer/surface or ground water.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

Package disposal

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Additional information

Suggested European waste catalogue 18 01 07 - chemicals other than those mentioned in 18 01 06. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

14.1 UN/ID number: Not regulated for transportation

14.2 UN proper shipping name: Not regulated for transportation

14.3 Transport hazard class(es): Not regulated for transportation

14.4 Packing group: Not regulated for transportation

14.5 Environmental hazards: Not regulated for transportation

14.6 Special precautions for user: None

14.7 Maritime transport in bulk according to IMO instruments: Not applicable

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 15 Regulatory information (Continued)

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

No ingredients listed.

New Jersey Dept. of Health Right To Know (RTK) List

No ingredients listed.

Pennsylvania Right To Know (RTK) List

No ingredients listed.

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany)

WGK 1, low water endangering

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 15 Regulatory information (Continued)

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

Not applicable.

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product is exempt from WHMIS label and SDS requirements.

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

No ingredients listed.

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Turkey-REACH - KKDIK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 0 Health: 2 Reactivity with water: 0 Physical contact: 2	Code 0=None 1=Slight 2=Caution 3=Severe
--------------------------------------	---	--

Revision changes

Updated sections 1, 2, 3, 4, 8 and 15

Document version and issue/revision date

Revision Date (year/month/day) 2023/12/23

Last Revision Date (year/month/day) 2023/05/23

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 16 Other information (Continued)

Hazard Classification Procedure

Document ID: OSR6119-75

Version: 09

This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.

Description of hazard class and hazard statements from Section 3

Aquatic Acute 1 - Aquatic Hazard Acute, Category 1

Acute Tox. Dermal 2 - Acute Toxicity Dermal, Category 2

Acute Tox. Inhal. 2 - Acute Toxicity Inhalation, Category 2

Acute Tox. Oral 3 - Acute Toxicity Oral, Category 3

Eye Dam. 1 - Eye Damage Category 1

Eye Irrit. 2 - Eye Irritation Category 2

Aquatic Longterm 1 - Aquatic Hazard Long term, Category 1

Skin Corr. 1C - Skin Corrosion Category 1C

Skin Irrit. 2 - Skin Irritation Category 2

Skin Sens. 1A - Skin Sensitization Category 1A

STOT SE 3 - Specific Target Organ Toxicity Single Exposure Category 3

EUH071 - Corrosive to the respiratory tract.

H301 - Toxic if swallowed.

H310 - Fatal in contact with skin.

H314 - Causes severe skin burns and eye damage.

H315 - Causes skin irritation.

H317 - May cause an allergic skin reaction.

H318 - Causes serious eye damage.

H319 - Causes serious eye irritation.

H330 - Fatal if inhaled.

H335 - May cause respiratory irritation.

H400 - Very toxic to aquatic life.

H410 - Very toxic to aquatic life with long lasting effects.

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)

ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail

CLP - Classification, Labeling and Packaging

DFGMAK - Republic Germany's maximum exposure limit

EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms

GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)

HCS - Hazard Communication Standard

IARC - International Agency for Research on Cancer

IATA DGR - International Air Transport Association Dangerous Goods Regulation

ICAO - International Civil Aviation Organization

IDLH - Immediately Dangerous to Life or Health

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 16 Other information (Continued)

IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL – Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA – Time Weighted Average
UN GHS - United Nations Globally Harmonized System
US DOT - United States Department of Transportation
US OSHA - United States Occupational Safety and Health Administration
vPvB - very Persistent and very Bioaccumulative substances
WHMIS - Workplace Hazardous Material Information System

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.



SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
 Revision Date (year/month/day) 2023/12/23
 Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name GGT R2
Part number Component of P/N OSR6119, OSR6219

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use For In Vitro Diagnostic Use. See product literature for details.

1.3 Details of the supplier of the safety data sheet

Manufacturer

Beckman Coulter, Inc.
 250 S. Kraemer Blvd
 Brea, CA 92821, U.S.A.
 Tel: 800-854-3633

Supplier

CANADA
 Beckman Coulter Canada LP
 7075 Financial Drive
 Mississauga, ON L5N 6V8
 Canada
 1-800-463-7828

UNITED KINGDOM
 Beckman Coulter (UK) Ltd.
 Oakley Court
 Kingsmead Business Park, London
 Road
 High Wycombe
 United Kingdom HP11 1JU
 01494 441181

AUSTRALIA
 Beckman Coulter Australia Pty Ltd
 23-27 Chaplin Drive
 Lane Cove NSW 2066
 Australia
 ABN 81 002 011 672
 24 Hour emergency contact phone
 number:
 1800 060 881

SWITZERLAND
 Beckman Coulter Eurocenter SA
 22, rue Juste-Olivier, Case Postale
 1044,
 CH-1260 Nyon 1, Switzerland.
 Telephone: +41 (0)22 365 36 11
 Monday through Friday, 9:00 am to
 7:00pm)

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

NEW ZEALAND
Beckman Coulter NZ
Unit J, 33 Walmsley Road, Otahuhu,
Auckland 1062, New Zealand
Hours available: 08:30 - 17:00

ICELAND / ÍSLAND
Beckman Coulter AB
Ekbacksvägen 28
168 69 Bromma
Sweden
Phone No.: +46 80564 85 900
Hours available: 08.00-16.30

MALTA
DX Distributor:
Cherubino Ltd
DELFI Building, Sliema Road, Gzira,
GZR 1637
Telephone: +356 21343270
Hours available: 08:30 – 17:00

e-mail address

SDSNT@beckman.com
Further information Contact:
Customer support Unit, Beckman Coulter Ireland Inc.
Technical Service Department Tel. +001-800-854-3633 (PST)
E-mail: Techsupportuk@beckman.com

1.4 Emergency telephone number

Telephone number (24H)

Chemtrec Emergency Tel No. U.S.A. 800-424-9300, International (001) 703-527-3887
Tel +353 (0)65 683 1170; 08:00 - 16:30 hrs Mon-Thur, 08:00 - 15:30 hrs Fri (GMT) Tel +001-800-223-0130 (PST)

Distributor and emergency phone no.

Refer to attached list, Document ID: [472050](#), for local distributor and emergency phone numbers.

UNITED STATES - Emergency Phone (24h): Chemtrec (800) 424-9300, International (001) 703-527-3887

CANADA - Poison Centre: 1-844-764-7669; Centre antipoison du Québec: 1-800-463-5060

UNITED KINGDOM - For UK and Scotland: Emergency Call 999

IRELAND - National Poisons Information Centre Phone No.: Members of Public: +353 (01) 809 2166 (8:00 am to 10:00 pm 7 days a week); Healthcare Professionals: +353 (01) 809 2566 (24 hour service)

AUSTRALIA - 24 Hour emergency contact phone number: 1800 060 881

NEW ZEALAND - 24 Hour emergency number: 0800 446 109

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 1 Identification of the substance/mixture and of the company/undertaking (Continued)

Section 2 Hazards identification

2.1 Classification of the substance or mixture

Product description

In vitro diagnostic reagent.
Light yellow; Liquid; Mild odor

Classification according to EC 1272/2008 (CLP/GHS)

Skin Sensitization Category 1, H317

Classification according to US-OSHA (HCS 29 CFR 1910.1200) and UN GHS

Not classified as hazardous per US-OSHA HCS 2012 and UN GHS

2.2 Label elements

According to EC 1272/2008 (CLP/GHS), US-OSHA and UN GHS

Hazardous ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and
2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1)

Pictogram



Signal word

WARNING

Hazard statements

H317 May cause an allergic skin reaction.

Precautionary statements

Prevention

P261 Avoid breathing vapours.
P272 Contaminated work clothing should not be allowed out of the workplace.
P280 Wear protective gloves, protective clothing and eye/face protection.

Response

P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 Take off contaminated clothing and wash it before use.

Storage

None

Disposal

P501 Dispose of contents/container in accordance with local/national regulations
Product label will display most significant precautionary statements.

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

See Section 11 Toxicological Information for more detailed health information.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 3 Composition and information on ingredients

3.2 Mixtures

Hazardous ingredients:		Hazard classification of pure ingredients		
Chemical name	% by wt.	EU 1272/2008 CLP/GHS	GHS	Note
reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) CAS # 55965-84-9 EINECS # Not available Index # 613-167-00-5	< 0.05	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 M-factor Acute=100 M-factor Chronic=100 Skin Corr. 1C, H314 Skin Sens. 1A, H317 EUH071 Specific Concentration Limit (SCL) Skin Irrit. 2 H315 >= 0.06% - < 0.6% Eye Dam. 1 H318 >= 0.6% Skin Corr. 1C H314 >= 0.6% Eye Irrit. 2 H319 >= 0.06% - < 0.6% Skin Sens. 1A H317 >= 0.0015% Acute Toxicity Estimates (ATE) ATE Dermal = 87.12 mg/kg ATE Inhalation - Vapors = 0.5 mg/L ATE Oral = 53 mg/kg	Acute Tox. Dermal 2, H310 Acute Tox. Inhal. 2, H330 Acute Tox. Oral 3, H301 Aquatic Acute 1, H400 Aquatic Longterm 1, H410 Eye Dam. 1, H318 Skin Corr. 1C, H314 Skin Sens. 1A, H317	9

9 - Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6] (3:1) is the active ingredient of ProClin 300.

See section 8 for available Occupational exposure limits

See Section 15 for additional regulatory information

See Section 16 for description of hazard class and hazard statements

Section 4 First aid measures

4.1 Description of first aid measures

Inhalation

If product is inhaled, move exposed individual to fresh air. If individual is not breathing, begin artificial respiration by trained personnel and obtain medical attention immediately.

Eye contact

If product enters eyes, rinse eyes gently with water as a precaution.

Skin contact

In case of skin contact, rinse with water as a precaution.

Ingestion

If product is ingested, rinse mouth with water. If irritation or discomfort occurs, obtain medical attention immediately.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 4 First aid measures (Continued)

4.2 Most important symptoms and effects, both acute and delayed

May cause an allergic skin reaction.
See Section 11 Toxicological Information for more detailed health information.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available. Refer to Section 4.1.

Section 5 Firefighting measures

5.1 Extinguishing media In case of fire use carbon dioxide (CO₂), dry chemical, water spray or foam.
For large fires use extinguishing media suitable for surrounding fire.

5.2 Special hazards arising from the substance or mixture Special fire and explosion hazards

No special hazards determined.

Hazardous combustion products

No combustion products posing significant hazards are expected from this product (an aqueous solution).

5.3 Advice for firefighters

Protective equipment

Self-contained breathing apparatus is recommended for firefighters in all chemical fire situations.

Additional information

No further relevant information available.

Section 6 Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions

Observe general safety guidelines for protection; avoid eye and skin contact.
Wear protective gloves, protective clothing and eye/face protection.

6.2 Environmental precautions

Contain spill to prevent migration.
Do not allow the undiluted product to enter sewers/surface or ground water.
Dispose of contents/container in accordance with local regulations

6.3 Methods and material for containment and cleaning up

Spill and leak procedures

Absorb spilled material with an appropriate inert, non-flammable absorbent and dispose according to local regulations.

6.4 Reference to other sections

Refer sections 8 and 13.

Section 7 Handling and storage

7.1 Precautions for safe handling Use good laboratory procedures; avoid eye and skin contact.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 7 Handling and storage (Continued)

7.2 Conditions for safe storage, including any incompatibilities

Store at 2 to 8°C, as directed on the product label.
To maintain product quality, store according to the instructions in the product labeling.
Store away from strong acids, strong bases, strong oxidizers and incompatible materials (section 10).

7.3 Specific end uses

No further relevant information available.

Section 8 Exposure controls and personal protection

8.1 Control parameters

Exposure limits

US OSHA None established

ACGIH None established

ACGIH Biological Exposure Indices (BEI)

None established

DFG MAK None established

Ireland None established

IOELVs

None established

NIOSH None established

China None established

Croatia None established

Japan None established

Sweden (AFS 2015:7 and amendments)

None established

Turkey None established

8.2 Exposure controls

Engineering controls

No special engineering controls are required. Use with good general ventilation.

Eye protection

Safety glasses or chemical goggles should be worn to prevent eye contact.

Refer U.S. OSHA 29 CFR 1910.133, European Standard EN166 or appropriate government standards.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 8 Exposure controls and personal protection (Continued)

Skin protection	Wear impervious gloves such as Nitrile or equivalent and protective clothing. Refer to U.S. OSHA 29 CFR 1910.138, European Standard EN 374, EN 14605:2005+A1:2009 or appropriate government standards.
Respiratory protection	Under normal conditions, the use of this product should not require respiratory protection. If overexposure should occur and ventilation is not adequate to maintain airborne concentrations at acceptable levels, the use of respiratory protection should be evaluated by a qualified professional.

Section 9 Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	Liquid	Density and/or relative density	1.05
Color	Light yellow	Solubility	
Odor	Mild odor	Water	Miscible
pH	6	Organic	Not determined
Freezing point	Similar to water, approximately 0°C	Partition coefficient n-octanol/water (log value)	Not determined
Boiling point or initial boiling point and boiling range	Similar to water, approximately 100°C	Auto-ignition temp.	Product is not selfigniting
Flash point	Not applicable	Decomposition temperature	Not determined
Flammability	Not applicable	Vapor pressure	Similar to water, approximately 23 hPa
		Kinematic viscosity	Not determined
Lower and upper explosion limit	Not applicable		
Relative vapor density	Not determined		
Particle characteristics	Not applicable		

9.2 Other information

Information with regard to physical hazard classes

Not applicable
No further relevant information available.

Other safety characteristics

No further relevant information available.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 10 Stability and reactivity

10.1 Reactivity	No further relevant information available.
10.2 Chemical stability	The product is stable in accordance with recommended storage conditions.
10.3 Possibility of hazardous reactions	No further relevant information available.
10.4 Conditions to avoid	To maintain product performance keep away from strong acids, strong bases, strong oxidizers. Avoid exposure to heat and direct sunlight.
10.5 Incompatible materials	No further relevant information available.
10.6 Hazardous decomposition products	No decomposition products posing significant hazards would be expected from this product (an aqueous solution).

Section 11 Toxicological information

11.1 Information on hazard classes

Toxicity data for hazardous ingredients

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) CAS # 55965-84-9	Dermal LD50 Rabbit 87.12 mg/kg (ECHA_API); Oral LD50 Rat 53 mg/kg (NLM_CIP)
---	---

Primary routes of exposure	Eye contact, ingestion, inhalation, and skin contact.
Acute toxicity	Not classified based on available data.
Skin corrosion/irritation	Not classified based on available data.
Serious eye damage/irritation	Not classified based on available data.
Respiratory or skin sensitisation	May cause sensitization by skin contact.
Germ cell mutagenicity	Not classified based on available data.
Carcinogenicity	No ingredients in this product are listed as carcinogens by ACGIH, IARC, NTP, OSHA or 1272/2008 EC regulation.
Reproductive toxicity	Not classified based on available data.
Specific target organ toxicity (STOT) – single exposure	Not classified based on available data.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 11 Toxicological information (Continued)

Specific target organ toxicity (STOT) – repeated exposure

Not classified based on available data.

Aspiration hazard

Not classified based on available data.

11.2 Information on other hazards

Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for health according to REACH Article 57(f).

Other information

No further relevant information available.

Section 12 Ecological information

12.1 Toxicity

Fresh water species

No information available.

Microtox/organisms

No information available.

Water flea

No information available.

Fresh water algae

No information available.

12.2 Persistence and degradability

Not determined for the product.

12.3 Bioaccumulative potential

Not determined for the product.

12.4 Mobility in soil

Not determined for the product.

12.5 Results of PBT and vPvB assessment

Not determined for the product. PBT: Not applicable, vPvB: Not applicable.

12.6 Endocrine disrupting properties

This product does not have substance(s) of endocrine disrupting properties for environment according to REACH Article 57(f).

12.7 Other adverse effects

This product contains environmentally hazardous substance below the cutoff level. Refer section 3 for ingredient information. Do not allow undiluted product to enter sewer/surface or ground water.

Section 13 Disposal considerations

13.1 Waste treatment methods

Product waste disposal

Chemical residues and remains should be routinely handled as special waste. This must be disposed of in compliance with anti-pollution and other laws of the country concerned. To ensure compliance we recommend that you contact the relevant (local) authorities and/or an approved waste-disposal company for information.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Package disposal

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations. If unsure of the applicable requirements, contact the authorities for information.

Additional information

Suggested European waste catalogue 18 01 07 - chemicals other than those mentioned in 18 01 06. Dispose in accordance with national, state and local waste regulations.

Section 14 Transport information

Transportation of this product is not regulated under ICAO, IATA DGR, IMDG, US DOT, European ADR and RID or Canadian TDG.

14.1 UN/ID number: Not regulated for transportation

14.2 UN proper shipping name: Not regulated for transportation

14.3 Transport hazard class(es): Not regulated for transportation

14.4 Packing group: Not regulated for transportation

14.5 Environmental hazards: Not regulated for transportation

14.6 Special precautions for user: None

14.7 Maritime transport in bulk according to IMO instruments: Not applicable

Section 15 Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal and State Regulations

SARA 313 (Section 313, Title III reporting requirements)

No ingredients listed.

CERCLA (The Comprehensive Environmental Response, Compensation, and Liability Act) 40 CFR 302.4

No ingredients listed.

California Proposition 65

Chemical which is known to the State of California to cause cancer

No ingredients listed.

Chemical which is known to the State of California to cause development toxicity

No ingredients listed.

Chemical which is known to the State of California to cause male reproductive toxicity

No ingredients listed.

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 15 Regulatory information (Continued)

Chemical which is known to the State of California to cause female reproductive toxicity

No ingredients listed.

Massachusetts Right To Know (RTK) List

No ingredients listed.

New Jersey Dept. of Health Right To Know (RTK) List

No ingredients listed.

Pennsylvania Right To Know (RTK) List

No ingredients listed.

EU Regulations

This SDS complies with EC Regulations 1907/2006 (REACH) and amendments.

Water Hazard Class (Germany)

WGK 1, low water endangering

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Substances Subject to Suspicious Transactions Reporting

No ingredients listed.

Regulation (EU) 2019/1148 on the marketing and use of explosives precursors - Restricted Explosives Precursors

No ingredients listed.

REACH 1907/2006 EC - Candidate List of Substances of Very High Concern (SVHC)

No ingredients listed.

REACH 1907/2006 EC - Annex XVII – Restrictions on Certain Dangerous Substances

Not applicable.

REACH 1907/2006 EC - Annex XIV - list of substances subject to authorisation

No ingredients listed.

Refer to Section 3

UK Regulations

UK REACH Regulation (as Amended) - List of substances subject to authorisation

Refer to Section 3

Canada

This product is exempt from WHMIS label and SDS requirements.

China

Catalog of Hazardous Chemicals - Hazardous Chemicals

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 15 Regulatory information (Continued)

No ingredients listed.

Inventory - China - Inventory of Existing Chemical Substances (IECSC)

All ingredients are listed or exempted.

Turkey

Turkey-REACH - KKDIK Regulation - Annex 17 – Restrictions

No ingredients listed.

International

UN/FAO/Rotterdam Convention - Chemicals Subject to Prior Informed Consent (PIC)

No ingredients listed.

15.2 Chemical safety assessment

A Chemical Safety Assessment has not been carried out.

Some hazardous ingredients listed in Section 15 are below the cutoff limits of 0.1% for carcinogen, mutagen and reproductive toxin and 1% for other health hazards required for reporting in Section 3.

Section 16 Other information

Beckman Coulter safety rating	Flammability: 0 Health: 2 Reactivity with water: 0 Physical contact: 2	Code 0=None 1=Slight 2=Caution 3=Severe
Revision changes	Updated sections 1, 2, 3, 4, 8 and 15	
Document version and issue/revision date	Revision Date (year/month/day) 2023/12/23 Last Revision Date (year/month/day) 2023/05/23 Document ID: OSR6119-75 Version: 09	
Hazard Classification Procedure	This mixture was classified using the calculation method for human health and environmental hazards. Physical hazards were determined based on the specification.	
Description of hazard class and hazard statements from Section 3	Aquatic Acute 1 - Aquatic Hazard Acute, Category 1 Acute Tox. Dermal 2 - Acute Toxicity Dermal, Category 2 Acute Tox. Inhal. 2 - Acute Toxicity Inhalation, Category 2 Acute Tox. Oral 3 - Acute Toxicity Oral, Category 3 Eye Dam. 1 - Eye Damage Category 1 Aquatic Longterm 1 - Aquatic Hazard Long term, Category 1 Skin Corr. 1C - Skin Corrosion Category 1C Skin Sens. 1A - Skin Sensitization Category 1A EUH071 - Corrosive to the respiratory tract. H301 - Toxic if swallowed.	

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 16 Other information (Continued)

Abbreviations and acronyms

H310 - Fatal in contact with skin.
H314 - Causes severe skin burns and eye damage.
H317 - May cause an allergic skin reaction.
H318 - Causes serious eye damage.
H330 - Fatal if inhaled.
H400 - Very toxic to aquatic life.
H410 - Very toxic to aquatic life with long lasting effects.
ACGIH - American Conference of Governmental Industrial Hygienists (ACGIH)
ADR and RID - European Agreement Concerning the International Carriage of Dangerous Goods by Road and Rail
CLP - Classification, Labeling and Packaging
DFGMAK - Republic Germany's maximum exposure limit
EC50 - Concentration of a substance in an environmental medium expected to produce a certain effect in 50% of test organisms
GHS - Globally Harmonized System of Classification and Labelling of Chemicals (GHS)
HCS - Hazard Communication Standard
IARC - International Agency for Research on Cancer
IATA DGR - International Air Transport Association Dangerous Goods Regulation
ICAO - International Civil Aviation Organization
IDLH - Immediately Dangerous to Life or Health
IMDG - International Maritime Dangerous Goods
IMO - International Maritime Organization
IOELVs - European Unions' Indicative Occupational Exposure Limit Values
LC50 - Concentration of a substance in water causing death (50% of the tested population) to aquatic life
LD50 - Lethal Dose 50%
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety and Health Administration
PBT - Persistent Bioaccumulative and Toxic substances
PEL - Permissible Exposure Limit
SARA - Superfund Amendments and Reauthorization Act
STEL - Short Term Exposure Limit
STLV - Short Term Limit Value
STV - Short Term Value
TDG - Canadian Transportation of Dangerous Goods Regulations
TLV - Threshold Limit Value
TWA - Time Weighted Average
UN GHS - United Nations Globally Harmonized System

SAFETY DATA SHEET

Document ID: OSR6119-75 Version 09
Revision Date (year/month/day) 2023/12/23
Last Revision Date (year/month/day) 2023/05/23

Section 16 Other information (Continued)

US DOT - United States Department of Transportation

US OSHA - United States Occupational Safety and Health Administration

vPvB - very Persistent and very Bioaccumulative substances

WHMIS - Workplace Hazardous Material Information System

For further information, please contact your local Beckman Coulter, Inc. representative.

WHILE BECKMAN COULTER, INC. BELIEVES THE INFORMATION CONTAINED HEREIN IS VALID AND ACCURATE, BECKMAN COULTER, INC. MAKES NO WARRANTY OR REPRESENTATION AS TO ITS VALIDITY, ACCURACY, OR CURRENCY. BECKMAN COULTER, INC. SHALL NOT BE LIABLE OR OTHERWISE RESPONSIBLE IN ANY WAY FOR USE OF EITHER THIS INFORMATION OR MATERIALS TO WHICH IT APPLIES. DISPOSAL OF HAZARDOUS MATERIALS MAY BE SUBJECT TO LOCAL LAWS OR REGULATIONS.