

Beckman Coulter
 MIAMI AND KENTUCKY REAGENTS MANUFACTURING
CERTIFICATE OF ANALYSIS

IDENTIFICATION:	Product Name:	6C® CELL CONTROLS
	Part Number:	628027
	Lot Number (Finished Product):	4220210K
	Lot Number (Sub-assembly):	123176110, 133186110, 143196110
	Date of Manufacture:	13 Aug 2025
	Expiration Date:	01 Nov 2025

TEST	SPECIFICATION	RESULTS	Pass/Fail
Microbial Integrity	No Growth	Meets Specification	Pass
Value Verification	Must meet established limits for appropriate parameters on designated instruments.	Meets Specification	Pass
pH	6.60 - 6.90	Meets Specification	Pass
Osmolality	Level 1: 301-339 Level 2: 318-353 Level 3: 323-358	Meets Specification Meets Specification Meets Specification	Pass Pass Pass
Hemolysis (if applicable)			
Level 1 / Level 2 Level 3	≤+2 hemolysis ≤+3 hemolysis	Meets Specification Meets Specification	Pass Pass
Homogeneity (if applicable)	≤ 2 SD limits	Meets Specification	Pass

The product was manufactured in accordance with the US FDA Quality System Regulation (GMP).

Signed by:	APPROVAL	TITLE	DATE
<i>Pedro Garcia Colon</i>		Manufacturing Quality Supervisor	25-Aug-2025

Signer Name: Pedro Garcia Colon
 Signing Reason: I am the author of this document
 Signing Time: 25-Aug-2025 | 12:42:12 PM PDT

FDB933D53ED64FE4A2CA0C1E014EFC91