



**Beckman Coulter, Inc.**  
Certificate of Analysis

Product Name DxH 500 Series Diluent

Reorder Part Number B36845

Lot Number 0153170

Storage Conditions 2 - 30°C

Date of Manufacture 2025-08-20  
(yyyy-mm-dd)

Expiration Date 2026-08-14  
(yyyy-mm-dd)

TEST	SPECIFICATION	RESULTS	PASS/FAIL
General Inspection	Components and Labeling must be complete and correct	Meets specification	Pass

PERFORMANCE TEST			
pH	Within Acceptance range: 6,80 - 7,20	Meets specification	Pass
Osmolality	Within Acceptance range: 325 - 335 mOsm/kg	Meets specification	Pass
Particle Count	Within Acceptance range: <= 250/mL	Meets specification	Pass

Approval	Title	DATE (yyyy-mm-dd)
O. Spangenberg	Quality Assurance	2025-08-22

B36845 AC