




**MIAMI AND KENTUCKY REAGENT MANUFACTURING
CERTIFICATE OF ANALYSIS**

IDENTIFICATION:	Product Name:	DXH DILUENT
	Part Number:	628017
	Lot Number:	3558650
	Date of Manufacture:	05-NOV-2025
	Expiration Date:	06-MAY-2027

TEST	SPECIFICATION	RESULTS	Pass/Fail
pH	7.0 – 7.2	Meets Specifications	Pass
Osmolality	331 – 341 mOsm/KG	Meets Specifications	Pass
Background Count	≤ 25/0.1 mL	Meets Specifications	Pass
Sodium	222 – 236 mEq/L	Meets Specifications	Pass
Conductivity	19.00 – 20.00 mS/cm @25 C	Meets Specifications	Pass

The product was manufactured in accordance with the US FDA Quality System Regulation (GMP).

APPROVAL	TITLE	DATE
Ramona Thompson 	Quality Control Manager	10-NOV-2025