






 Beckman Coulter Ireland, Inc.
 Lismeehan
 O'Callaghan's Mills
 Co. Clare, Ireland
 www.beckmancoulter.com

 Quality Control
 Diagnostics

Commision No. : 180153740 Customer : US Only

	:	800-7212	
	:	7212755	M
	:	2026-03-23	
	:	2027-03-23	
		30°C	
		2°C	

Container, Test Strips and Label

Visual control of container, test strips, label and residual moisture




Parameter	: Visual Control of container, test strips and label
Results	: Meets Requirement

Parameter	: Residual Moisture of Stopper	Parameter	: Residual Moisture of Desiccant Pack
Results	: 3.37 %	Results	: 5.37 %
Upper Limit	: 8.00 %	Upper Limit	: 8.00 %
Lower Limit	: - %	Lower Limit	: - %

Functional Test with Test Urine

Test Strips must meet the criteria as per IRL-UAQC-PCD-0005

Parameter	: Functional Testing with Test Urine
Results	: Meets Requirement

Signature:	 <u>EAMONN LENIHAN</u>	Date:	<u>22 Apr 2026</u>
Verified By:	 <u>KEVIN O'SULLIVAN</u>	Date:	<u>22 Apr 2026</u>
Released by:	 <u>EMMA KING</u>	Date:	<u>22 Apr 2026</u>