



**KREFELD REAGENT MANUFACTURING  
CERTIFICATE OF ANALYSIS**

IDENTIFICATION:	Product Name:	DXH DILUENT
	Part Number:	628017
	Lot Number :	256514K
	Date of Manufacture:	2026-04-28 (yyyy-mm-dd)
	Expiration Date:	2027-10-26 (yyyy-mm-dd)

TEST	SPECIFICATION	RESULTS	PASS/FAIL
pH	7.00 - 7.20	Meets Specification	Pass
Osmolality	331 - 341 mOsm/kg	Meets Specification	Pass
Background Count	≤ 250/mL	Meets Specification	Pass
Sodium	222 - 236 mEq/L	Meets Specification	Pass
Conductivity	19.00 - 20.00 mS/cm @ 25°C	Meets Specification	Pass

The product was manufactured in accordance with the EN ISO 13485:2016.

Approval	Title	DATE (yyyy-mm-dd)
O. Spangenberg	Quality Assurance	2026-04-30